



Mar 12, 2021

Making Quenching and Blocking Solution

Kenneth Schackart¹, Kattika Kaarj¹¹University of Arizona

1

Works for me

dx.doi.org/10.17504/protocols.io.bmnfk5bn

Yoon Lab

Kenneth Schackart
University of Arizona

SUBMIT TO PLOS ONE

ABSTRACT

This protocol describes how to make Quenching solution and Blocking solution used for protein EDAC particle coupling. Glycine serves as a primary amine source for quenching the coupling reaction. Bovine serum albumin (BSA) is used to coat particle surface to prevent nonspecific (e.g. hydrophobic) interaction. Adapted from Bangs Laboratories Inc. TechNote 205 Covalent Coupling.

ATTACHMENTS

[bangs_laboratories_technote205_covalent_coupling.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.bmnfk5bn

PROTOCOL CITATION

Kenneth Schackart, Kattika Kaarj 2021. Making Quenching and Blocking Solution. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bmnfk5bn>



KEYWORDS

quenching solution, glycine, edac, protein coupling

LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 23, 2020

LAST MODIFIED

Mar 12, 2021

PROTOCOL INTEGER ID

42407

MATERIALS TEXT

MATERIALS

 Glycine Fisher

Scientific Catalog #G46


 BSA Sigma

Aldrich Catalog ##A8806


Reagents:

- Bovine serum albumin (BSA)
- Glycine

Make Quenching Solution



- 1 Add  **100 mg** BSA to  **100 mL** DW to make **[M]0.1 % w/v** solution.
Swirl to dissolve BSA.

Recommended concentration from Tech Note 205 is **[M]0.05 % w/v** - **[M]1 % w/v** BSA. Volume and concentration can be adjusted as necessary.

- 2 Add  **263 mg** glycine to make **[M]35 Milimolar (mM)** solution.

Recommended concentration from Tech Note 205 is **[M]30 Milimolar (mM)** - **[M]40 Milimolar (mM)** glycine. Volume and concentration can be adjusted as necessary.

Make Blocking Solution

- 3 Add  **1 g** BSA to  **100 mL** DW to make **[M]1 % w/v** solution.
Swirl to dissolve BSA.

Store Solutions

- 4 Store solutions at  **4 °C**.