



JAN 31, 2024

# GFP-Clu-tails purification from Escherichia coli cells

Andreas Bracher<sup>1</sup>, F Ulrich Hartl<sup>1</sup>

<sup>1</sup>Department of Cellular Biochemistry, Max Planck Institute of Biochemistry



Patricia Yuste-Checa

## ABSTRACT

This protocol details how to purify the fusion protein GFP-Clu-tails from Escherichia coli.

## ATTACHMENTS

[GFP-Clu-tails purification protocol\\_protocols.io.docx](#)

## MATERIALS

### Buffers

Binding buffer: PBS + [M] 20 millimolar (mM) imidazole

Elution buffer: PBS + [M] 250 millimolar (mM) imidazole

OPEN ACCESS



DOI:

[dx.doi.org/10.17504/protocols.io.n2bvj3w55lk5/v1](https://dx.doi.org/10.17504/protocols.io.n2bvj3w55lk5/v1)

**Protocol Citation:** Andreas Bracher, F Ulrich Hartl 2024. GFP-Clu-tails purification from Escherichia coli cells. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n2bvj3w55lk5/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Jan 25, 2024

Last Modified: Jan 31, 2024

PROTOCOL integer ID: 94463

**Funders Acknowledgement:**

Aligning Science Across  
Parkinson's  
Grant ID: ASAP-000282

## His<sub>6</sub>-Ubiquitin-GFP-Clu-tails expression and cell lysis

1 Express His<sub>6</sub>-Ubiquitin-GFP-Clu-tails in E. coli BL21 (DE3) codon+RIL cells cultured in 4 L LB Medium



Overnight at 18 °C with 0.25 millimolar (mM) IPTG.

2 Centrifuge culture and keep pellet.



3 Re-suspend pellet in 70 mL volume of ice-chilled Binding buffer, add Complete protease inhibitor cocktail (Roche) and 1 millimolar (mM) phenylmethylsulfonyl fluoride (PMSF).

4 Lyse cells by ultrasonication in ice bath (15 cycles of 00:00:25 ultrasonication with 00:01:35 intermittent cooling). 2m

5 Clear lysate by centrifugation at 22000 rpm in a JA25.50 rotor at 4 °C.

## Ni<sup>2+</sup>-chelating Sepharose metal affinity chromatography


6 Load supernatant onto a 8 mL Ni-chelating Sepharose (Cytiva 17-0575-01) column previously equilibrated with binding buffer by gravity flow at 4 °C .

7 Wash the column with 5 CV of ice-chilled Binding buffer.






8 Wash the column with 2 CV of ice-chilled PBS + 50 Molarity (m) imidazole.

9 Elute His6-Ubiquitin-GFP-Clu-tails protein with 6x 5 mL of ice-chilled Elution buffer. Collect fractions of 5 mL volume. Bright green color indicates presence of His6-Ubiquitin-GFP-Clu-tails. Store fractions On ice .



10 Analyze eluted fractions by SDS-PAGE and Coomassie blue staining.




## Protease cleavage and removal of protease

- 11 Pool fractions containing His6-Ubiquitin-GFP-Clu-tails peak (here fractions 4 and 5). Cleave off His6-Ubiquitin with His-tagged Usp2 during  04:00:00  On ice . 4h
- 12 Exchange protein buffer to Binding buffer with a HiPrep™ 26/10 Desalting column (Cytiva 17-5087-01) equilibrated with Binding buffer at  4 °C using bright green fluorescence as a marker for protein containing fractions.
- 13 Remove uncleaved material and His-tagged Usp2 by metal affinity chromatography as above. Collect flow-through and wash fractions with bright green fluorescence.
- 14 Concentrate to  1.5 mL volume by ultrafiltration using 10 kDa cut-off spin concentrator at  4 °C .

## Size exclusion chromatography on HiLoad 16/600 Superdex-200 (Cytiva 28...

- 15 Apply concentrate on a HiLoad 16/600 Superdex-200 (Cytiva 28-9893-35) column equilibrated with PBS. Develop the column at  4 °C and collect  3 mL fractions.
- 16 Analyze eluted fractions by SDS-PAGE and Coomassie blue staining.



- 17 Merge fractions 16 and 17. Concentrate to 0.7 mL volume by ultrafiltration using 10 kDa cut-off spin concentrator at 4 °C , aliquot and flash-freeze purified GFP-Clu-tails in liquid nitrogen for storage at -70 °C .

#### Note

Concentrations were determined by absorbance at 488 nm (eGFP!) using absorbance coefficients of 61,000 M<sup>-1</sup> cm<sup>-1</sup> or 1.45 L g<sup>-1</sup> cm<sup>-1</sup>.

**Approximate yield:** From 4 L of culture around 2.5 mg of GFP-Clu-tails were obtained.