



🌐 CHEM 584--Cloning sgRNA Sequences into lentiCRISPRv2

Ken Christensen¹¹Brigham Young University

protocol .



Ken Christensen
Brigham Young University

This protocol describes how to clone your chosen sgRNA sequence into the lentiCRISPRv2 plasmid. A method to validate your cloning is also included here.

Ken Christensen 2021. CHEM 584--Cloning sgRNA Sequences into lentiCRISPRv2. **protocols.io**
<https://protocols.io/view/chem-584-cloning-sgrna-sequences-into-lenticrisprv-bznyp5fw>



 protocol ,

Nov 01, 2021

Nov 02, 2021

54712

Design and Order Guide Sequences

- 1 Design and order sgRNA sequence oligonucleotides based on your previously identified target sequences.

Oligo 1: 5'-CACCGNNNNNNNNNNNNNNNNNNNN

Oligo 2: 3'- C N N N N N N N N N N N N N N N N N N N C A A A

Important Note: Do not include the NGG PAM in your designed oligonucleotides

lentiCRISPRv2 Digestion

- 2** Digest 2 μ g of lentiCRISPRv2 plasmid with BsmBI overnight at 55°C:

Set up the following reaction:

2.5 μ l 10X NEB r3.1 Buffer

x μ l ddH₂O to make up 25 μ l final volume

1 μ l BsmBI

2 μ g lentiCRISPRv2 plasmid*

*Add the plasmid last, into the reaction mixture, pipette gently to mix and spin briefly.

Gel Purify the Digested Plasmid

- 3 Use the Zymoclean Gel DNA Recovery kit to gel purify the digested plasmid.

A ~2 kb filler piece should be present on the gel. **Only purify the larger band.** Leave the ~2 kb band.

Anneal the oligonucleotides

- 4 Suspend the oligonucleotides at 100 μ M in autoclaved ddH₂O

- 5 Set up annealing reaction:

2 μ l Oligo 1

2 μ L Oligo 2

2 μ l 10X T4 Ligation Buffer (NEB)

14 μ l ddH₂O

- 6 Anneal the oligonucleotides in a thermocycler:

95°C for 5 min

Ramp temperature down to 25°C at 5°C/min

- 7 Dilute annealed oligonucleotides 1:200 with autoclaved ddH₂O (1 μ l annealed oligo pair + 199 μ l water).

Ligation into lentiCRISPRv2

- 8 Set up the ligation reaction and incubate at 16°C overnight (with the 2X Instant Sticky End Master Mix, you may be able to significantly shorten this incubation).

x μ l digested plasmid (50 ng)

1 μ l diluted oligonucleotide duplex

6 μ l 2X Instant Sticky End Master Mix

x μ l ddH₂O to make a final volume of 12 μ l

Transform into NEB STBL competent cells

- 9 Add up to 5 μ l of your ligation reaction into NEB STBL Mix & Go component cells following the

Mix & Go Competent Cell protocol. Plate your cells on Amp/Carb plates.

Use Colony PCR to screen for the presence of cloned oligo pairs

- 10 Use the Colony PCR protocol to screen for the presence of cloned oligo pairs using the 2X TaqDog Master Mix.

Run 10-20 µl of your PCR products on a 0.8% agarose gel. Use the 100 bp ladder for your gel.

Primers (provided):

lentiCRISPRv2 For: GTGGAAAGGACGAAACACCG

lentiCRISPRv2 Rev: CTAGGCACCGGATCAATTGC

Expected amplicon for positive clones = 248 bp

Isolate plasmid for positive clones

- 11 Prepare overnight cultures for positive clones.

Purify plasmid from overnight cultures using the Zymo Plasmid Miniprep-Classical protocol.

Sequence positive plasmids

- 12 Send plasmid for Sanger Sequencing (Eton Biotechnology).

Sequencing primer:

hU6-F: GAGGGCCTATTTCCCATGATT