



Jul 15, 2020

Huh7.5_SOP

Melissa Teusel¹, Irina Titkova¹, Ina Schmitt¹, Lena Postawa¹, **Marcel Schilling¹**, Ursula Klingmüller¹

¹Division Systems Biology of Signal Transduction, German Cancer Research Center (DKFZ), INF 280, Heidelberg, Germany

1 Works for me dx.doi.org/10.17504/protocols.io.biapkadn

Marcel Schilling

ABSTRACT

SOP for thawing, cultivation and freezing of the Huh7.5 cell line by the Division Systems Biology of Signal Transduction, German Cancer Research Center (DKFZ).

DOI

dx.doi.org/10.17504/protocols.io.biapkadn

PROTOCOL CITATION

Melissa Teusel, Irina Titkova, Ina Schmitt, Lena Postawa, Marcel Schilling, Ursula Klingmüller 2020.
Huh7.5_SOP. **protocols.io**
dx.doi.org/10.17504/protocols.io.biapkadn

KEYWORDS

Huh7.5 cell lines, cell culture, medium, thawing, freezing

LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 07, 2020

LAST MODIFIED

Jul 15, 2020

PROTOCOL INTEGER ID

38959

MATERIALS

NAME	CATALOG #	VENDOR
Fetal Bovine Serum	10270106	Gibco - Thermo Fischer
Neubauer Improved (NI) Hemocytometer	22-600- 100	Life Technologies
0.05% Trypsin-EDTA, phenol red	25300054	Invitrogen - Thermo Fisher
Glutamax (100x)	35050-061	Gibco - Thermo Fischer
DMEM, high glucose, no glutamine, no phenol red	31053028	Thermo Fisher
Bovine Serum Albumin	A9418	Sigma Aldrich
Penicillin-streptomycin (P/S)	15140122	Gibco - Thermo Fisher
Dimethylsulfoxide – DMSO	#41639-100ML	Merck Millipore Sigma
DPBS w/o: Ca and Mg	P04-36500	PAN Biotech
15 ml centrifuge tubes	188271	greiner bio-one
50 ml centrifuge tubes	227261	greiner bio-one
Tissue Culture Dish D:150 mm	93150	Techno Plastic Products (tpp)
Tissue Culture Test Plates 6-well	92006	Techno Plastic Products (tpp)

NAME	CATALOG #	VENDOR
CRYO.S 2 ML PP ROUND BOTTOM INTERNALTHREAD	122263-2DG	greiner bio-one
Cool Cell LX	210004	Biozym
Mr. Frosty	5100-0001	Thermo Fisher Scientific
Trypan Blue solution 0.4% for microscopy	935395	Merck Millipore Sigma

Cells

- 1 Huh7.5 cells were kindly provided by Charles M. Rice (The Rockefeller University, NY, RRID:CVCL_7927)

- 1.1 Cells are cultivated in incubator at 37°C and 5 % CO₂ incubation and 95 % relative humidity

Medium

- 2 **Growth Medium**
DMEM
1% Glutamax
10% FCS
1% P/S
 - 3 **Growth factor-depleted Medium**
DMEM
1% Glutamax
1% P/S
1 mg/ml BSA
 - 4 **Freezing Medium**
70% Growth-Medium
20% FCS
10% DMSO

Thawing of cells

- 5 Pre-warm growth medium in **37 °C** water bath
 - 6 Thaw the cells in a **37 °C** water bath (not completely, there should be a small visible ice clump inside the tube)
 - 7 Transfer the cells into a 15 ml centrifuge tube containing 9 ml pre-warmed growth medium
 - 7.1 To this aim, use 1000 µl pipette and add some medium to cryotube, and mix. With this step the little ice clump should disappear
 - 8 Transfer all liquid from the cryotube to the centrifuge tube, rinse cryotube again and transfer to centrifuge tube

9  **1000 rpm, Room temperature 00:03:00 , 130 x g**

10 Aspirate supernatant

11 Resuspend cells in 10 ml fresh growth medium (~7 times)

12 Transfer the cells to a 150-mm dish (~ 2×10^6) in 25 ml growth medium

13 Let cells adhere to the surface of the dish over night in incubator

14 Replace the medium with 25 ml fresh growth medium the next day

15 Split after 2-3 days as described in "Cultivation of cells" (check for confluency to be around 80% to 90%; for example: thawed on Thursday, change medium on Friday and split on Monday)

Cultivation of cells

16 Cells are used for experiments until passage 25 (total passage number starting from the stocks provided by Charles M. Rice, The Rockefeller University)


17 Cells are split every 3 to 4 days (e.g. Monday to Friday and Friday to Monday)

18 On one ~80-90% confluent 150-mm plate, $7-10 \times 10^6$ cells can be expected

19 Aspirate the growth medium

20 Wash the cells with 10 ml DPBS once

21 Detach the cells from surface by adding 0.05% Trypsin-EDTA (3 ml for 150-mm dish)

- 22 Place the cells with Trypsin-EDTA in an incubator and incubate for 3 to 4 minutes by which the cells should be detached from the surface
- 23 Tap the dish carefully and stop the enzyme reaction of Trypsin-EDTA by adding growth medium containing FCS (7 ml for 150-mm dish)
- 24 Resuspend the cells and transfer to centrifuge tube
- 25  **1000 rpm, Room temperature 00:03:00 , 130 x g**
- 26 Aspirate supernatant
- 27 Resuspend in 10 ml fresh growth medium (~7 times) to remove all clumps
- 28 Dilute 50 µl of cells 1:2 using medium, count the cells applying a 1:2 dilution with Trypan blue using a Neubauer improved hemocytometer
- 29 Plate the cells as required for your planned experiment or for maintaining the cell culture
 - 29.1 Keeping cells in culture:
4 days: 2×10^6 cells in 25 ml growth medium / 150-mm dish
3 days: 3×10^6 cells in 25 ml growth medium / 150-mm dish
 - 29.2 For experiments: usually 0.6×10^6 cells / well of a 6-well plate in 1.5 ml of medium is used
 - 29.3 After 24 h, wash cells 3 times with DPBS , add 1.5 ml of growth-factor depleted medium and incubate for 3 h
 - 29.4 Stimulate cells and perform experiment

Freezing of cells

- 30 Count the cells as described above [go to step #19](#)
- 31 Transfer cell suspension to centrifuge tube
- 32  **1000 rpm, Room temperature 00:03:00 , 130 x g**
- 33 Label cryotubes (name of cell line, number of cells, passage, date and your initials)
- 34 Resuspend cells to a density of 2×10^6 cells / ml in freezing medium
 - 34.1 Calculate needed amount of freezing medium (E.g. if freezing 6×10^6 cells, 3 ml are needed)
 - 34.2 Prepare freezing medium, mix well
 - 34.3 Resuspend pellet in freezing medium using a serological pipette and transfer to cryotubes (1 ml/tube)
- 35 Transfer cryotube to Mr. Frosty or Cool Cell LX, only use room temperature Mr. Frosty or Cool Cell LX, never cold ones
- 36 Transfer cells in cryotubes in Mr. Frosty or Cool Cell LX to -80°C
- 37 Transfer tubes to liquid nitrogen tank on the next day
- 38 Record name of cell line, number of cells, passage, date, your name and location of the tube in liquid nitrogen tank