



1 ▼

Mar 22, 2022

# A three-dimensional immune-oncology model for studying in vitro human NK cell cytotoxic activity V.1

PLOS One ✓

Nontaphat Thongsin<sup>1,2</sup>, Methichit Wattanapanitch<sup>1</sup><sup>1</sup>Siriraj Center for Regenerative Medicine, Research Department, Faculty of Medicine Siriraj Hospital, Mahidol University, Bangkok, Thailand;<sup>2</sup>Department of Immunology, Faculty of Medicine Siriraj Hospital, Mahidol University, Bangkok, Thailand and

1

[dx.doi.org/10.17504/protocols.io.b2ufqetn](https://doi.org/10.17504/protocols.io.b2ufqetn)

PLOS ONE Lab Protocols

Tech. support email: [plosone@plos.org](mailto:plosone@plos.org) methichit.wat

Immunotherapy has emerged as a promising therapeutic approach for treating several forms of cancer. Adoptive cell transfer of immune cells, such as natural killer (NK) cells, provides a powerful therapeutic potential against tumor cells. In the past decades, two-dimensional (2D) tumor models have been used to investigate the effectiveness of immune cell killing. However, the 2D tumor models exhibit less structural complexity and cannot recapitulate the physiological condition of the tumor microenvironment. Thus, the effectiveness of immune cells against tumor cells using these models cannot fully be translated to clinical studies. In order to gain a deeper insight into immune cell-tumor interaction, more physiologically relevant in vivo-like three-dimensional (3D) tumor models have been developed. These 3D tumor models can mimic the dynamic cellular activities, making them much closer to the in vivo tumor profiles. Here, we describe a simple and effective protocol to study the cytotoxic activity of primary human NK cells toward the 3D tumor spheroids. Our protocol includes isolation and expansion of human NK cells, labeling and formation of tumor spheroids, co-culture of NK cells and tumor spheroids, and evaluation of cytotoxic activity using a confocal microscope. This protocol is also applicable to other types of tumors and immune cells.

DOI

[dx.doi.org/10.17504/protocols.io.b2ufqetn](https://doi.org/10.17504/protocols.io.b2ufqetn)<https://doi.org/10.1371/journal.pone.0264366>

Nontaphat Thongsin, Methichit Wattanapanitch 2022. A three-dimensional immune-oncology model for studying in vitro human NK cell cytotoxic activity.  
**protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.b2ufqetn>



#### Thailand Research Fund

Grant ID: RSA6280090

#### Siriraj Research Fund, Faculty of Medicine Siriraj Hospital, Mahidol University

Grant ID: R016333015

#### The National Research Council of Thailand (NRCT)

Grant ID: NRCT5-RGJ63012-126

#### The Distinguished Research Professor Grant

Grant ID: NRCT 808/2563

#### The Program Management Unit for Human Resources & Institutional Development, Research and Innovation

Grant ID: B05F630080

#### The New Discovery and Frontier Research Grant, Mahidol University

Grant ID: R016420010

protocol

Thongsin N, Wattanapanitch M (2022) A three-dimensional immune-oncology model for studying in vitro primary human NK cell cytotoxic activity. PLOS ONE 17(3): e0264366. <https://doi.org/10.1371/journal.pone.0264366>

NK cells, cytotoxicity assay, tumor spheroid formation, immune-oncology model, confocal microscope, killing assay, cancer spheroids, 3D tumor spheroids

protocol ,

Dec 14, 2021

Dec 23, 2021

55911

#### Isolation of peripheral blood mononuclear cells (PBMCs)

1. 30-ml peripheral blood in vacutainer EDTA blood collection tubes
2. Vacutainer® EDTA blood collection tubes (BD Biosciences)

3. IsoPrep (Robbins Scientific)
4. Rosewell Park Memorial Institute (RPMI) 1640 (RPMI-1640) medium (Gibco™)
5. CryoStor® CS10 Cell Freezing Medium (StemCell™ Technologies)
6. 1 × phosphate-buffered saline (PBS)
7. 50-ml conical tubes
8. Hemocytometer
9. Slow cooling cryo-container
10. Cryogenic vials
11. Sterile seropipettes and pipette tips (1000 µl)
12. Transfer pipettes
13. Inverted microscope
14. Tabletop centrifuge
15. Tissue culture hood and humidified incubator (37 °C, 5% CO<sub>2</sub>)
16. Liquid nitrogen tank

### **Negative selection of peripheral blood NK (PB-NK) cells**

1. MojoSort™ Human NK cell Isolation Kit (Biolegend)
2. MojoSort™ Magnet (Biolegend)
3. Sorting buffer: 1 × PBS, 0.5% BSA, and 2 mM EDTA
4. PB-NK culture medium: RPMI-1640, 10% fetal bovine serum (FBS), 2 mM GlutaMAX™, 1% non-essential amino acid, 1% penicillin/streptomycin, and 100 U/ml hIL-2
5. 25-cm<sup>2</sup> cell culture flask
6. Hemocytometer
7. 5-ml round-bottom tube with cap
8. 50-ml conical tubes
9. 70-µm cell strainer
10. Slow cooling cryo-container
11. Cryogenic vials
12. CryoStor® CS10 Cell Freezing Medium (StemCell™ Technologies)
13. Sterile seropipettes and pipette tips (200 and 1000 µl)
14. Inverted microscope
15. Tabletop centrifuge
16. Vortex
17. Tissue culture hood and humidified incubator (37 °C, 5% CO<sub>2</sub>)
18. Liquid nitrogen tank

### **Expansion of PB-NK cells**

1. Irradiated genetically modified membrane-bound IL-21-K562 (mIL-21-K562) cell line
2. mIL-21-K562 cryopreservation medium: 90% FBS and 10% DMSO
3. PB-NK culture medium: RPMI-1640, 10% FBS, 2 mM GlutaMAX™, 1% non-essential amino acid, 1% penicillin/streptomycin, and 100 U/ml hIL-2
4. 25-cm<sup>2</sup> cell culture flask
5. Hemocytometer
6. Slow cooling cryo-container
7. Cryogenic vials
8. Sterile seropipettes and pipette tips (200 and 1000 µl)
9. Inverted microscope
10. Tabletop centrifuge

11. Tissue culture hood and humidified incubator (37 °C, 5% CO<sub>2</sub>)
12. Liquid nitrogen tank

#### **Evaluation of PB-NK cells using flow cytometry**

1. FAC buffer: 3% FBS in 1 × PBS
2. Fluorescence-conjugated anti-human CD45, CD3, CD56 (Biolegend)
3. Zombie Violet™ Fixable Viability Kit (Biolegend)
4. Human AB serum (Sigma-Aldrich)
5. 1% paraformaldehyde (Sigma-Aldrich)
6. Hemocytometer
7. 5-ml round-bottom tube with cap
8. Sterile seropipettes and pipette tips (200 and 1000 µl)
9. Inverted microscope
10. Tabletop centrifuge
11. Flow cytometer (LSRFortessa™)

#### **Culture of KKU213A cell line.**

1. KKU213A culture medium: DMEM/F12, 10% FBS, 2 mM GlutaMAX™, and 1% penicillin/streptomycin
2. 0.1% Trypsin-EDTA
3. 1 × PBS
4. 25-cm<sup>2</sup> cell culture flask
5. 15-ml conical tubes
6. Hemocytometer
7. Sterile seropipettes and pipette tips (200 and 1000 µl)
8. Inverted microscope
9. Tabletop centrifuge
10. Tissue culture hood and humidified incubator (37 °C, 5% CO<sub>2</sub>)

#### **KKU213A labeling using CFSE**

1. Carboxyfluorescein succinimidyl ester (CFSE) (Biolegend)
2. KKU213A culture medium: DMEM/F12, 10% FBS, 2 mM GlutaMAX™, and 1% penicillin/streptomycin
3. 1 × PBS
4. Sterile seropipettes and pipette tips (200 and 1000 µl)
5. Hemocytometer
6. Tabletop centrifuge
7. Tissue culture hood and humidified incubator (37 °C, 5% CO<sub>2</sub>)

#### **Formation of three-dimensional (3D) tumor spheroid**

1. CFSE labeled-KKU213A cell line
2. KKU213A culture medium: DMEM/F12, 10% FBS, 2 mM GlutaMAX™, and 1% penicillin/streptomycin
3. Matrigel® matrix (Corning®)
4. Ultra-low attachment round-bottom 96-well plate (Corning®)
5. Pre-chilled sterile pipette tips (200 and 1000 µl)
6. Tabletop centrifuge

7. Tissue culture hood and humidified incubator (37 °C, 5% CO<sub>2</sub>)
8. Icebox

### **Cytotoxic activity assay**

1. Expanded PB-NK cells
2. PB-NK culture medium: RPMI-1640, 10% FBS, 2 mM GlutaMAX™, 1% non-essential amino acid, 1% penicillin/streptomycin
3. CFSE labeled-KKU213A tumor spheroids
4. Propidium iodide (PI) (Thermo Fisher Scientific)
5. Sterile seropipettes and pipette tips (200 and 1000 µl)
6. Hemocytometer
7. 15-ml conical tubes
8. Tabletop centrifuge
9. Inverted microscope
10. Confocal microscope (Nikon Instruments Inc., Melville, NY, USA)
11. Tissue culture hood and humidified incubator (37 °C, 5% CO<sub>2</sub>)
12. NIS-Elements software

## Isolation and expansion of peripheral blood NK (PB-NK) cells

2w

4h

### **1 Isolation of PB-NK cells by negative selection using NK cell isolation kit**

We isolate peripheral blood mononuclear cells (PBMCs) from 30 ml of peripheral blood collected in vacutainer® EDTA blood collection tubes using density gradient centrifugation.

1. Collect 30 ml of peripheral blood to vacutainer® EDTA blood collection tube and mix the blood thoroughly by inverting.
2. Transfer 15 ml of blood to two 50-ml conical tubes.
3. Add 15 ml of 1× PBS to the blood and mix thoroughly using a sterile serological pipette.
4. Prepare two 50-ml conical tubes containing 9 ml of IsoPrep.
5. Gently overlay 30 ml of diluted blood onto the IsoPrep layer using a sterile serological pipette and centrifuge at 800 × g for 30 min, acceleration = 5 and break = off at room temperature.
6. Remove plasma layer using a transfer pipette.
7. Collect the PBMCs to a 50-ml conical tube using a transfer pipette.
8. Wash the cells with 20 ml of 1× PBS and centrifuge at 400 × g for 10 min at room temperature.
9. Repeat the washing step (step 8).
10. After centrifugation, discard the supernatant and resuspend the cell pellet with 10 ml of basal RPMI-1640 medium.
11. Perform cell counting using trypan blue exclusion assay to determine the cell number and viability.
12. Perform PB-NK cell isolation using the MojoSort™ Human NK cell Isolation Kit following the manufacturer's instruction.
13. Resuspend the cell pellet with the PB-NK culture medium and centrifuge the cells at 500 × g for 10 min at room temperature.
14. During centrifugation, perform cell counting using trypan blue exclusion assay to determine

the cell number and viability.

15. After centrifugation, resuspend the cell pellet with the PB-NK culture medium and collect  $1 \times 10^6$  viable cells for PB-NK cell expansion,  $4 \times 10^5$  viable cells for flow cytometric analysis, and cryopreserve using CryoStor<sup>®</sup> at a density of  $2 \times 10^6$  viable cells/ml/vial.

## 2 Expansion of PB-NK cells

1w 5d

Expansion of the PB-NK cells is performed using the genetically modified mL-21-K562 cell line irradiated with  $\gamma$ -radiation at 10,000 rad. One-day post-irradiation, the mL-21-K562 cell line is cryopreserved in the mL-21-K562 cryopreservation medium.

1. Quickly thaw the cryopreserved irradiated mL-21-K562 cells in a 37 °C water bath until only a small ice crystal remains.
2. Wipe the outside of the cryogenic vial with 70% ethanol.
3. Transfer the cell suspension using a transfer pipette to a 15-ml conical tube containing 9 ml of basal RPMI-1640 medium and take out 10  $\mu$ l of the cell suspension for counting.
4. Centrifuge the cells at  $500 \times g$  for 5 min at room temperature.
5. During centrifugation, perform cell counting using trypan blue exclusion assay to determine the cell number and viability.
6. After centrifugation, discard the supernatant and resuspend the cell pellet with 1 ml of PB-NK culture medium.
7. Co-culture the PB-NK cells (Heading 1, step 15) with the irradiated mL-21-K562 cells at a 1:2 ratio at a seeding density of  $3 \times 10^5$  PB-NK cells/ml in the PB-NK culture medium and incubate at 37 °C, 5% CO<sub>2</sub>.
8. Change the PB-NK culture medium every three days and add freshly thawed irradiated mL-21-K562 cells at 1:1 ratio every six days.

## 3 Evaluation of the sorted PB-NK cells using flow cytometry

2h

We perform flow cytometric analysis to determine the purity of the PB-NK cells after cell sorting and expansion by staining with the common NK cell markers including CD45, CD3, and CD56. NK cells are identified by CD45<sup>+</sup> CD3<sup>-</sup> CD56<sup>+</sup> population.

1. Centrifuge the cell suspension (Heading 1, step 15) at  $500 \times g$  for 5 min at room temperature.
2. Discard the supernatant and block the potential non-specific binding using 100  $\mu$ l of 10% human AB serum in FAC buffer for 30 min at 4 °C.
3. Add 0.2  $\mu$ l of Zombie Violet™ Fixable Viability Kit, vortex briefly, and transfer 50  $\mu$ l of the cell suspension into two 5-ml round-bottom tubes for the NK cell marker and the isotype control staining.
4. Add 1  $\mu$ l of the fluorescence-conjugated anti-human CD45, CD3, and CD56 or the fluorescence-conjugated isotype control to the cell suspension. Incubate for 15 min in the dark.
5. Add 3 ml of FAC buffer, vortex briefly, and centrifuge at  $500 \times g$  for 5 min at room temperature.
6. Discard the supernatant and fix the cells using 500  $\mu$ l of 1% paraformaldehyde. The stained

cells can be stored at 4 °C for up to 2 weeks.

7. Determine the percentage of NK cell population using a flow cytometer (LSRFortessa™) and the FlowJo software. Exclude the dead cells, which stained positive for Zombie Violet™.

#### Formation of three-dimensional (3D) tumor spheroids

3d

- 4 For the 3D spheroid formation, the cholangiocarcinoma cell line (KKU213A) is cultured in the medium containing 2.5% Matrigel® matrix on an ultra-low attachment round-bottom 96-well plate for two days. 30m

#### KKU213A cell line passaging

Passaging is performed when the cells reach 70-80% confluence, usually 3-4 days, and the cells are seeded onto a 25-cm<sup>2</sup> cell culture flask.

1. Aspirate the culture medium and wash the cells with 3 ml of 1 × PBS.
2. Add 1 ml of 0.1% trypsin-EDTA and incubate the cells at 37 °C for 5 min.
3. After incubation, add 4 ml of DMEM/F-12 medium and resuspend the cells vigorously using a serological pipette.
4. Transfer the cell suspension to a 15-ml conical tube, take out 10 µl of the cell suspension for counting, and centrifuge the cells at 500 × g for 5 min at room temperature.
5. During centrifugation, perform cell counting using the trypan blue exclusion assay to determine the cell number and viability.
6. After centrifugation, discard the supernatant and resuspend the cell pellet with 1 ml culture medium.
7. Transfer 1.5 × 10<sup>5</sup> viable cells to a 25-cm<sup>2</sup> culture flask containing 5 ml culture medium and incubate at 37 °C, 5% CO<sub>2</sub>.

#### 5 KKU213A cell line labeling using CFSE

1h

1. Aspirate the culture medium and dissociate the cell using 0.1% trypsin-EDTA (Heading 4, steps 1-6).
2. Transfer 2 × 10<sup>5</sup> viable cells to a 15-ml conical tube, wash with 1 ml of 1 × PBS and centrifuge at 500 × g for 5 min at room temperature.
3. Discard the supernatant and resuspend the cell pellet with 1 ml of 1 × PBS.
4. Add 2.5 µM CFSE solution to the cell suspension, vortex briefly, and incubate at 37 °C, 5% CO<sub>2</sub> for 20 min in the dark.
5. After incubation, add 5 ml of the KKU213A culture medium to the cell suspension and centrifuge at 500 × g for 5 min at room temperature.
6. Discard the supernatant and resuspend the cell pellet with 1 ml of the culture medium.

#### 6 Formation of 3D tumor spheroids

40m

1. Prepare the KKU213A culture medium containing 2.5% Matrigel® matrix and chill the



medium on ice.

2. Aliquot appropriate cell number of the CFSE-labeled KKU213A cells (Heading 5) to a 15-ml conical tube and centrifuge at  $500 \times g$  for 5 min at room temperature.
3. Resuspend the CFSE-labeled KKU213A cells with the pre-chilled KKU213A culture medium containing 2.5% Matrigel<sup>®</sup> matrix.
4. Seed the CFSE-labeled KKU213A cells onto an ultra-low attachment round-bottom 96-well plate at a seeding density of  $2 \times 10^3$  cells/100  $\mu$ l culture medium.
5. Centrifuge the 96-well plate at  $1,000 \times g$  for 10 min at 4 °C and incubate at 37 °C, 5% CO<sub>2</sub> for two days.

Thaw Matrigel<sup>®</sup> at 4 °C overnight and always keep Matrigel<sup>®</sup> under 4 °C until use to prevent gelation. Plasticware that is in contact with Matrigel<sup>®</sup> should be pre-chilled before use.

#### Cytotoxicity assay

4d

- 7 Cytotoxic activity is performed by co-culturing the expanded PB-NK cells with the tumor spheroids at 1:1, 2.5:1, 5:1, and 10:1 (effector: target) ratios for three days. Cytotoxic activity can be determined by the fluorescence intensity of propidium iodide (PI), which stained positive for dead cells, using a confocal microscope. 40m

#### Co-culture of PB-NK cells with KKU213A tumor spheroids

1. Change the medium of the PB-NK cells 24 h prior to co-culture.
2. On day 2 of spheroid formation, transfer the PB-NK cells to a 15-ml conical tube, take out 10  $\mu$ l of the cell suspension for counting, and centrifuge at  $500 \times g$  for 5 min at room temperature.
3. During centrifugation, perform cell counting using trypan blue exclusion assay to determine the cell number and viability.
4. After centrifugation, discard the supernatant, resuspend the cell pellet with 1 ml of PB-NK culture medium, and aliquot the appropriate number of PB-NK cells into a 15-ml conical tube for co-culture.
5. Centrifuge at  $500 \times g$  for 5 min at room temperature.
6. Discard the supernatant and resuspend the cell pellet with an appropriate volume of PB-NK culture medium, without hIL-2 for co-culture, add 2  $\mu$ g/ml of propidium iodide (PI) and mix vigorously using a sterile pipette.

The volume of PB-NK culture medium is 100  $\mu$ l/well.

7. Add 100  $\mu$ l of the suspended PB-NK cells to a well of the ultra-low attachment round-bottom



96-well plate.

1. Gently add the PB-NK cells to avoid disturbing the spheroids. Total volume of culture medium is 200 µl/well and the final concentration of PI is 1 µg/ml.

8. Centrifuge 96-well plate at 1,000 × g for 10 min at 4 °C.
9. Incubate the cells at 37 °C, 5% CO<sub>2</sub> for three days.

We recommend performing 3 replicates for each ratio. Prepare control wells containing only the tumor spheroids for determining the spontaneous tumor cell death, and the tumor spheroids with 0.1% Triton X-100 for determining the maximum cell death.

## 8 Evaluation of cytotoxic activity using a live imaging confocal microscope

3h

1. After three days of co-culture, determine the mean fluorescence intensity (MFI) of PI under a confocal microscope.
2. Cytotoxic activity of the PB-NK cells can be calculated using the following formula:

$$\% \text{Specific killing} = \frac{[(\text{Experimental MFI} - \text{Spontaneous MFI})]}{[(\text{Maximum MFI} - \text{Spontaneous MFI})]} \times 100$$

The experimental MFI represents the mean fluorescence intensity of PI in the presence of effectors at a given effector to target (E: T) ratio, whereas the spontaneous MFI represents the mean fluorescence intensity of PI in the absence of effector cells. Maximum MFI represents the mean fluorescence intensity of PI from tumor spheroids that are treated with 0.1% Triton™ X-100.