



Aug 17, 2022

# CRISPR puromycin

Rebecca Berrens<sup>1</sup><sup>1</sup>IDRM, Paediatrics, Oxford University

1 Works for me

 Share[dx.doi.org/10.17504/protocols.io.yxmvmndpog3p/v1](https://dx.doi.org/10.17504/protocols.io.yxmvmndpog3p/v1)

BerrensLab



Rebecca Berrens

## ABSTRACT

This protocol is for CRISPR knock-in generation in mouse embryonic stem cells

DOI

[dx.doi.org/10.17504/protocols.io.yxmvmndpog3p/v1](https://dx.doi.org/10.17504/protocols.io.yxmvmndpog3p/v1)

## PROTOCOL CITATION

Rebecca Berrens 2022. CRISPR puromycin . **protocols.io**  
<https://protocols.io/view/crispr-puromycin-cfeqtjdw>



## LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Aug 17, 2022

## LAST MODIFIED

Aug 17, 2022

## PROTOCOL INTEGER ID

68784

## BEFORE STARTING

Have ES cells ready growing happily. Mostly 2 passages after thawing

- 1 Plate cells on 6-well plate on SNL a day before transfection (number of cells shown in the table below).<sup>30m</sup>

- 2 Transfect cells using Lipofectamine 2000 protocol and plasmid DNA amounts as shown in the following table:

E14 1x10 <sup>6</sup> cells	E14 1x10 <sup>6</sup> cells
6:1 Molar Ratio dTAG-SETDB1 to MA299/300	No DNA Control
100uL OPTIMEM + 6.46 uL (2.5 ug) dTAG-SETDB1 + 0.44 uL (869.3ng) MA299/300	100ul OPTIMEM
100uL OPTIMEM + 10 uL LF2000	100uL OPTIMEM + 10uL LF2000

Transfection guideline

- 3 A day after transfection cells were passaged from 6-well plate into 10cm Petri dishes with<sup>30m</sup> SNL. Cells from each well were split into 3 Petri dishes in 3 different densities (1/2, 1/3 and 1/5 of the cells).
- 4 24 hours later apply puromycin selection (at the adequate concentration for your cell line of<sup>10m</sup> interest) for 48hours to select for cells that contain the Cas9 plasmid (there is no selection in the targeting vector!).
- 5 After 48 hours of puromycin selection feed cells with regular ES medium every day until ready<sup>10m</sup> to pick colonies (usually 12-14 days after transfection).
- 6 pick colonies (usually 12-14 days after transfection).
  - follow ES cell picking protocol