

Aug 29, 2024

## Startle induced negative geotaxis (SING)

DOI

**[dx.doi.org/10.17504/protocols.io.5jyl82637l2w/v1](https://dx.doi.org/10.17504/protocols.io.5jyl82637l2w/v1)**

Natalie Kaempf<sup>1</sup>, Uli Pech<sup>1</sup>, Patrik Verstreken<sup>1</sup>

<sup>1</sup>VIB-KU Leuven Center for Brain & Disease Research, 3000 Leuven, Belgium.

ASAP Collaborative Rese...



Natalie Kaempf

VIB KU Leuven

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.5jyl82637l2w/v1](https://dx.doi.org/10.17504/protocols.io.5jyl82637l2w/v1)

**Protocol Citation:** Natalie Kaempf, Uli Pech, Patrik Verstreken 2024. Startle induced negative geotaxis (SING). **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.5jyl82637l2w/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** July 29, 2024

**Last Modified:** August 29, 2024

**Protocol Integer ID:** 104252

**Keywords:** ASAPCRN, climbing assay, behavior, dopaminergic dysfunction, SING



**Funders Acknowledgement:**

**Aligning Science Across**

**Parkinson's**

**Grant ID: ASAP-000430**

**EMBO long-term postdoctoral**

**fellowship**

**Grant ID: ALTF\_299-2019**

**Research project, FWO**

**Vlaanderen**

**Grant ID: G0A5219N**

**Research project, FWO**

**Vlaanderen**

**Grant ID: G0B8119N**

**Methusalem project**

**Grant ID: METH/21/05**

**(3M210778)**

**Research project, KU Leuven**

**Parkinson Fonds**

**Grant ID: EQZ-PARFON-O2010**

**Opening the Future grant,**

**Leuven Universiteitsfonds**

**(LUF)**

**Grant ID: EQZ-OPTFUP-O2010**

**Research project, FWO**

**Vlaanderen**

**Grant ID: G031324N**


## Abstract

This protocol describes how to perform the Startle induced negative geotaxis (SING) assay. SING defects in parkinsonism fly models have previously been shown to be caused by dopaminergic impairments.



## Startle induced negative geotaxis (SING)

30s

- 1 for the description of the SING apparatus please refer to:  
Benzer, Seymour. 1967. "BEHAVIORAL MUTANTS OF *Drosophila* ISOLATED BY COUNTERCURRENT DISTRIBUTION." *Proceedings of the National Academy of Sciences* 58 (3): 1112–19. <https://doi.org/10.1073/pnas.58.3.1112>.  
  
Inagaki, Hidehiko K, Azusa Kamikouchi, and Kei Ito. 2010. "Methods for Quantifying Simple Gravity Sensing in *Drosophila Melanogaster*." *Nature Protocols* 5 (1): 20–25. <https://doi.org/10.1038/nprot.2009.196>.
- 2 set a specific time period during the day (e.g. Zeitgeber 3-5) to assess negative geotaxis behavior, which will be kept the same for all following experiments
- 3 flies are transferred into the first tube of the apparatus without prior anesthesia and allowed to adjust  
for some minutes to experimental environment (24±1°C, 40-60% humidity in behavior chamber, note down the environmental condition and time for every genotype)
- 4 tap the flies down five times and allowed to climb up for  00:00:30
- 5 after these 30 s the flies that reached the upper tube are moved to the next tube
- 6 repeat the procedure 4 times
- 7 count the number N of (male/female) flies with the corresponding genotype per tube to calculate the SING score, which is normalized to the mean SING score of the control.
- 8 
$$\text{SING score} = ((N_1 * 0) + (N_2 * 1) + (N_3 * 2) + (N_4 * 3) + (N_5 * 4)) / 4(N_1 + N_2 + N_3 + N_4 + N_5)$$
$$N_k = \text{number of flies in the } k\text{th tube}$$

30s



## Protocol references

Benzer, Seymour. 1967. "BEHAVIORAL MUTANTS OF *Drosophila* ISOLATED BY COUNTERCURRENT DISTRIBUTION." *Proceedings of the National Academy of Sciences* 58 (3): 1112–19. <https://doi.org/10.1073/pnas.58.3.1112>.

Inagaki, Hidehiko K, Azusa Kamikouchi, and Kei Ito. 2010. "Methods for Quantifying Simple Gravity Sensing in *Drosophila Melanogaster*." *Nature Protocols* 5 (1): 20–25. <https://doi.org/10.1038/nprot.2009.196>.