



VERSION 1  
JAN 19, 2023

OPEN  ACCESS

**DOI:**  
[dx.doi.org/10.17504/protocols.io.5qpvp6bv4o/v1](https://dx.doi.org/10.17504/protocols.io.5qpvp6bv4o/v1)

**Protocol Citation:** Rebecca Andrews 2023. Single-molecule Immunofluorescence Tissue Staining Protocol for Oligomer Imaging. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5qpvp6bv4o/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Jan 16, 2023

**Last Modified:** Jan 19, 2023

**PROTOCOL integer ID:**  
75340

**Keywords:**  
immunofluorescences,  
oligomer imaging

# Single-molecule Immunofluorescence Tissue Staining Protocol for Oligomer Imaging V.1

Rebecca Andrews<sup>1</sup>

<sup>1</sup>Univeristy of Cambridge



Rebecca Andrews  
Univeristy of Cambridge

## ABSTRACT

This protocol details about immunofluorescence staining for oligomer imaging.

## ATTACHMENTS

[kb3ib25np.pdf](#)

## GUIDELINES


- Use only clean bottles, flasks, magnetic stirrers, tweezers, weighing spatulas, measuring cylinders – everything should be cleaned, dried and covered if left on the side before next use.
- Everything should be handled with clean tweezers – gloves should not touch the samples, solutions and ideally anything placed into the solutions where the slides are.

## MATERIALS

### Materials and Reagents

- Microtome
- Glass slides
- Xylene solution
- 100% alcohol
- Methanol
- Hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) solution
- Citrate buffer pH6
- Milli Q water
- Pressure cooker
- PBS
- Goat Serum 10%
- AlexaFluor antibody
- 0.1% Sudan black solution
- Vectashield
- Overslip



**Immunofluorescences staining protocol for oligomer imaging**

1 Cut  8 µm tissue sections on a microtome and load onto glass slides.

2 Dry slides  Overnight at  37 °C – **cover over the top.**

10m




3 Before staining commences keep slides for a few hours but ideally  Overnight at  60 °C .


10m




4 De-wax sections through three pots of xylene solution. Use each fresh pots of xylene each time.

4.1 De-wax section through pot of xylene solution for  00:02:00 . (1/3)

2m

4.2 De-wax section through pot of xylene solution for  00:02:00 . (/23)

2m


4.3 De-wax section through pot of xylene solution for  00:02:00 . (3/3)

2m


5 Take sections through two pots of 100% alcohol.

**Note**



Use fresh pots each time – methylated spirits.

5.1 Take sections through pot of 100% alcohol for  00:02:00 . (1/2)

2m

5.2 Take sections through pot of 100% alcohol for  00:02:00 . (2/2)

2m



6 Put slides into methanol + hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) solution (100 ml: 1 ml) for  00:10:00 at  Room temperature **fresh pot each time.**

10m

Note

This process will block any staining of endogenous peroxidase in the tissue sections.


7 Perform necessary antigen retrieval pretreatments by pressure cooking in citrate buffer.

7.1 Pressure cook sections in citrate buffer  6 for  00:10:00 at pressure (wait for it to release high pressure air) in **cleaned pressure cooker.**

10m

7.2 Cool under running **Milli Q** water – *never* under tap water.



8 Block non-specific antigen/antibody binding by placing sections in **PBS and Goat Serum 10% (taken from the DNA PAINT protocol)** for  00:30:00 .


30m

9 Apply primary antibody for  01:00:00 at  Room temperature .

1h


10 Wash in PBS with fresh buffer (**at least filtered if not cell culture grade**).



10.1 Wash  00:05:00 in PBS clean squirty bottle with fresh buffer. (1/3)

5m


10.2 Wash  00:05:00 in PBS clean squirty bottle with fresh buffer. (2/3) 5m

10.3 Wash  00:05:00 in PBS clean squirty bottle with fresh buffer. (3/3) 5m

11 Apply secondary AlexaFluor antibody for  01:00:00 at  Room temperature in the dark. 1h



12 Wash in PBS.



12.1 Wash  00:05:00 in PBS in the dark. (1/3) 5m

12.2 Wash  00:05:00 in PBS in the dark. (2/3) 5m

12.3 Wash  00:05:00 in PBS in the dark. (3/3) 5m

13 Add **filtered (0.22  $\mu$ m)** 0.1% Sudan black solution (0.1% sudan black/70% ethanol) for  00:10:00 at  Room temperature in the dark. 10m



14 Wash 2-3 times in 30% ethanol.



15 Mount section with Vectashield and coverslip (**Plasma cleaned slides**).

16 Take for imaging.

