

Oct 01, 2024

## pMOA 189F-682R

DOI

**dx.doi.org/10.17504/protocols.io.rm7vzz1rvx1w/v1**

Roey Angel<sup>1</sup>, Eva Petrova<sup>1</sup>

<sup>1</sup>Institute of Soil Biology and Biogeochemistry, Biology Centre CAS

Anaerobic and Molecular Microbiology Lab, Biology Centre CAS

Tech. support email: [eva.petrova@bc.cas.cz](mailto:eva.petrova@bc.cas.cz)



**Eva Petrova**

Soil and Water Research Infrastructure

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.rm7vzz1rvx1w/v1](https://dx.doi.org/10.17504/protocols.io.rm7vzz1rvx1w/v1)

**Protocol Citation:** Roey Angel, Eva Petrova 2024. pMOA 189F-682R. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.rm7vzz1rvx1w/v1>

**Manuscript citation:**

Holmes AJ, Costello A, Lidstrom ME, Murrell JC. Evidence that particulate methane monooxygenase and ammonia monooxygenase may be evolutionarily related. FEMS Microbiol Lett. 1995;132:203–8.

<http://www.sciencedirect.com/science/article/pii/037810979500311R>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** May 24, 2018

**Last Modified:** October 01, 2024

**Protocol Integer ID:** 12426

**Keywords:** PCR, pmoA, PMMO, Methanotrophs, Methane oxidising bacteria, Particulate methane monooxygenase



## Abstract

A PCR assay targeting the alpha subunit of the particulate methane monooxygenase gene (*pmoA*). This is a general assay for methane oxidising bacteria using primers 189F–682R.

**Note:** these primers might occasionally also amplify the *amoA* gene of ammonia-oxidising bacteria.

189F GGN GAC TGG GAC TTC TGG (target seq. 189 - 207) Holmes et al. (1995),

682R GAA SGC NGA GAA GAA SGC (target seq. 664 - 682) Holmes et al. (1995).

Fragment size: 508bp.



## PCR mixture

1	Reagent	Final. conc.	1 tube (50 $\mu$ l)	1 tube(25 $\mu$ l)	96 tubes (25 $\mu$ l x100)	96 tubes (12.5 $\mu$ l x100)
	PCR H <sub>2</sub> O		37.9	18.95	1255	627.5
	10X Buffer (I for RT and II for DNA)	1x	5	2.5	500 (Promega)	250 (Promega)
	dNTP mixture (2.0 mM each)	0.2mM	5	2.5	250	125
	MgCl <sub>2</sub> (50 mM)*	3mM	1.5*	0.75*	250	125
	BSA (20 $\mu$ g/ $\mu$ l)	0.8 $\mu$ g/ $\mu$ l	2	1	100	50
	<b>189f (25<math>\mu</math>M)</b>	0.4 $\mu$ M	0.8	0.4	40	20
	<b>682r (25 <math>\mu</math>M)</b>	0.4 $\mu$ M	0.8	0.4	40	20
	AccuPrime™ Taq	?	1	0.5	15	7.5
	Template		1	0.5	0.5 x 100	0.25 x 100
	<b>Final volume</b>		<b>50</b>	<b>25</b>	<b>2500</b>	<b>1250</b>

\* Buffer contains 1.5 mM MgCl<sub>2</sub> at final concentration

## PCR program

- 2
1. 94°C – 4'
2. x **10** {
  - a. 94°C – 1'
  - b. 62°C – 1°C – 1'
  - c. 68°C – 1'}
3. x **25** {
  - a. 94°C – 1'
  - b. 52°C – 1'
  - c. 68°C – 1'}
4. 68°C – 10'