




MAR 07, 2024

HaloTag autophagy flux assay

 In 1 collection

Dan Tudorica¹

¹University of California, Berkeley



Dan Tudorica
Hurley Lab, QB3, UC Berkeley

ABSTRACT

Pule-chaseable method for quantifying autophagy flux

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.x54v9pwkmg3e/v1

Protocol Citation: Dan Tudorica 2024. HaloTag autophagy flux assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.x54v9pwkmg3e/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Nov 09, 2023

Last Modified: Mar 07, 2024

Day 1: Seeding plates

- 1 Seed 50,000 HeLa cells in each chamber of a 12-well treated tissue culture plate. Maintain cells in DMEM + 10% FBS + 1x Pen-strep.
- 2 Grow for 2 days until cells are ~80% confluent

Day 3: Transient Transfection

- 3 If cells do not need to be transfected, skip to "day 5"
- 4 Wash cells into DMEM + 10% FBS without antibiotic prior to transfection.
- 5 Prepare reagents for lipofection with lipofectamine 3000. Each well requires 1.2 µg of pmCherry vector or construct. Prepare a master mix consisting of 1.2 µg of DNA, 125 µL of opti-mem low serum medium, and 2.4 µL of P3000 reagent per well. Prepare an additional mix consisting of 125 µL of opti-mem and 2.4 µL of lipofectamine 3000 per well.
- 6 Combine the DNA + P3000 and lipofectamine solutions together, and allow to incubate for 10 minutes at room temperature.
- 7 After 10 min, dispense 250 µL of combined DNA/P3000/lipofectamine mixture into each well to be transfected.

- 8 Return plates to incubator, and incubate for two days before conducting experiment

Day 4: Swap medium

- 9 Swap medium of transfected cells with fresh DMEM + 10% FBS + Pen-strep

Day 5: Autophagy Induction and Cell Collection

- 10 Replace growth medium with labeling medium consisting of growth medium + 100 nM Janelia Fluor 549-conjugated HaloLigand.
- 11 Incubate at 37 C for 20 min.
- 12 Aspirate off labeling medium and wash with PBS or growth medium.
- 13 Wash cells into either growth medium + DMSO, growth medium + 10 μ M Oligomycin/ 5 μ M Antimycin, or Earle's Balanced Salt solution to start the experiment.
- 14 Return cells to incubator for either 3 H (for starvation autophagy measurements via EBSS) or 6 H (for mitochondrial depolarization via OA).

- 15 Collect cells via scraping into ice-cold PBS.
- 16 Pellet cells via tabletop centrifuge, aspirate off supernatant, and snap freeze in liquid nitrogen. Store at -80 C until next day.

Day 6: HaloTag gel shift measurement and quantification

- 17 Lyse cells via resuspension in 30 μ L lysis buffer and 30 minute incubation at 4 C (50 mM HEPES pH 7.4, 150 mM NaCl, 0.5% NP-40, 1 mM TCEP, cOmplete protease inhibition tablet (Roche)).
- 18 Measure protein concentration via a 96-well plate BCA assay. Briefly, pipette 200 μ L of reconstituted BCA reagent (Thermo Fisher) in each well of a 96-well plate. Include 5 additional wells for a 3 mg/mL, 2 mg/mL, 1 mg/mL, PBS blank, and Lysis buffer blank. Into each 200 μ L well, pipette 1 μ L of sample, standard, or blank. Cover plate and incubate at 37 C for 30 min before measuring the absorbance of each well at 560 nm. If protein concentration is especially high, it might be necessary to make a dilute sample of each well sample prior to performing BCA assay.
- 19 In excel, subtract blank readings from measured absorbances, and plot a standard curve using the BSA standards. Using the standard curve, determine the concentration of each protein samples, and the volume necessary to load 20 μ g each.
- 20 Prepare samples for SDS-PAGE via addition of loading buffer and 25 mM TCEP. Heat samples at 95 C for 2 min before loading onto gel. Omit a ladder.
- 21 Run gel at 120 V until the dye band has reached the bottom of the gel. Gently cut off dye band at the bottom of the gel, and transfer the gel to a clean container. Rinse in deionized water before imaging on a ChemiDoc using the AlexaFluor 548 channel. Use AutoExpose protocol unless trying to image particularly faint bands.

- 22** Open image file in Fiji, define rectangular selection that encompasses the first lane, and use the Gel Analyzer tool to highlight each lane before plotting the trace.
- 23** In the trace window, define a baseline, and then measure the area under each section of the curve that describes the band of interest. The highest MW band should represent the Su9-halo or the LC3-halo reporter, while the lowest MW band should represent processed halo.
- 24** Calculate the fraction of total Halo signal constituted by the lowest band over the total halo signal to quantify autophagy flux.