

Oct 24, 2021

DNA extraction from mouthwash samples

V.1

Ahmed A Shibl¹, Anique Ahmad¹, Tsedenia Denekeew¹, Aashish Jha¹¹New York University, Abu Dhabi

protocol .

Ahmed Shibl
New York University, Abu Dhabi

DNA extraction

Ahmed A Shibl, Anique Ahmad, Tsedenia Denekeew, Aashish Jha 2021. DNA extraction from mouthwash samples . **protocols.io**
<https://protocols.io/view/dna-extraction-from-mouthwash-samples-by9mpz46>



protocol ,

Oct 20, 2021

Oct 24, 2021

54285

QIAGEN - DNeasy PowerSoil Pro Kit

Ice

Vortex

BeadBeater

Centrifuge

Liquid N2

Set centrifuge to 4°C

Keep CD2 solution on ice

Prepare and label collection tubes, microcentrifuge tubes, and MB spin columns

1 Thaw mouthwash samples on ice for  00:30:00

30m

Transfer desired volume into 1.5 mL or 2 mL eppendorfs

1.1 1-2 mL


2 Centrifuge transferred samples at maximum speed at 4 °C 10m


2.1 Repeat step 1.1 if more volume is needed to see a pellet or a high yield of DNA is required

3 Discard the supernatant carefully without disturbing the pellet

4 Add  800 µL CD1 and vortex to resuspend pellet


5 Transfer entire eppendorf content to PowerBead Pro tubes











6 Secure PowerBead Pro tubes onto the bead beater and run at maximum speed for  00:05:00 5m

7 Centrifuge the PowerBead Pro tubes at maximum speed for  00:01:00 1m

8 Transfer the supernatant carefully without disturbing the pellet, into a clean 2 mL microcentrifuge tube

9 Add  200 µL CD2 into the 2 mL microcentrifuge tube and vortex for  00:00:10 10s

10 Centrifuge the 2 mL microcentrifuge tube at maximum speed for  00:01:00 1m

- 11 Transfer the supernatant carefully without disturbing the pellet, into a clean 2 mL microcentrifuge tube
- 12 Add  **600 μ L** CD3 into the 2 mL microcentrifuge tube and vortex for  **00:00:10** ^{10s}
- 13 Load  **650 μ L** of the lysate onto MB spin columns and centrifuge at maximum speed for  **00:01:00** ^{1m}
- 13.1 Discard flow-through and repeat step 13 to consume all the lysate
- 14 Place spin column onto new collection tube, add  **500 μ L** EA, and centrifuge at maximum speed for  **00:01:00** ^{1m}
- 15 Discard flow-through and place spin column back into the collection tube
- 16 Add  **500 μ L** CD5 onto spin columns and centrifuge at maximum speed for  **00:01:00** ^{1m}
- 17 Discard flow-through and place spin column back into the collection tube
- 18 Centrifuge at maximum speed for  **00:02:00** and place spin column into 1.5 mL eppendorf tube ^{2m}
- 19 Add  **50-100 μ L** nuclease-free water to the center of the spin column and leave at room temperature for around  **00:07:00** ^{7m}
- 20 Centrifuge at maximum speed for  **00:01:00** , measure with Qubit and flash freeze DNA ^{1m}

with liquid N₂