



VERSION 2

JAN 31, 2024

OPEN ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.36wgq4enyvk5/v2

External link:

<https://doi.org/10.1073/pnas.2121362119>

Protocol Citation: Diep R Ganguly 2024. Plant RNA purification using TRIzol (TRI reagent). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.36wgq4enyvk5/v2> Version created by Diep R Ganguly

MANUSCRIPT CITATION:

🌐 Plant RNA purification using TRIzol (TRI reagent) V.2

Diep R Ganguly¹

¹University of Pennsylvania

Gregory lab



Diep R Ganguly

University of Pennsylvania, The Australian National Universi...

ABSTRACT

Extraction of total RNA from plant tissue using TRIzol (or TRI reagent) followed by DNA nuclease treatment.

ATTACHMENTS

[TRIreagent_SDS.pdf](#) [chloroform_SDS.pdf](#)

GUIDELINES

TRI reagent and chloroform are hazardous - handle with care, operate in fume hood, wear vinyl gloves and safety glasses.

MATERIALS

- TRIzol (or TRI reagent)
- Chloroform (or chloroform : IAA [24:1])
- Isopropanol
- 80% Ethanol
- Nuclease-free H₂O (e.g. 0.01% DEPC-treated H₂O) or Tris-EDTA (10 mM Tris-Cl, pH 6.5, 0.1 mM EDTA)
- 2 mL safe-lock Eppendorf tubes
- 1.5 mL microcentrifuge tubes
- Tissue lyser or mortar and pestle
- RNase-Free DNase Set (Qiagen)
- RNA loading dye, 2X (NEB)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jan 30, 2024

Last Modified: Jan 31, 2024

PROTOCOL integer ID: 94427

SAFETY WARNINGS



TRI reagent - hazardous

Chloroform - hazardous


Ensure you read SDS documents (attached) and organise appropriate waste vessels (fume hood).

BEFORE START INSTRUCTIONS

Ensure benches and equipment are RNase free.

For RNase-Free DNase Set (Qiagen): Prepare DNase I stock solution by dissolving the lyophilized DNase I (1500 Kunitz units) in 550 µl of the RNase-free water. Mix gently by inverting the vial. Divide into single use aliquots and store at -20 °C. Thawed aliquots can be stored at 2–8°C for up to 6 weeks. Do not refreeze the aliquots after thawing.

RNA purification

- 1 Collect 50-100 mg of plant tissue and freeze immediately in liquid N₂.
- 2 Grind tissue to fine powder under liquid N₂ using tissue lyser or mortar + pestle, then immediately add 1 mL TRI reagent (1 mL per 100 mg tissue).

Note, achieving a fine grind is critical to high yields of intact RNA.
- 3 Invert each tube by hand ~20x and incubate at room temperature for 5 minutes (DO NOT vortex samples as it may result in RNA degradation).
- 4 Add 1/5 volume of pre-mixed chloroform:isoamyl alcohol (24:1), cap tubes, shake vigorously (by hand) for 15 seconds (solution should become cloudy), then incubate at room temperature for 3 minutes.

5 Centrifuge at 14,000 rcf for 10 minutes at 4°C.

6 Transfer the upper aqueous phase to a new microfuge tube (approx. 400-600 µL).

7 Repeat steps 4 and 5 (approx. 300-400 µL).

Note, if you are observing buffer and salt carryover in your purified RNA (high 230 nm absorbance), reduce volume of upper-phase recovered.

8 (Optional) If the expected RNA concentration is ≤ 10 µg/mL, add 1/10 volume of 3M NaOAc (pH 5.5) and/or 100 µg/mL GlycoBlue (or glycogen).



9 Add equal volume of 100% isopropanol, then mix by inverting tubes ~20x by hand.

10 Incubate at -20°C for 1 hour. Alternatively, incubate overnight to capture small RNAs.

11 Centrifuge samples at 14,000 rcf for 20 minutes at 4°C.

12 Remove the supernatant, you should observe a white pellet.

- 13 Add 1 mL of 80 % ethanol and invert tube ~10x.
- 14 Centrifuge samples at 10,000 rcf for 5 minutes at room temperature.
- 15 Remove supernatant, carefully since the pellet often becomes dislodged at this step.
- 16 Air-dry pellet at room temperature for 5 minutes.
- 17 Resuspend pellet in RNase-free water (e.g. 0.01% DEPC-treated water) or Tris-EDTA (10 mM Tris-Cl, pH 6.5, 0.1 mM EDTA).

DNA nuclease treatment and ethanol precipitation

- 18 Make up volume of RNA solution to 87.5 μ L with nuclease-free water. This can be performed with an aliquot or total sample from the previous step.
- 19 Add 10 μ L Buffer RDD and 2.5 μ L DNase I stock solution (Qiagen RNase Free DNase Set) and mix with gentle pipetting.
- 20 Incubate at room temperature for 5-10 minutes.

- 21** Add 500 μ L of 100% ethanol.
- 22** (Optional) Add 3 μ L glycogen or GlycoBlue, and 10 μ L NaOAc (pH 5.5) to aid RNA precipitation.



- 23** Mix by gentle inversion.
- 24** Incubate samples for at least 1 hour at -20 °C (overnight, if purifying small RNAs).
- 25** Centrifuge at 14,000 rcf at 4 °C for 20 minutes.
- 26** Remove supernatant and rinse pellet with 1 mL of 80% ethanol.
- 27** Centrifuge at 10,000 rcf for 5 minutes at 4 °C.
- 28** Remove supernatant without disturbing pellet and air-dry for 2 minutes.

- 29** Resuspend pellet in RNase-free water (e.g. 0.01% DEPC-treated water) or Tris-EDTA (10 mM Tris-Cl, pH 6.5, 0.1 mM EDTA).

Quality control

- 30** Take 50-100 ng aliquot of RNA and mix 1:1 with 2X RNA loading dye (NEB).
- 31** Incubate RNA at 65 °C for 5 minutes.
- 32** Load and run samples on a 1% agarose TBE gel.
- 33** Nanodrop RNA to check for purity and quantity.