



Ovary tissue dissociation V.2

Regina Hoo¹, Carmen Sancho¹, Roser Vento-Tormo¹

¹Wellcome Sanger Institute

Version 2 ▾

Mar 27, 2021

1

Works for me

dx.doi.org/10.17504/protocols.io.656hg9e

Vento-Tormo

Tech. support email: rv4@sanger.ac.uk

Regina Hoo

SUBMIT TO PLOS ONE

ABSTRACT

Single cell suspension from ovary tissue and digestion of oocytes

DOI

dx.doi.org/10.17504/protocols.io.656hg9e

PROTOCOL CITATION

Regina Hoo, Carmen Sancho, Roser Vento-Tormo 2021. Ovary tissue dissociation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.656hg9e>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 06, 2019

LAST MODIFIED

Mar 27, 2021

PROTOCOL INTEGER ID

27550

MATERIALS TEXT

MATERIALS

[Beta-mercaptoethanol](#) Contributed by users

[Falcon® Conical Tubes, 50 mL 500 Tubes](#) Stemcell

Technologies Catalog #38010

[DNAse](#)

Sigma Catalog #4716728001

[Buffer TCL](#)

Qiagen Catalog #1031576

[PBS Invitrogen - Thermo Fisher](#)

[HypoThermosol® FRS Preservation solution](#) Sigma

Aldrich Catalog #H4416

[100 µm Cell](#)

Strainer Falcon Catalog #352360

[Liberase™ TM Research Grade](#) Sigma

Aldrich Catalog #5401119001

[Collagenase IA](#) Sigma

Aldrich Catalog #C2674-100MG

[RPMI 1640 Medium](#) Thermo Fisher

Scientific Catalog #11875093

[10X RBC Lysis Buffer \(Multi-](#)

species) **eBioscience Catalog #00-4300-54**

[HBSS calcium magnesium no phenol red](#) Thermo Fisher

Scientific Catalog #14025050

[DMEM/F-12 HEPES](#) Thermo Fisher

Scientific Catalog #11330032

Prepare collagenase mix

1

Heat-inactivate FBS at **56 °C** for **00:30:00** before use.

Collagenase mix recipe:

| Product | Stock | Final volume (10 ml) | Concentration |
|----------------------------|---------------------------------|----------------------|---------------|
| RPMI or Hams F12 + 10% FBS | 9 ml RPMI or Hams F12 + 1ml FBS | 8.8 ml | |
| Collagenase IA | 10 mg/ml | 1 ml | 1 mg/ml |
| Liberase TM | 5 mg/ml | 100 µl | 50 µg/ml |
| DNase I | 10 mg/ml | 100 µl | 0.1 mg/ml |

Prepare oocyte lysis buffer

2

| Product | Stock | Final volume (1 ml) | Concentration |
|----------------------|-------|---------------------|---------------|
| TCL lysis buffer | | 990 ul | |
| beta-mercaptoethanol | | 10 ul | 1% |

*aliquote in low bind PCR tube 10ul/tube

Tissue dissociation and digestion

- 3 Wash tissue with PBS.
- 4 Place wet tissue under a petri dish. Take 2 scalpels and roughly mince up the tissue. This step is crucial to increase the efficiency of the digestion.
- 5 Transfer contents to 50ml falcon containing the collagenase mix (~ **10 mL** /tissue but it will depend on the size of the tissue.
- 6 Tighten lid and then seal with parafilm.
- 7 Incubate at **37 °C** for **00:45:00**. Shacking every 10 min during the incubation is recommended.
- 8 Filter sample through small strain (100 um) and discard the filter.

*Note: If collecting oocytes, proceed to step 18.
- 9 Centrifuge 450 g, **00:05:00** (0.5 rcf, 5mins).
- 10 Remove carefully 90% of the media by decanting the tube (leave between 0.5 ml and 1 ml).
- 11 Resuspend the cell pellet in 5 ml of PBS.
- 12 Check cell viability with trypan blue staining and on a with disposable C-chip haematocytometer.

If cells is viable proceed to next step.
- 13 Centrifuge 450 g for 5min. Discard supernatant.

- 14 Resuspend the cell pellet in 5 ml of PBS. Centrifuge 450 g for 5min. Discard supernatant.
- 15 Resuspend cells with 5ml of PBS. Proceed to cell count with disposable C-chip haematocytometer.
- 16 Re-count cells to obtain accurate cell suspension concentration for loading 10X instrument.
- 17 Freeze down excess cells, re-spin down any remaining cells, at 4°C, 450g, 5mins.

Proceed to "Cryopreservation of single cell suspension from tissue"

****If collecting for oocytes**

- 18 Filter sample through small strain (100 um) – do not discard retained tissue that will be put it back in RPMI 10% and observed under the microscope.
- 19 Filtered material: add ~ **5 mL** of RPMI 10%.
- 20 Transfer oocytes manually using high resolution microscope into drops of oocyte media.
- 21 Clean oocytes.
- 22 1. Digest oocytes (to be optimised).
Options:
 - A) Trypsin-EDTA (0.25%) phenol red (cat. #: 25200072, LifeTech).
 - B) Accutase
 - C) Collagenase Mix
- 23 Transfer oocytes and granulosa cells into **10 µl** of lysis buffer in a PCR tube.
- 24 Centrifuge the other non-oocyte cells at 450 g, **00:05:00** (0.5 rcf, 5mins)
- 25 Cells are resuspended in 2 ml of freezing media and aliquot in 2 cryopreservation tubes (1ml/each). Refer to "Cryopreservation of single cell suspension from tissue" protocol.