



1 ▼

May 30, 2022

VOC and VOI (SARS-Cov-2) identification by Sanger Sequencing V.1

Laís Ceschini¹, Matheus Filgueira Bezerra²,Viviane do Carmo Vasconcelos de Carvalho³, Cássia Docena³,Gabriel Luz Wallau⁴, Marcelo Henrique Santos Paiva⁵¹Departamento de Entomologia, Instituto Aggeu Magalhães (IAM);²Departamento de Microbiologia, Instituto Aggeu Magalhães (IAM);³Núcleo de Plataforma Tecnológica (NPT);⁴Núcleo de Bioinformática, Instituto Aggeu Magalhães, Fiocruz, Departamento de Entomologia, Instituto Aggeu Magalhães (IAM);⁵Departamento de Entomologia, Instituto Aggeu Magalhães (IAM), Núcleo de Ciências da Vida, Universidade Federal de Pernambuco

1

dx.doi.org/10.17504/protocols.io.ewov1nxqkgr2/v1

Rede Covid FIOCRUZ-PE



Laís Ceschini

The method hereby described is an update of a rapid and accessible protocol based on Sanger sequencing that is able to discriminate the main SARS-CoV-2 VOCs (Variants of Concern) and VOIs (Variants of Interest), according to each characteristic mutational signature at the Spike receptor binding domain (RBD) and an additional mutational profile of the N-terminal domain (NTD) of the Spike protein. Although this approach does not substitute whole-genome sequencing, in a scenario that combines the rapid spread of new VOCs around the world with supply shortages and lack of technical infrastructure, it represents a powerful tool that allows a broader network of laboratories to perform molecular surveillance of SARS-CoV-2 VOCs, improving its capacity to report more results within in a timely manner.

DOI

dx.doi.org/10.17504/protocols.io.ewov1nxqkgr2/v1

Laís Ceschini, Matheus Filgueira Bezerra, Viviane do Carmo Vasconcelos de Carvalho, Cássia Docena, Gabriel Luz Wallau, Marcelo Henrique Santos Paiva 2022. VOC and VOI (SARS-Cov-2) identification by Sanger Sequencing .

protocols.io<https://dx.doi.org/10.17504/protocols.io.ewov1nxqkgr2/v1>

CNPq

Grant ID: 303902/2019-1

genotyping, RNA virus, Spike protein, molecular assay, screening

protocol ,

May 30, 2022

May 30, 2022

63521

cDNA Synthesis

2h 15m

1

The cDNA was prepared according to the manufacturer's instructions: High-Capacity cDNA Reverse Transcription Kit

[Applied Biosystems™ High-Capacity cDNA Reverse Transcription Kit](#) **Applied Biosystems Catalog #4368814**

Mix the following components in an 0.2mL 8-strip tube or 96 well PCR plate;

Component	Value
10X RT Buffer	2.0 µL
dNTP Mix (100mM)	0.8 µL
10X RT Random primers	2.0 µL
MultiScribe Reverse Transcriptase	1.0 µL
H ₂ O	4.2 µL
Template RNA	10.0 µL
Total	20.0 µL

2

2h 15m

Incubate the reaction as follows:

Time **Temperature**

🕒 00:10:00 at 🌡 25 °C

🕒 02:00:00 at 🌡 37 °C ,

🕒 00:05:00 at 🌡 85 °C

Hold at 🌡 4 °C

Primers sequences

3

Primer sets targeting the Spike residue binding domain (RBD).

A	B	C	D
Primer set	Flanked region*	Amplicon	Covered mutations
Artic primers 76 Left 5'-AGGGCAAACCTGGAAAGATTGCT-3' 77 Right 5'-CAGCCCCTATTAAACAGCCTGC-3'	22819-23500	725 bp	N439K, L452R/Q, Y453F, S477N, T478K, E484K/Q, S494P, N501Y/T/S, A570D, D614G
In house 1MS Fw 5'-TAACGCCACCAGATTTGCAT-3' 2MS Rv 5'-ACACGCCAAGTAGGAGTAAGT-3'	22607-23446	878 bp	K417T/N, N439K, L452R/Q, Y453F, S477N, T478K, E484K/Q, S494P, N501Y/T/S, A570D, D614G

*not including the primer binding site.

Primer set targeting the Spike N-terminal domain (NTD).

A	B	C	D
Primer set	Flanked region*	Amplicon	Covered mutations
Artic primers 71 Left 5'-ACAAATCCAATTCAGTTGTCTTCCTATTC-3' 73 Right 5'-CACCAGCTGTCCAACCTGAAGA-3'	21386-22324	989 bp	Δ69-70/144-145, Δ157-158, ✖241-243, S13I, L18F, T19R, T20N, P26S, Q52R, A67V, V70F, G75V, T76I, D80A, T95I, D138Y, W152C, E156G, R190S, D215G, A222V, W258L.

*not including the primer binding site.

PCR amplification 11m 55s

4

The PCR was performed under conditions standardized using the

[Taq Platinum DNA Polymerase Invitrogen - Thermo Fisher](#)




Mix the following components in an 0.2mL 8-strip tube or 96 well PCR plate;

Component	Value
10x Buffer	2.5 µL
MgCl ₂	0.5 µL
dNTP (10 mM)	1.0 µL
Forward primer (10uM)	0.5 µL
Reverse primer (10uM)	0.5 µL
Taq Polymerase	0.25 µL
H ₂ O	18.25 µL
cDNA input	1.5 µL
Total	25 µL

5 Incubate the reaction as follows:

11m 55s

Step	Time	Temperature	Cycle
Initial denaturation	00:05:00	98 °C	1x
Denaturation	00:00:30	98 °C	35x
Annealing	00:00:35	59 °C	35x
Extension	00:00:50	72 °C	35x

Final extension	 00:05:00	 72 °C	1x
Hold	Indefinite	 4 °C	

Electrophoresis, quantification and Sequencing

6

- Agarose gel was prepared at **0.1 mg/mL** and stained with Sybr Safe (Sigma-Aldrich).

- PCR products was quantified using

Nanodrop 2000C

Thermo Scientific TSC-ND2000C

(1uL per sample) and diluted to a final concentration of 30 ng/uL.

- Sequencing reaction is performed with BigDye Terminator v3.1 (Applied Biosystems) and run in capillary electrophoresis (ABI 3500, Applied Biosystems), according to the manufacturer's instructions.