

Aug 14, 2020

Human Islet Microvasculature Immunofluorescence in Optically Cleared Samples

Martha Campbell Thompson¹, Elizabeth Butterworth Hosaka², Katelyn N Carty²¹University of Florida; ²co-author**1** Works for me dx.doi.org/10.17504/protocols.io.y3tfynnJesus Peñaloza
University of Florida

ABSTRACT

This protocol describes the immunostaining performed on five human control pancreas samples, cleared using a modified passive CLARITY (PACT) method according to our published methods. Islets were identified based on alpha-cells stained with glucagon while collagen-containing basement membranes in the extracellular matrix were stained using anti-collagen IV. Smooth muscle cells surrounding arteries and arterioles were identified using anti-SMA primary antibody directly conjugated with Cy3 fluorophore. An accompanying protocol describes using Vesselucida360 to contour islets and determine islet basement membrane and SMA densities and morphometric variables (<https://dx.doi.org/10.17504/protocols.io.bjfkj6>).

EXTERNAL LINK

<https://www.protocols.io/view/human-pancreas-pact-optical-clearing-and-high-reso-9gbh3sn/materials>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Butterworth E, Dickerson WD, Vijay V, Weitzel K, Cooper J, Atkinson EW, Coleman JE, Otto KJ, Campbell-Thompson M. High resolution 3D Imaging of the Human Pancreas Neuro-Insular Network. J Vis Exp. 2018 Jan 29;(131). doi: 10.3791/56859.2018. PMID: 29443037.

ATTACHMENTS

[jove-protocol-56859-high-resolution-3d-imaging-of-the-human-pancreas-neuro-insular-network.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.y3tfynn

PROTOCOL CITATION

Martha Campbell Thompson, Elizabeth Butterworth Hosaka, Katelyn N Carty 2020. Human Islet Microvasculature Immunofluorescence in Optically Cleared Samples . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.y3tfynn>

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Butterworth E, Dickerson WD, Vijay V, Weitzel K, Cooper J, Atkinson EW, Coleman JE, Otto KJ, Campbell-Thompson M. High resolution 3D Imaging of the Human Pancreas Neuro-Insular Network. J Vis Exp. 2018 Jan 29;(131). doi: 10.3791/56859.2018. PMID: 29443037.

EXTERNAL LINK

<https://www.protocols.io/view/human-pancreas-pact-optical-clearing-and-high-reso-9gbh3sn/materials>

KEYWORDS

human islet, glucagon, optical clearing, confocal microscopy, 3D, PACT, collagen IV, SMA, Vesselucida360, alpha-cells, pancreas, neuro-insular network, passive CLARITY

LICENSE

———— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

IMAGE ATTRIBUTION

An islet is shown from case 6232 with all 3 fluorescent channels visible.

CREATED

Mar 12, 2019


LAST MODIFIED

Aug 14, 2020

PROTOCOL INTEGER ID

21331

GUIDELINES

This protocol is for PACT-cleared human pancreas samples, fixed with 4% paraformaldehyde for at minimum  **48:00:00**, 4% agarose embedded and sectioned at 400 um or greater thickness. The Triton-X concentration in the antibody diluent was increased from 0.1% to 0.5% based on experiments that showed antibody penetration into the middle of the sample was improved without loss of signal.


MATERIALS

NAME	CATALOG #	VENDOR
Glucagon (GCG) Mouse anti-human	ab10988 (Lot# GR290488-2)	Abcam
Goat anti-ms Dylight 405 hi-cross	35500BID	Thermo Fisher Scientific
PBS Phosphate Buffered Saline 10X Solution	BP399-1	Fisher Scientific
Goat normal serum	S-1000	Vector Laboratories
Triton X-100	T8787-100	Sigma Aldrich
Collagen IV (Col14) Rabbit anti-human	ab6586	Abcam
Smooth muscle actin- alpha (a-SMA) Mouse anti-human-Cy3 [clone 1A4]	C6198	Sigma – Aldrich
Goat anti-rb AF488 hi cross	A-11034	Thermo Fisher Scientific

STEPS MATERIALS

NAME	CATALOG #	VENDOR
PBS Phosphate Buffered Saline 10X Solution	BP399-1	Fisher Scientific
Goat normal serum	S-1000	Vector Laboratories
Triton X-100	T8787-100	Sigma Aldrich

MATERIALS TEXT





Vesselucida Explorer

v2019.1.1 [↗](#)

by MBF Bioscience

EQUIPMENT

NAME	CATALOG #	VENDOR
LSM 710	LSM 710	

SAFETY WARNINGS

Follow all laboratory safety procedures when handling hazardous chemicals and materials.

DISCLAIMER:

Different primary and secondary antibody lots may differ in affinities and should be independently optimized.

BEFORE STARTING

Ensure samples are cleared and the SDS is well-removed by extensive washing in PBS at [↗ Room temperature](#) on a rocker plate.

- 1 For these samples, optical clearing was performed using 4% SDS at [↗ 60 °C](#) with rocking for up to 3 weeks until the samples reached transparency. Our original protocol used 8% SDS. We have subsequently tested both 4% and 8% SDS at [↗ Room temperature](#) and [↗ 60 °C](#). Different tissues should be tested independently for optimal results.

Prepare the Antibodies

- 2 Make antibody dilution buffer (ABD): 1x PBS containing 2% normal goat serum (NGS) and 0.5% Triton-X 100.



PBS Phosphate Buffered Saline 10X Solution

by Fisher Scientific

Catalog #: BP399-1



Goat normal serum

by Vector Laboratories

Catalog #: S-1000



Triton X-100

by Sigma Aldrich

Catalog #: T8787-100

- 3 Dilute the primary antibodies to 1:200 in ABD and prepare at least 1 mL per sample.

Stain the Tissue

- 4 Incubate the samples with primary antibodies on a rocker plate for 48:00:00 1-5 days at Room temperature .

- 5 Wash the tissue sections with 1xPBS for five times, minimum 01:00:00 at Room temperature .



PBS Phosphate Buffered Saline 10X Solution

by Fisher Scientific

Catalog #: BP399-1


- 6 Dilute the secondary antibodies to 1:500 and prepare at least 1 mL per sample. Incubate for 48:00:00 or more on a rocker at Room temperature .

- 7 Wash the tissue with 1 x PBS 5 times 01:00:00 at Room temperature .

- 8 Equilibrate the samples in RIMS containing 0.5% sodium azide for at least 12:00:00 before imaging.
Transfer to an imaging dish or chamber slide using RIMS as mounting media.



Imaging

- 9 Perform confocal microscopy. Each microscope will require optimization of imaging parameters.



LSM 710
confocal laser microscope

Zeiss LSM 710 [↗](#)




- 10 The following lasers were used for imaging each antigen in these images:

Track 1- 488- Collagen IV
Track 2- 405- GCG, 561- SMA

Results

11



Islet alpha-cells- GCG Stain - Specific Stain, Intensity Good, Background Low
Collagen IV- basement membranes- Specific Stain, Intensity Good, Background medium
SMA- smooth muscle cells surrounding vessels - Specific Stain, Intensity Good, Background Low

- 12 Analyze islet microvascular density by contouring the GCG+ area for islet volume followed by collagen IV and SMA morphometry using Vesselucida360.