



Cell seeding V.2

PMAT0001¹

¹Nanyang Technological University

Version 2

Oct 21, 2020

1

Works for me

This protocol is published without a DOI.

PMAT0001

PROTOCOL CITATION

PMAT0001 2020. Cell seeding. **protocols.io**
<https://protocols.io/view/cell-seeding-bnp8mdrw>
 Version created by PMAT0001

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 21, 2020

LAST MODIFIED

Oct 21, 2020

PROTOCOL INTEGER ID

43488

MATERIALS

NAME	CATALOG #	VENDOR
1X Dulbecco's Phosphate Buffered Saline (DPBS)	14190094	Thermo Fisher Scientific
Penicillin Streptomycin	15140 122	Invitrogen - Thermo Fisher
DMEM	11885	Invitrogen - Thermo Fisher
Foetal Bovine Serum (FBS) Triple 0.1 µm Sterile Filtered 500 ml Australian Origin	FBS-AU-015	Serana

BEFORE STARTING

- Aseptic techniques***
- Wipe down hood and every item introduced into the hood with 70% ethanol

- 1 Take cells from the freezer and thaw for 00:02:00 in the 37 °C water bath. 2m
- 2 Take a T75 plate and place all needed items in ethanol-sterilized laminar flow hood.
- 3 Draw the needed amount of cells into a 50mL centrifuge tube.
- 4 Add about 7 mL of PBS into the centrifuge tube.

- 5 Spin the mixture in a centrifuge for 🕒 **00:03:00** at 1500rpm.
- 6 Decant supernatant and add PBS to the cell pellet once again.
- 7 Resuspend cell pellet with the PBS using a micropipette.
- 8 Repeat steps 5 to 7.
- 9 Add DMEM (with FBS and antibiotics) into centrifuge tube and resuspend cell pellet.
- 10 Add DMEM to the T75 plate (sufficient to cover the entire surface area, about 🧴 **5 mL**).
- 11 Add 🧴 **5 mL** of the cell mixture into the T75 flask.
- 12 Incubate T75 flask at 🌡 **37 °C** , 5% CO₂.