



MAR 07, 2024

## 🌐 Stable cell line generation via retrovirus

📁 In 1 collection

Dan Tudorica<sup>1</sup>

<sup>1</sup>University of California, Berkeley



Dan Tudorica  
Hurley Lab, QB3, UC Berkeley

### ABSTRACT

Production of pseudotyped virus an subsequent transduction.

OPEN  ACCESS



#### DOI:

[dx.doi.org/10.17504/protocols.io.81wgbxq2qlpk/v1](https://dx.doi.org/10.17504/protocols.io.81wgbxq2qlpk/v1)

**Protocol Citation:** Dan Tudorica 2024. Stable cell line generation via retrovirus. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.81wgbxq2qlpk/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**Created:** Nov 09, 2023

**Last Modified:** Mar 07, 2024

# 1

## **Day 0: Thaw fresh cryovial of cells to be transduced**

Thaw ~ 1 week prior to planned transduction day, use low passage number if possible. Aim for cells to be 80% confluent on day of transduction. Recommend seeding four wells of a 12-well plate (100,000 cells/well) the day before transduction (after collecting virus).

## **Day 1: HEK293T Transfection**

1.  
Prepare 2 80% confluent 10 cm plates of HEK293T for transfection
  
1.  
Prepare 3 mL of warm Opti-mem solution.  
Add 10 µg of retro/lentiviral transfection plasmid, 10 µg VSV-G plasmid, and 10 µg pCMV-MLV (if retroviral)[\[DT1\]](#) [\[DT2\]](#) or 10 µg pCMV R8.74 (if lentiviral)
  
1.  
Add 90 µL of LT-1 reagent and swirl.
  
1.  
Incubate at room temperature for 15 minutes.

1.

Add 1.5 mL dropwise to each 10 cm HEK293T plate.

1.

Incubate for 3 days in incubator before harvesting viral particles. If a fluorescent marker was included, it should begin to be visible on days 1,2 post transfection.

**SAFETY:** After

this point all media and plasticware is to be treated as potentially containing viral particles. Decontaminate using 10% bleach for 10 minutes before discarding as red bag waste. Rinse pipettes with 10% bleach before discarding as biohazard waste.

**Day 4: Harvest viral particles**

1.

Collect media from both HEK293T plates, combine into one 20 mL portion in a 50 mL falcon tube. Bleach and discard used HEK293T plate.

1.

Spin 50 mL tube containing harvested viral particles in Sorvall for 2 min at 2,000 rpm[DT3] [DT4] to clarify media and pellet stray cells

1.

Aspirate off media without disturbing pellet.

1.  
Add 6 mL Lenti-X concentrator solution  
to 18 mL of recovered media (concentrator is 4x). Invert to mix.
  
1.  
Incubate in refrigerator 1 H – O/N<sup>[DT5]</sup>
  
1.  
Spin 50 mL conical containing viral  
particles in centrifuge at 1,500 rcf for 45 min
  
1.  
Remove 2 mL of supernatant from top of  
pelleted solution and save.
  
1.  
Aspirate off remaining media to recover  
viral pellet. Bleach excess media.
  
1.  
Resuspend retro/lenti pellet in 2 mL of  
saved medium.
  
1.  
Titrate concentrated retro/lenti  
solution into prepared 12-well containing target cells. Recommend  
100-200-400-800 µL titration for 4 total attempts at transfection. Typically  
400 or 800 uL will yield ~100% transfection efficiency
  
1.  
Return 12-well plate containing  
transduced cells to incubator.

### **Day 5: Check on transduction**

1.  
May be able to see fluorescence today.

1.  
Swap media to fresh.

### **Day 6-7: Check on transduction**

### **Day 8-:**

1.  
Allow cells to grow until almost confluent on 12-well plate. Trypsinize and replate in a 6-well, allow to grow until confluent, then replate in a 10 cm plate.
1.  
Split into two 10 cm plates, Freeze one plate down into 4 cryovials the next day to save progress. Label with estimated transduction efficiency if able to measure.

1.

Take other plate and select with antibiotics if able to/necessary, or clonally isolate populations if able to/necessary.