



Dec 24, 2020

CalculatingGenerationTime

Elizabeth Fozo¹¹In-house protocol

1

Works for me

This protocol is published without a DOI.

Eadewunm

ABSTRACT

Calculation of generation times (doubling times)

PROTOCOL CITATION

Elizabeth Fozo 2020. CalculatingGenerationTime. **protocols.io**
<https://protocols.io/view/calculatinggenerationtime-bq2fmybn>

KEYWORDS

calculation , generation times

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 24, 2020

LAST MODIFIED

Dec 24, 2020

PROTOCOL INTEGER ID

45863

DISCLAIMER:

DISCLAIMER: THIS WORK IS IN PROGRESS. IT IS FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer-reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.


ABSTRACT

Calculation of generation times (doubling times)

Steps

- 1 Grow ON of strain of interest in appropriate media
- 2 In the morning, check the culture(s) under the scope to verify pure culture. Dilute the ON into the appropriate media.

- 3 At hourly intervals, record OD 600 nm. Make sure you record the time in minutes! For example, 2 hours is 120 minutes; 3 hours and 30 minutes is 210 minutes, etc.
- 4 Plot your OD 600 versus time so that you know that you are in EXPONENTIAL PHASE to calculate GENERATION TIMES



- 5 In order to calculate the generation/doubling time, calculate the log of the OD 600 values you have obtained. This will give you negative numbers unless the OD 600 value is above 1.0 (which in that case you are likely NOT in exponential phase).
- 6 Calculate the slope (K) of the line by selecting points in mid to late exponential phase. From the example above:
Slope = $(N - N_0) / (T - T_0) = K$ (slope is called K in these calculations)
- 7 Calculate the generation time using the following formula:
 $t_g = 0.3 / K = \text{minutes to double}$
N.B 0.3 is a constant
- 8 Here are some values as an example:

Time (minutes)	OD 600	log OD 600
180 (T_0)	0.3576 (N_0)	-0.447
291 (T)	0.9525 (N)	-0.021

$$\text{Slope (K)} = (N - N_0) / (T - T_0) = ((-0.021) - (-0.447)) / (291 - 180) = 0.0038$$

$$t_g = 0.3 / K = 0.3 / 0.0038$$

$$t_g = 78.3 \text{ minute}$$