



APR 16, 2024

🌐 siRNA Transfection of Dispersed Islet Cells

Amanda Gomes¹, Xiong Liu¹, Xiaoqing Dai¹

¹University of Alberta



Aliya F Spigelman
University of Alberta

ABSTRACT

siRNA Transfection of Dispersed Islet Cells

MATERIALS

Silencer™ FAM-labeled Negative Control No. 1 siRNA

Catalog number: AM4620 (ThermoFisher)

Opti-MEM™ I Reduced Serum Medium

Catalog number: 31985070 (Gibco)

Lipofectamine™ RNAiMAX Transfection Reagent

Catalog number: 13778075 (ThermoFisher)

Accell Non-targeting Control Pool

Catalog ID: D-001910-10-05

siRNA of your target gene

BEFORE START INSTRUCTIONS

Reconstitute your siRNA to a 20 µM working solution with nuclease-free water.

Work in a cell culture hood, in an RNase-free environment, and keep the hood light off.

Clean the cell culture hood, turn the UV light on, and thaw reagents on ice.

OPEN  ACCESS



Protocol Citation: Amanda Gomes, Xiong Liu, Xiaoqing Dai 2024. siRNA Transfection of Dispersed Islet Cells.

protocols.io

<https://protocols.io/view/sirna-transfection-of-dispersed-islet-cells-czbnx2me>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Aug 30, 2023

Last Modified: Apr 16, 2024

PROTOCOL integer ID: 87118

- 1** Disperse islets, and resuspend cells in Opti-Mem media after centrifugation. Distribute 100 μ L of media with cells per 35 mm dish.

Remember to make extra dishes for qPCR to test the efficiency of the transfection.

- 2** For each control dish (calculate the volume of reagents according to the number of dishes you wish to make):

Tube 1: 50 μ L OptiMem media + 3 μ L Lipofectamine RNAiMAX Transfection Reagent

Tube 2: 50 μ L OptiMem media + 0.6 μ L control siRNA + 0.4 μ L dye indicator (SilencerTM FAM-labeled Negative Control)

For each dish with the target gene silenced:

Tube 3: 50 μ L OptiMem media + 3 μ L Lipofectamine RNAiMAX Transfection Reagent

Tube 4: 50 μ L OptiMem media + 0.6 μ L target siRNA + 0.4 μ L dye indicator (SilencerTM FAM-labeled Negative Control)

Wait 5 minutes after making the 4 tubes.

- 3** Mix tubes 1 and 2, and tubes 3 and 4. Wait 15 minutes.

- 4** Add 103 μ L of transfection reagents to each dish with 100 μ L of media with dispersed cells.

- 5** Wait 4-6h, and add 2 ml of human media (DMEM) to the dishes.

- 6** After 2-3 days, perform desired experiment and add 0.5 -1 mL of Trizol to the extra dishes. Freeze cells with Trizol in the -80°C freezer until RNA extraction.