



Apr 18, 2022

Cell culture, transfection, and imaging

Will Hancock-Cerutti^{1,2,3}, Pietro De Camilli^{1,3}

¹Departments of Neuroscience and of Cell Biology, Howard Hughes Medical Institute, Program in Cellular Neuroscience, Neurodegeneration and Repair, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

²Interdisciplinary Neuroscience Program and MD-PhD Program, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

³Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815

1



dx.doi.org/10.17504/protocols.io.36wgq41j3vk5/v1

William Hancock-Cerutti

This protocol describes general procedures for culturing HeLa cells, transient transfection, and imaging using an Andor Dragonfly spinning disk confocal system.

[dxwubkjz7.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.36wgq41j3vk5/v1

Will Hancock-Cerutti, Pietro De Camilli 2022. Cell culture, transfection, and imaging. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.36wgq41j3vk5/v1>



Cell culture, Transfection, Spinning-disk confocal, Fluorescent microscopy, ASAPCRN

protocol ,

Jul 08, 2021

Apr 18, 2022

Jul 08, 2021

Urmilas

Jul 14, 2021

William Hancock-Cerutti

51390

DMEM Solution:

A	B
FBS	10%
Penicillin	100 U/ml
Streptomycin	100 mg/ml
L-glutamine	2 mM

* (all from Gibco)

General preparation

- 1 Culture the HeLa-M cells at **37 °C** in 5% CO₂ and DMEM containing 10% FBS, **100 U/ml** penicillin, **100 mg/mL** streptomycin, and **2 Milimolar (mM)** L-glutamine (all from Gibco).

Note: For general maintenance, when cells reached 80-90% confluency, they were deattached from the dish with Trypsin and diluted 1:20 in a new dish.

- 2 For live-cell imaging experiments, seed the cells on glass-bottomed dishes (MatTek) at a concentration of 35,000 cells per dish and transfect after **24:00:00** using FuGene HD (Promega) in Opti-MEM (Gibco). 1d

- 3  1d

Image the cells **24:00:00** after transfection.

- 4 Just before imaging, remove the growth medium and replace with prewarmed live-cell imaging solution (Life Technologies).

- 5    30m

For lysotracker experiments, incubate the cells in **[M]50 Nanomolar (nM)** LysoTracker Red DND-99 (ThermoFisher) in complete DMEM for **🕒00:30:00**, wash twice with media, then image in live-cell imaging solution.

6 Perform all live-cell imaging at **🌡37 °C** and 5% CO₂.

7 

Perform spinning-disk confocal microscopy using an Andor Dragonfly system equipped with a plan apochromat objective (63×, 1.4 NA, oil) and a Zyla scientific CMOS camera.

8 For any given experiment, use the same exposure time, laser power, and gain for image acquisition to allow for quantitative comparison.

Imaging of cells stably expressing STING-GFP

3d 14h

9 Generate the cells stably expressing STING-GFP as described elsewhere.

10 Culture the stable STING-GFP HeLa-M cells at **🌡37 °C** in 5% CO₂ and DMEM containing 10% FBS, **📦100 U/ml** penicillin, **[M]100 mg/mL** streptomycin, and **[M]2 Milimolar (mM)** Lglutamine (all from Gibco).

11 

3d

For experiments using siRNA, transfect 60 pmols of the indicated siRNA using **📦6 µL** Lipofectamine RNAiMax (ThermoFisher) in Opti-MEM (Gibco) per dish according to manufacturer protocol. Image the cells **🕒72:00:00** after siRNA transfection.

12 

14h

For experiments using cGAMP, transfect **📦50 µg/L** of cGAMP using **📦18 µL** Lipofectamine RNAiMax (ThermoFisher) in Opti-MEM (Gibco) per dish according to manufacturer protocol. Image the cells **🕒14:00:00** after transfection.

