

May 14, 2024

Phage screening - colorimetric test

DOI

dx.doi.org/10.17504/protocols.io.5jyl82jmdl2w/v1

Kesia da Silva¹, Jason Andrews¹

¹Stanford University



Kesia da Silva

Stanford University

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.5jyl82jmdl2w/v1

Protocol Citation: Kesia da Silva, Jason Andrews 2024. Phage screening - colorimetric test. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.5jyl82jmdl2w/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 13, 2024

Last Modified: May 14, 2024

Protocol Integer ID: 99704

Keywords: Bacteriophage, detection, colorimetric assay, environmental surveillance

Abstract

Phage Screening - Colorimetric Test protocol involves a series of steps designed to detect the presence of specific bacteriophages through a colorimetric assay. This comprehensive procedure aims to provide a robust method for screening and assessing bacteriophage activity within microbial cultures, essential for various research and diagnostic applications.



Attachments




Colorimetric Phage a...

822KB

Materials

1. 15ml Falcon tube
2. 0.5% sodium hypochlorite
3. 70% ethanol
4. Cotton/Tissue paper
5. Gloves
6. Tube holder
7. Discard bag
8. Permanent marker
9. 0.22µm syringe filter
10. 5ml Syringe
11. LB broth
12. Pipette sets
13. Sterile tips - 200µl 1000µl
14. Disposable Serological Pipets, Sterile, 25 mL

Protocol materials

 50ml AlamarBlue Cell Viability Assay Reagent **G-Biosciences Catalog #786-923**

Step 10



Day 1

- 1 Prepare LB-agar plates.
- 2 Prepare LB-broth.
- 3 To recover bacteria from stock open the tube and use a sterile loop, toothpick or pipette tip to scrap off some bacteria and streak on appropriate media plates using the overnight culture and grow at 37 °C .
- 4 If using lyophilized strain skip step 3 and 5 and [go to step #6.1](#)

Day 2

- 5 Next, start overnight cultures by inoculating fresh liquid media (5 mL) with a single colony and grown for 14-16 hours at 37 °C in a shaker incubator.


Note







Note: Make sure you are using only one colony to start the experiment.

Day 3


- 6 Inoculate 5 mL of fresh LB (1% inoculum) and grow at 37 °C until the cells reach a density of $OD_{600}=0.2$. The amount of host culture to add is not critical, but it should be enough that would normally produce a saturated culture over the period of the enrichment (0.2 OD is enough or 0.5 McFarland Standard i.e. 1.5×10^8 CFU/ml).










**Note**

If using lyophilized strain  [go to step #6.1](#)


- 6.1 Using a pipette, aseptically add  500 μL of the recommended growth medium to the freeze-dried material and mix well.
- 6.2 Transfer the entire suspension to a test tube containing  5 mL of the recommended media.
- 7 In a new falcon tube add  1 mL of filtered  water sample to  900 μL of LB-broth and inoculate with  100 μL of a fresh host culture (0.2 OD) and mix gently.


Note

If using lyophilized strain take  100 μL of host culture previously diluted in step 6.2.

- 7.1 In another falcon tube inoculate a positive control adding 10-50uL (this amount will depend of the concentration of your phage stock) of phage stock to  1900 μL of LB-broth with  100 μL of a fresh host culture and mix gently.
- 8 Incubate this enrichment culture at  37 °C for at least 2 hours (or under whatever conditions the host favors). At the same time inoculate a new tube with  5 mL of fresh LB (1% inoculum) and grow at  37 °C for 2 hours.
- 9 Filter the culture through a 0.22 μm syringe filter. The filtrate is now ready for testing.
- 10 In a new falcon tube add  100 μL of sample (filtrate from previous step) to  1800 μL of LB-broth and inoculate with  100 μL of a fresh host culture (0.2 OD) and add  50 μL



of  50ml AlamarBlue Cell Viability Assay Reagent **G-Biosciences Catalog #786-923** (resazourin) and mix gently.

- 11 Incubate the culture at  37 °C for 3 hours. Check the culture every hour to monitor color change.
- 12 If the sample is positive for phage culture will remain blue. If the sample is negative culture will change to pink color.