




Jun 23, 2022

In vitro germination of *Austropuccinia psidii* urediniospores

Alyssa M Martino¹, Rebecca M Degnan²¹School of Life and Environmental Sciences, The University of Sydney, Camperdown NSW 2006, Australia;²School of Chemistry and Molecular Biosciences, The University of Queensland, Brisbane City QLD 4072, Australia

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.3byl4b73jvo5/v1 Alyssa Martino

ABSTRACT

Optimisation of *Austropuccinia psidii* urediniospore germination for use in RNA extraction and cytogenetics.

DOI

dx.doi.org/10.17504/protocols.io.3byl4b73jvo5/v1

PROTOCOL CITATION

Alyssa M Martino, Rebecca M Degnan 2022. In vitro germination of *Austropuccinia psidii* urediniospores . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.3byl4b73jvo5/v1>



FUNDERS ACKNOWLEDGEMENT

Australian Research Council
Grant ID: LP190100093

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jan 09, 2022

LAST MODIFIED

Jun 23, 2022

Spore collection, desiccation, and storage

2d

- 1 Harvest fresh spores from heavily infected leaves by shaking into a paper bag.
- 2 Move fresh spores to a glass petri dish with no lid. Transfer petri dish to a desiccator with silica gel beads for 24 - 48 hours to dry the spores.
- 3 Sieve spores to remove any plant matter or dirt.
- 4 Aliquot spores into Nunc CryoTubes or similar for short or long term storage.

Spores can be frozen at -80°C or in liquid nitrogen for long-term storage but may reduce germination rates. Use of spores immediately after desiccation is recommended.

Inocula preparation and plating

2h

- 5 Create 2% water agar plates.
- 6 In a 15 mL centrifuge tube, make up inocula to a concentration of approximately 1 mg/mL in sterile distilled water (SDW) with 5ul/mL Tween 20 (0.05%). Mix thoroughly by gently inverting the tube several times.

A haemocytometer can be used to count and confirm spore concentration.

- 7 Rehydrate spores in inocula for 30 minutes prior to use.

- 8 Pipette 2 uL aliquots of inocula onto agar plate, covering the entire plate with distinct 2 uL droplets and seal the plates with parafilm.

Incubation 8h

- 9 Germinate in incubator for a minimum of 8 hours in the dark at 18°C and 75% relative humidity.

Determine Germination Percentage 1h

- 10 For each aliquot, count the number of germinated spores under a dissecting microscope. A germinated spore is considered to be a spore with a germination tube at least double the length of the urediniospore.

Germination rates of approximately 70 - 80% are expected when using fresh spores.