



Oct 18, 2020

## PCR

Jiaxin Li<sup>1</sup><sup>1</sup>South China University of Technology

1

Works for me

[dx.doi.org/10.17504/protocols.io.bnienmcbe](https://dx.doi.org/10.17504/protocols.io.bnienmcbe)

Jiaxin Li

DOI

[dx.doi.org/10.17504/protocols.io.bnienmcbe](https://dx.doi.org/10.17504/protocols.io.bnienmcbe)

PROTOCOL CITATION

Jiaxin Li 2020. PCR. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bnienmcbe>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 16, 2020

LAST MODIFIED

Oct 18, 2020

PROTOCOL INTEGER ID

43302

- 1 Set up a small box with ice, put DNA and 2×high Taq Master Mix into it before going into the Bio-rad S1000TM thermocycler cycler.
- 2 Add the following reagent to a PCR tube (50 µl) and program the thermocycler as follows:

## Q5 polymerase

3

component	25ul	50ul
5XQ5 buffer	5ul	10ul
10mM dNTPs	0.5ul	1ul
Forward Primer 10µM	1.25ul	2.5ul
Reverse Primer 10µM	1.25ul	2.5ul
Template (<1000ng)	According to different template	According to different template
Q5 polymerase	0.25ul	0.5ul
ddWater	Add to 25ul	Add to 50ul

4

Temperature	time
98°C	03:00
98°C	00:05-00:10
Tm-5°C	00:10-00:30
72°C	20-30s/kb
72°C	02:00
4°C	∞

KOD FX

5

component	volumn/ul	Final Concentration
2xPCR buffer	25uL	1X
2mM dNTPs	10uL	400umol/L
10umol/L Forward Primer	1.5uL	0.3umol/L
10umol/L Reverse Primer	1.5uL	0.3umol/L
Template DNA	2uL	Genomic DNA≤200ng/50ul
		Plasmid DNA10ng/ul
		cDNA≤200ng
		Crude sample≤0.5-4ul
ddH2O	Add to 50uL	10ul
KOD FX	1uL	1unit/50μl
Total	50μl	

6

Temperature	time
98°C	T1= 03:00
98°C	T2= 00:10
Tm-5°C	T3=00:30
68°C	T4 : 1min/kb
68°C	T6>T2+T3+T4
4°C	∞

T5 Taq Polymerase

7

Component	Volume/ul	Final Concentration
2XT5 Super PCR Mix	5ul	1x
10uM Forward Primer	0.4ul	0.4uM
10uM Reverse Primer	0.4ul	0.4uM
Template DNA	1ul	<1ug
ddH2O	Add to 10ul	
Total	10ul	

Temperature	time
98°C	03:00
98°C	00:10
Tm-5°C	00:10
72°C	20s/kb
72°C	02:00
4°C	∞