

Sep 19, 2024

## Multiplex Nested PCR for Vibrio Cholera

DOI

**[dx.doi.org/10.17504/protocols.io.36wgqznzbkgk5/v1](https://dx.doi.org/10.17504/protocols.io.36wgqznzbkgk5/v1)**

Shannon Fitz<sup>1</sup>, Alex Shaw<sup>1</sup>, Dilip Abraham<sup>2</sup>, michael Owusu<sup>3</sup>

<sup>1</sup>Imperial College London; <sup>2</sup>CMC Vellore; <sup>3</sup>Kwame Nkrumah University of Science and Technology

GEMS - Genomic Environ...



Shannon Fitz

Imperial College London

OPEN  ACCESS



DOI: **[dx.doi.org/10.17504/protocols.io.36wgqznzbkgk5/v1](https://dx.doi.org/10.17504/protocols.io.36wgqznzbkgk5/v1)**

**Protocol Citation:** Shannon Fitz, Alex Shaw, Dilip Abraham, michael Owusu 2024. Multiplex Nested PCR for Vibrio Cholera. protocols.io **<https://dx.doi.org/10.17504/protocols.io.36wgqznzbkgk5/v1>**

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** September 19, 2024

**Last Modified:** September 19, 2024

**Protocol Integer ID:** 108018

**Keywords:** PCR, Amplicon, vibrio cholera, multiplex PCR

**Funders Acknowledgement:**

**The Gates Foundation**

**Grant ID:** INV-049092 PA4273



## Abstract

Primers targeting various regions of the V. Cholera genome were obtained from the literature or designed, and optimized for a multiplex nested PCR assay.

The following primers were extracted from Hoshino et al., 1998: O1\_rfbFor\_Inner (modified here by Dr Alex Shaw), O1\_rfbRev\_Inner, O139\_rfbFor\_Inner, O139\_rfbRev\_Outter.


The following primers were obtained from Theron et al., 2000: ctxA Inner Rev, ctxA\_rev.

The following primers were obtained from Nandi et al., 2000: ctxA\_For, ompW\_For, ompW\_Rev, ompW\_Inner\_Rev.

The following primers were designed by Dr Alex Shaw and are unpublished:

ToxR\_For\_Outter, Tox\_Rev\_Outter, ToxR\_For\_Inner, O1\_rfbFor\_Outter, O1\_rfbRev\_Outter, O139\_rfbFor\_Outter, O139\_rfbRev\_Outter.

## Materials

 DreamTaq PCR Master Mix (2X) **Thermo Fisher Catalog #K1072**

### Approximate total cost per sample: £0.9


17 primers required (cost excluded from estimate as primers do not need to be ordered each time)

Dreamtaq cost per sample: ~£0.45 per sample, per PCR run

### Extra equipment required:

Vortex, mini centrifuge, thermocycler

## Protocol materials

 DreamTaq PCR Master Mix (2X) **Thermo Fisher Catalog #K1072**

Materials, Step 1.1

## Multiplex nested PCR for *V.Cholera*

### 1 Assemble primer pool:

A	B	C	D
Primer name	Primer	Pool C	Pool D
ToxR_For_Outer	AGATGTTTCGGATTAGGACAC	C	
ToxR_Rev_Outer	ATGGCATCGTTAGGGTTAGCA	C	D
ToxR_For_Inner	GCAGCAACGAAAGCCGAATT		D
ctxA_For	CTCAGACGGGATTTGTTAGGCACG	C	D
ctxA_rev	CGATGATCTTGGAGCATTCCCAC	C	
ctxA_InnerRev	GAGTATGGAATCCACCTAAAGC		D
O1_rfbFor_Outer	CCCGACAGCCAGTGAGATAC	C	
O1_rfbRev_Outer	CGTATTGCGGCGGTAAAAGG	C	
O1_rfbFor_Inner	GGTTTCACTGAACAGATGGG		D
O1_rfbRev_Inner	GGTCATCTGTAAGTACAAC		D
O139_rfbFor_Outer	AACGTAGGCACTTGAGAGGC	C	
O139_rfbRev_Outer	TCGCCGGTCGACTGTTTAAC	C	
O139_rfbFor_Inner	AGCCTCTTTATTACGGGTGG		D
O139_rfbRev_Inner	GTCAAACCCGATCGTAAAGG		D
ompW_For	CACCAAGAAGGTGACTTTATTGTG	C	D
ompW_Rev	GAACTTATAACCACCCGCG	C	
ompW_Inner_Rev	GGTTTGTCTGAATTAGCTTCACC		D

Bands for the *V. cholera* reactions:




A	B	C
	First round amplicon length (bp)	Second round amplicon length (bp)
ToxR	875	723
ctxA	444	223
O1_rfbF	1085	196
O139_rfb	1018	451
ompW	586	303

Reconstitute primers to 100  $\mu$ M using nuclease-free water (or if primer manufacturer recommends otherwise, follow their recommendations for reconstituting primers).

Create working stocks of 10  $\mu$ M using nuclease-free water. Create 10  $\mu$ M primer pools C and D by mixing together 10  $\mu$ L of each primer marked as being part of the pool. Scale up as needed.

### 1.1 *V. cholera* first round PCR reaction

 DreamTaq PCR Master Mix (2X) **Thermo Fisher Catalog #K1072**

Prepare the following Master mix on ice in a 1.5ml Eppendorf Lobind tube per number of samples/controls + 10%:

A	B	C
	1 Reaction ( $\mu$ L)	Reactions
DreamTaq 2x master mix	12.5	$\mu$ L
Water	6.5	$\mu$ L
Primer pool C	1	
Total volume	20	

Briefly vortex and centrifuge down the master mix and aliquot 20  $\mu$ L into each PCR tube.

Add 5  $\mu$ L of extracted DNA from each sample to a tube.

Briefly vortex, and centrifuge down the PCR mixes.

### 1.2 Amplify using the following cycling conditions:

A	B	C	D
CYCLE	STEP	TEMP ( $^{\circ}$ C)	TIME



A	B	C	D
1	Initial Denaturation	95	2 minutes
35	Denaturation	95	30 seconds
	Annealing	56	30 seconds
	Extension	72	3 minutes
1	Final Extension	72	10 minutes
-	Hold	10	-

### 1.3 *V. Cholera* second round PCR reaction

Prepare the following Master mix on ice in a 1.5ml Eppendorf Lobind tube per number of samples/controls + 10%:

A	B	C
	1 reaction (µL)	Reactions
DreamTaq 2x master mix	12.5	µL
Water	6.5	µL
Primer pool D	1	µL
Total volume	20	

Briefly vortex and centrifuge down the master mix and aliquot 20 µL into each PCR tube.

Add 5 µL of first round amplicon and vortex briefly to mix. Briefly centrifuge down the PCR mixes.

### 1.4 Amplify using the following cycling conditions:

A	B	C	D
<b>CYCLE</b>	<b>STEP</b>	<b>TEMP (°C)</b>	<b>TIME</b>
1	Initial Denaturation	95	2 minutes
35	Denaturation	95	30 seconds
	Annealing	55	30 seconds
	Extension	72	1 minute



A	B	C	D
1	Final Extension	72	10 minutes
-	Hold	10	-

1.5 Check amplicons using an agarose gel or tapestation.

## Protocol references

Katsuaki Hoshino, Shinji Yamasaki, Asish K. Mukhopadhyay, Soumen Chakraborty, Arnab Basu, Sujit K. Bhattacharya, G. Balakrish Nair, Toshio Shimada, Yoshifumi Takeda, Development and evaluation of a multiplex PCR assay for rapid detection of toxigenic *Vibrio cholerae* O1 and O139, *FEMS Immunology & Medical Microbiology*, Volume 20, Issue 3, March 1998, Pages 201–207, <https://doi.org/10.1111/j.1574-695X.1998.tb01128.x>

Nandi B, Nandy RK, Mukhopadhyay S, Nair GB, Shimada T, Ghose AC. Rapid method for species-specific identification of *Vibrio cholerae* using primers targeted to the gene of outer membrane protein OmpW. *J Clin Microbiol*. 2000 Nov;38(11):4145-51. doi: 10.1128/JCM.38.11.4145-4151.2000. PMID: 11060082; PMCID: PMC87555.

Theron, J., Cilliers, J., Du Preez, M., Brözel, V.S. and Venter, S.N. (2000), Detection of toxigenic *Vibrio cholerae* from environmental water samples by an enrichment broth cultivation–pit-stop semi-nested PCR procedure. *Journal of Applied Microbiology*, 89: 539-546. <https://doi.org/10.1046/j.1365-2672.2000.01140.x>