



Nucleic acid extraction - RNeasy PowerMicrobiome Kit (Qiagen) V.2

Ad  la  de Roguet¹, Shuchen Feng¹

¹UWM

Version 2
Jan 05, 2021

In Development dx.doi.org/10.17504/protocols.io.bpg5mjy6

mclellan lab



McLellan Lab

DOI

dx.doi.org/10.17504/protocols.io.bpg5mjy6

PROTOCOL CITATION

Ad  la  de Roguet, Shuchen Feng 2021. Nucleic acid extraction - RNeasy PowerMicrobiome Kit (Qiagen).
protocols.io
<https://dx.doi.org/10.17504/protocols.io.bpg5mjy6>
Version created by [McLellan Lab](#)

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 07, 2020

LAST MODIFIED

Jan 05, 2021

PROTOCOL INTEGER ID

44285

GUIDELINES

When work is completed, remove equipment and supplies from the cabinet. Wipe the work area with 10% bleach, let stand for 10 min, rinse with water, then with 70% ethanol, and finally with RNAase AWAY.

MATERIALS TEXT

MATERIALS

 RNeasy PowerMicrobiome

Kit **Qiagen Catalog #26000-50**

Mini-Beadbeater-16
high-energy cell disrupter

BioSpec 607 
1 speed

- β -mercaptoethanol
- 0.5 mL free-standing microcentrifuge tubes (low binding)

SAFETY WARNINGS

The nucleic acid extraction has to be performed in a chemical safety to avoid any inhalation of beta-mercaptoethanol.

BEFORE STARTING

1. Clean the working area and all equipment: wipe down with 70% ethanol and let dry. Then, wipe down using RNase AWAY and let dry.
2. Warm RNase-Free Water at 55C.

Nucleic acid extraction (in the chemical safety cabinet)

7m

- 1 For **HA filter** extraction, add 6.5 μ L of beta-Mercaptoethanol to each of the 2-mL tubes. Go to **step 2**.

For **BCoV** extraction, add 5 μ L of beta-Mercaptoethanol to each of the 2-mL tubes. Warm the tube at 55C for 5 min and vortex for 15 seconds. Let sit for 10 min. Skip step 2.



Add the beta-Mercaptoethanol while samples are still frozen. Let thaw at room temperature.

2

For **HA filter** extraction, place the 2-mL tubes in the bead beater.

Mini-Beadbeater-16
high-energy cell disrupter

BioSpec 607 
1 speed

2.1 Bead beat for 00:02:30

2m







Start the bead beating when the beads start to be loose in the tubes.

2.2 Cooldown the samples on ice for ⌚00:05:00 .

5m

2.3 Repeat Steps 9.1 and 9.2 once ↺ .

- 3 Centrifuge at maximum speed for 1 min at room temperature. ⚙️150000 rpm, Room temperature , 00:01:00
- 4 For **HA filter** extraction, transfer 450 µL of supernatant to a Collection Tube (provided in the RNeasy PowerMicrobiome kit).
For **BCoV** extraction, transfer all supernatant.
- 5 Add 150 µL of Solution IRS and vortex briefly to mix. Incubate at 2–8°C for 5 min.
- 6 Centrifuge at maximum speed for 1 min ⚙️150000 rpm, Room temperature , 00:01:00 . Avoiding the pellet and transfer the supernatant to a new Collection Tube.
- 7 Add 650 µL each of Solution PM3 and Solution PM4. Vortex briefly to mix.
- 8 Load 650 µL of the mixture into an MB Spin Column.
- 9 Centrifuge at max speed for 1 min ⚙️150000 rpm, Room temperature , 00:01:00 . Discard the flow-through and repeat until all the mixture has been loaded onto the MB Spin Column.
- 10 Shake to mix Solution PM5 and add 650 µL to the MB Spin Column.
- 11 Centrifuge at max speed for 1 min ⚙️150000 rpm, Room temperature , 00:01:00 .
- 12 Discard flow-through. Add 600 µL of Solution PM4.

- 13 Centrifuge at max speed for 1 min  **150000 rpm, Room temperature , 00:01:00** .
- 14 Discard flow-through and centrifuge filter at max speed for an additional 2 min
 **150000 rpm, Room temperature , 00:02:00** .
- 15 Place the MB Spin Column in a clean 2-ml Collection Tube (provided in the RNeasy PowerMicrobiome kit).
- 16 Add 60 µL of RNase-Free Water (warmed to 55°C) to the center of the MB Spin Column membrane. Incubate at room temperature for at least 1 min  **00:01:00** .
- 17 Centrifuge at max speed for 1 min  **150000 rpm, Room temperature , 00:01:00** . Discard the MB Spin Column.



The DNA/RNA is now ready for downstream applications. RNA extract may be stored in RNase-free water at -80°C for 1 year.