



Mar 20,  
2020

# Generation of SARS-COV-2 RNA transcript standards for qRT-PCR detection assays

Chantal Vogels<sup>1</sup>, Joseph Fauver<sup>1</sup>, Isabel Ott<sup>2</sup>, Nathan Grubaugh<sup>1</sup>

<sup>1</sup>Department of Epidemiology of Microbial Diseases, Yale School of Public Health, <sup>2</sup>Department of Ecology and Evolutionary Biology, Yale University

1

Works for me

dx.doi.org/10.17504/protocols.io.bdv6i69e

Coronavirus Method Development Community



Nathan Grubaugh ⚡

## ABSTRACT

The protocol describes how to generate high-quality single-stranded RNA transcript standards (starting from a virus RNA stock) targeting the nsp10, RdRp, nsp14, envelope (E), and nucleocapsid (N) coding regions for use with [China CDC](#), [Hong Kong University \(HKU\)](#), [Corman et al. \(Berlin\)](#), and [US CDC SARS-CoV-2 primer and probe sets for qRT-PCR](#). (Sequences for transcripts generated, along with their corresponding assays, are provided under 'Guidelines'.)

## GUIDELINES

### RNA transcript sequences

Transcript name	Length (nt)	Genome position	Assays	Sequence
SARS-COV-2 nsp10 RNA	704	13,122 - 13,825	China CDC-ORF1	GUGGGGGACA ACCAAUCACUA AUUGUGUUAAG AUGUUGUGUAC ACACACUGGUA CUGGUCAGGCA AUAACAGUAC ACCGGAAGCCA AUAUGGAUCAA GAAUCCUUUGG UGGUGCAUCGU GUUGUCUGUAC UGCCGUUGCCA CAUAGAUAUC CAAAUCCUAAA GGAUUUUGUGA CUUAAAAGGUA AGUAUGUACAA AUACCUACAAC UUGUGCUAUG ACCUGUGGGU UUUACACUUA AAACACAGUCU GUACCGUCUGC GGUAUGUGGAA AGGUUAUGGCU GUAGUUGUGAU CAACUCCGCGA

				ACCCAUGCUUC AGUCAGCUGAU GCACAAUCGUU UUUAAACGGGU UUGCGGUGUAA GUGCAGCCCGU CUUACACCGUG CGGCACAGGCA CUAGUACUGAU GUCGUUACAG GGCUUUUGACA UCUACAAUGAU AAAGUAGCUGG UUUUGCJAAAU UCCUAAAAACU AAUUGUUGUCG CUUCCAAGAAA AGGACGAAGAU GACAAUUUAU UGAUUCUACU UUGUAGUUAAG AGACACACUUU CUCUACUACC AACAUGAAGAA ACAAUUUAUAA UUUACUUAAGG AUUGUCCAGCU GUUGCUAAACA UGACUUCUUA AGUUUAGAAUA GACGGUGACAU GGUACCACUA UAUCACGUCAA CGUCUACUAA AUACACAAUGG CAGACCUCGUC U
<b>SARS-COV-2 RdRp RNA</b>	883	15,094 - 15,976	Berlin-RdRp	AAUAGAGCUCG CACCGUAGCUG GUGUCUCUAC UGUAGUACUUA GACCAAUAGAC AGUUUCAUCAA AAAUUAUUGAA AUCAAUAGCCG CCACUAGAGGA GCUACUGUAGU AAUUGGAACAA GCAAAUUCUUA GGUGGUUGGC ACAACAUUUUA AAAACUGUUUA UAGUGAUGUAG AAAACCCUCAC CUUAUGGGUUG GGAUUAUCCUA

AAUGUGAUAGA  
 GCCAUGCCUAA  
 CAUGCUUAGAA  
 UUAUGGCCUCA  
 CUUGUUCUUGC  
 UCGCAAACAU  
 CAACGUGUUGU  
 AGCUUGUCACA  
 CCGUUUCUUA  
 GAUUAGCUAAU  
 GAGUGUCUCA  
 AGUAUUGAGUG  
 AAAUGGUCAUG  
 UGUGGCGGUUC  
 ACUAUAUGUUA  
 AACCAGGUGGA  
 ACCUCAUCAGG  
 AGAUGCCACAA  
 CUGCUUAUGCU  
 AAUAGUGUUUU  
 UAACAUUUGUC  
 AAGCUGUCACG  
 GCCAAUGUUAA  
 UGCACUUUUAU  
 CUACUGAUGGU  
 AACAAAAUUGC  
 CGAUAAGUAUG  
 UCCGCAUUUA  
 CAACACAGACU  
 UUAUGAGUGUC  
 UCUAUAGAAU  
 AGAGAUGUUGA  
 CACAGACUUUG  
 UGAAUGAGUUU  
 UACGCAUAUUU  
 GCGUAAACAUU  
 UCUCAUGAUG  
 AUACUCUCUGA  
 CGAUGCUGUUG  
 UGUGUUUCAU  
 AGCACUUAUGC  
 AUCUCAAGGUC  
 UAGUGGCUAGC  
 AUAAAGAACUU  
 UAAGUCAGUUC  
 UUUAUUAUCAA  
 AACAAUGUUUU  
 UAUGUCUGAAG  
 CAAAUGUUGG  
 ACUGAGACUGA  
 CCUUACUAAAG  
 GACCUCAGAA  
 UUUUGCUCUCA  
 ACAUACAAUGC  
 UAGUUAAACAG

				GGUGAUGAUUA UGUGUACCUUC CUUACCCAGAU CCAUCAAGAAU CCUAGGGGCCG GCUGUUUUGUA CAUUG
SARS-COV-2 nsp14 RNA	848	18,447 - 19,294	HKU-ORF1	UAGUGCUAAAC CACCGCCUGGA GAUCAAUUUAA ACACCUCAUAC CACUUAUGUAC AAAGGACUCC UUGGAAUGUAG UGCGUAUAAAG AUUGUACAAAU GUUAAGUGACA CACUUAUUUUU CUCUCUGACAG AGUCGUUUUG UCUUAUGGGCA CAUGGCUUUGA GUUGACAUCUA UGAAGUAUUUU GUGAAAUUAGG ACCUAGGCGCA CCUGUUGUCUA UGUGAUAGACG UGCCACAUGCU UUUCCACUGCU UCAGACACUUA UGCCUGUUGGC AUCAUUCUAAU GGAUUUGAUUA CGUCUAUAAUC CGUUUAUGAUU GAUGUUAACA AUGGGGUUUUA CAGGUAACCUA CAAAGCAACCA UGAUCUGUAUU GUCAAGUCCAU GGUAAUGCACA UGUAGCUAGUU GUGAUGCAAUC AUGACUAGGUG UCUAGCUGUCC ACGAGUGCUUU GUUAAGCGUGU UGACUGGACUA UUGAAUAUCCU AUAAUUGGUGA UGAACUGAAGA UUAUUGCGGCU UGUAGAAAGGU UCAACACAUGG

				UUGUUAAGCU GCAUUAUAGC AGACAAUCC CAGUUCUAC GACAUUGGUA CCCUAAAGCUA UUAAGUGUGUA CCUCAAGCUGA UGUAGAAUGGA AGUUCUAUGAU GCACAGCCUUG UAGUGACAAAG CUUAUAAAAUA GAAGAAUUAU CUAUUCUUAUG CCACACAUUCU GACAAUUCAC AGAUGGUGUAU GCCUAUUUUGG AAUUGCAAUGU CGAUAGAUUUC CUGCUAAUCC AUUGUUUGUAG AUUUGACACUA GAGUGCUAUCU AACCUAACUU GCCUGGUUGUG AUGGUGGCAGU II
<b>SARS-COV-2 envelope (E) RNA</b>	808	26,207 - 27,116	Berlin-E	GCGUGCCUUG UAAGCACAAGC UGAUGAGUACG AACUUAUGUAC UCAUUCGUUUC GGAAGAGACAG GUACGUUAAUA GUUAAUAGCGU ACUUCUUUUC UUGCUUUCGUG GUUUUCUUGCU AGUUACACUAG CCAUCCUACU GCGCUUCGAUU GUGUGCGUACU GCUGCAAUUAU GUUACGUGAG UCUUGUAAAAC CUUCUUUUUAC GUUUACUCUCG UGUUAAAAAUC UGAAUUCUUCU AGAGUCCUGA UCUUCUGGUCU AAACGAACUAA AUUUUAUUA GUUUUUCUGUU

				UGGAACUUUAA UUUUAGCCAUG GCAGAUUCCAA CGGUACUUAUA CCGUUGAAGAG CUUAAAAAGCU CCUUGAACAAU GGAACCUAGUA AUAGGUUUCCU AUUCCUUACAU GGAUUUGUCUU CUACAAUUUGC CUAUGCCAACA GGAAUAGGUUU UUGUAUAUAU UAAGUUAUUU UCCUCUGGCUG UUAUGGCCAGU AACUUUAGCUU GUUUUGUCUU GCUGCUGUUUA CAGAAUAAUU GGAUCACCGGU GGAAUUGCUAU CGCAAUGGCUU GUCUUGUAGGC UUGAUGUGGCU CAGCUACUUA UUGCUUCUUC AGACUGUUUGC GCGUACGCGUU CCAUGUGGUCA UUCAAUCCAGA AACUAACAUUC UUCUCAACGUG CCACUCCAUGG CACUAUUCUGA CCAGACCGCUU CUAGAAAGUGA ACUCGUAAUCG GAGCUGUGAUC CUUCGUGGACA UCUUCGUUUUG CUGGACACCAU CUAGGACGCUG UGACAUCAAGG ACCUGCC
<b>SARS-COV-2 nucleocapsid (N) RNA</b>	1363	28,068 - 29,430	China CDC-N; HKU-N; US CDC-N1; US CDC-N2; US CDC-N3	GAAUUGUGCGU GGAUGAGGCU GGUUCUAAUUC ACCCAUUCAGU ACAUCGAUAUC GGUAUUUAUAC AGUUUCCUGUU UACCUUUUACA

AUUAAUUGCCA  
GGAACCUAAAU  
UGGGUAGUCUU  
GUAGUGCGUUG  
UUCGUUCUAUG  
AAGACUUUUUA  
GAGUAUCAUGA  
CGUUCGUGUUG  
UUUUAGAUUUC  
AUCUAAACGAA  
CAAACUAAAAU  
GUCUGAUAAUG  
GACCCCAAAAU  
CAGCGAAAUGC  
ACCCCGCAUUA  
CGUUUGGUGGA  
CCCUCAGAUUC  
AACUGGCAGUA  
ACCAGAAUGGA  
GAACGCAGUGG  
GGCGCGAUCAA  
ACAACGUCGG  
CCCCAAGGUUU  
ACCCAAUAUA  
CUGCGUCUUGG  
UUCACCGCUCU  
CACUCAACAUG  
GCAAGGAAGAC  
CUUAAAUUCCC  
UCGAGGACAAG  
GCGUCCAAUU  
AACACCAUAG  
CAGUCCAGAUG  
ACCAAAUUGGC  
UACUACCGAAG  
AGCUACCAGAC  
GAAUUCGUGGU  
GGUGACGGUAA  
AAUGAAAGAUC  
UCAGUCCAAGA  
UGGUUUUUCUA  
CUACCUAGGAA  
CUGGGCCAGAA  
GCUGGACUCC  
CUAUGGUGCUA  
ACAAAGACGGC  
AUCAUUUGGGU  
UGCAACUGAGG  
GAGCCUUGAAU  
ACACCAAAAGA  
UCACAUUGGCA  
CCCGCAAUCCU  
GCUAACAAUGC  
UGCAUUCGUGC  
UACAACUCCU  
CAAGGAACAAC

AUUGCCAAAAG  
GCUUCUACGCA  
GAAGGGAGCA  
GAGGCGGCAG  
UCAAGCCUCUU  
CUCGUUCCUCA  
UCACGUAGUCG  
CAACAGUUCAA  
GAAAUUCAACU  
CCAGGCAGCAG  
UAGGGGAACUU  
CUCCUGCUAGA  
AUGGCUGGCAA  
UGGCGGUGAU  
GCUGCUUUGC  
UUUGCUGCUGC  
UUGACAGAUUG  
AACCAGCUUGA  
GAGCAAAUUGU  
CUGGUAAGGC  
CAACAACAACA  
AGGCCAACUG  
UCACUAAGAAA  
UCUGCUGCUGA  
GGCUUCUAAGA  
AGCCUCGGCAA  
AAACGUACUGC  
CACUAAAGCAU  
ACAAUGUAACA  
CAAGCUUUCGG  
CAGACGUGGUC  
CAGACAAACC  
CAAGGAAUUU  
UGGGGACCAG  
GAACUAAUCAG  
ACAAGGAACUG  
AUUACAAACAU  
UGGCCGCAAAU  
UGCACAAUUUG  
CCCCAGCGCU  
UCAGCGUUCUU  
CGGAAUGUCGC  
GCAUUGGCAUG  
GAAGUCACACC  
UUCGGGAACGU  
GGUUGACCUAC  
ACAGGUGCCAU  
CAAAUUGGAUG  
ACAAAGAUCCA  
AAUUUCAAGA  
UCAAGUCAUUU  
UGCUGAAUAAG  
CAUAUUGACGC  
AUACAAAACAU  
UCCCACCAACA  
GAGCCUAAAAA



				GGACAAAAGA AGAAGGCUGAU GAAACUCAAGC CUUACCGCAGA GACA
--	--	--	--	---

RNA transcripts generated by this protocol, along with their length (in nucleotides), their position on the SARS-CoV-2 reference genome, and the assays they serve as positive controls for (since all assays are multiplexed, the target is specified).

#### MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
Q5 High-Fidelity 2X Master Mix - 100 rxns	M0492S	New England Biolabs
Agilent High Sensitivity DNA Kit	5067-4626	Agilent Technologies
Bioanalyzer RNA 6000 Nano Kit	5067	Agilent Technologies
SuperScript™ IV VIL0™ Master Mix	11756500	Thermo Fisher
MEGAscript™ T7 Transcription Kit	AM1334	Thermo Fisher
Qubit™ dsDNA HS Assay Kit	Q32851	Thermo Fisher
Qubit™ RNA HS Assay Kit	Q32852	Thermo Fisher
Mag-Bind® TotalPure NGS	M1378-01	Omega Biotek

#### STEPS MATERIALS

NAME ▾	CATALOG # ▾	VENDOR ▾
SuperScript™ IV VIL0™ Master Mix	11756500	Thermo Fisher
Q5 High-Fidelity 2X Master Mix - 100 rxns	M0492S	New England Biolabs
Mag-Bind® TotalPure NGS	M1378-01	Omega Biotek
Qubit™ dsDNA HS Assay Kit	Q32851	Thermo Fisher
Agilent High Sensitivity DNA Kit	5067-4626	Agilent Technologies
MEGAscript™ T7 Transcription Kit	AM1334	Thermo Fisher
Qubit™ RNA HS Assay Kit	Q32852	Thermo Fisher
Bioanalyzer RNA 6000 Nano Kit	5067	Agilent Technologies

#### SAFETY WARNINGS

Please be aware that full length single-stranded positive sense SARS-CoV-2 RNA is considered infectious and must be handled in BSL-2 or above conditions. Please consult with your institutional biosafety team before handling SARS-CoV-2 RNA.

#### BEFORE STARTING

**Basic outline of steps:** cDNA synthesis > PCR > Purification > Quality + Quantity check > T7 RNA transcription > Purification > Quality + Quantity check

## Preparation of cDNA

- 1 Isolate viral RNA using Omega Viral DNA/RNA kit, Trizol, or equivalent.
- 2 Many different cDNA synthesis kits can be used, but choose something that is relatively high-fidelity. The current protocol uses SuperScript IV VILO Master Mix because the enzyme has low error rates and the protocol is fast and easy.



### SuperScript™ IV VILO™ Master Mix

by Thermo Fisher

Catalog #: 11756500


Component	Volume in 20 µL reaction
SSIV VILO Master Mix	4 µL
Nuclease-free water	11 µL
Virus RNA	5 µL

- 3 Run the following cycles on a thermocycler:

Temperature	Time
25°C	10 minutes
50°C	10 minutes
85°C	5 minutes
4°C	∞

- 4 Store samples at **4 °C (same day)** or **-20 °C (up to a week)** until ready for PCR.

- 5 Prepare the following primers at a concentration of 10  $\mu$ M, depending on which assay is needed. The nucleocapsid (N) transcript is for the US CDC-N1, N2, and N3 (N3 is not currently in use), HKU-N, and China CDC-N assays. The envelope (E) transcript is for use with the Berlin-E assay, nsp14 transcript is for use with the HKU-ORF1 assay, the RdRp transcript is for use with the Berlin RdRp assay (low sensitivity), and the nsp10 transcript is for use with the China CDC-ORF1 assay. [See here for assay primer sequences and references.](#)



Target	Primer	Sequence	Position	Amplicon size
N	N-Std-T7-Fwd	TAATACGACTCACTAT AGGGGAATTGTGCGT GGATGAGGC	28068	1363
	N-Std-Rev	TGTCTCTGCGGTAAG GCTTG	29411	
E	E-Std-T7-Fwd	TAATACGACTCACTAT AGGGGCGTGCCTTTG TAAGCACAA	26207	808
	E-Std-Rev	GGCAGGTCCTTGATG TCACA	27016	
nsp14	nsp14-Std-T7-Fwd	TAATACGACTCACTAT AGGGTAGTGCTAAAC CACCGCCTG	18447	848
	nsp14-Std-Rev	AACTGCCACCATCACAA ACCA	19275	
RdRp	RdRp-Std-T7-Fwd	TAATACGACTCACTAT AGGGAATAGAGCTCG CACCGTAGC	15094	883
	RdRp-Std-Rev	CATCTACAAAACAGCC GGCC	15957	
nsp10	nsp10-Std-T7-Fwd	TAATACGACTCACTAT AGGGGTGGGGGACAA CCAATCACT	13122	704
	nsp10-Std-Rev	AGACGAGGTCTGCCA TTGTG	13806	

Primers for PCR amplification. Note that forward primers include a T7 promoter sequence.

- 6 Prepare the following PCR reaction for each sample:



**Q5 High-Fidelity 2X Master Mix - 100**

rxns

by New England Biolabs

Catalog #: M0492S

Component	Volume in 25 µL reaction
Q5 2x Master Mix	12.5 µL
10 µM Primer F-T7	1.25 µL
10 µM Primer R	1.25 µL
Nuclease-free water	8 µL
cDNA	2 µL

- 7 Run the following cycles on a thermocycler:

Temperature	Time
98°C	30 seconds
98°C	10 seconds
55°C	30 seconds
72°C	90 seconds
<b>Repeat steps 2-4 for a total of 35 cycles</b>	
72°C	2 minutes
4°C	∞

Post-PCR cleanup

- 8 Allow Mag-Bind TotalPure NGS beads to equilibrate to room temperature, vortex until homogenous.



**Mag-Bind® TotalPure NGS**

by Omega Biotek

Catalog #: M1378-01

- 9 Bring PCR product volume up to 25 µL with water (if not at volume already).

- 10 Add 20  $\mu$ L of beads to 25  $\mu$ L of PCR product, mix well, and incubate at room temperature for 10 minutes.



If your lab has a KingFisher extractor to use for post-PCR cleanup steps, our automated protocols are available here: [https://github.com/grubauglab/Kingfisher\\_protocols](https://github.com/grubauglab/Kingfisher_protocols) (use 'purification.bd2' for this cleanup)

- 11 Place tubes on a magnetic stand and incubate until solution appears clear.
- 12 Discard supernatant without disturbing the beads.
- 13 While tubes are on the magnet, add 200  $\mu$ L of 80% ethanol, incubate for 30 seconds, and discard the ethanol wash.
- 14 Repeat previous 80% ethanol wash and remove as much ethanol as possible.
- 15 Leave tubes on magnet and air dry for 5 minutes.
- 16 Remove tubes from magnet and add 20  $\mu$ L of nuclease-free water. Mix well by pipetting.
- 17 Place tubes on magnet stand. When solution appears clear, remove supernatant without disturbing the beads and place into new tubes.

#### Intermediate quantification

- 18 Quantify the DNA concentration using the Qubit High Sensitivity DNA kit (or equivalent) from 1  $\mu$ L of each product, following manufacturer's protocol.



#### Qubit™ dsDNA HS Assay Kit

by Thermo Fisher

Catalog #: Q32851

- 19 Check DNA fragment distributions of the samples using the BioAnalyzer DNA 1000 kit, following manufacturer's protocol.



#### Agilent High Sensitivity DNA Kit

by Agilent Technologies

Catalog #: 5067-4626

## Reverse Transcription (Megascript T7 kit)

- 20 Important: thaw reagents on ice, but keep the 10X reaction buffer at room temperature while assembling the reaction. Add the 10X reaction buffer after water and NTPs have been added to the tube.



### MEGAscript™ T7 Transcription Kit

by Thermo Fisher

Catalog #: AM1334

- 21 Mix equal volumes (2 µL/sample) of the four ribonucleotide solutions together and add 8 µL of the mixture for each sample.
- 22 Use 2 µL of pTRI-Xef in a standard reaction as a transcription control. This should yield a transcript of 1.89 kb.

Component	Volume in 20 µL reaction
Nuclease-free water	Add to 20 µL
NTP mixture (ATP + CTP + GTP + UTP)	8 µL
10X reaction buffer	2 µL
Purified PCR product	0.1-0.2 µg PCR product
Enzyme mix	2 µL

- 23 Mix by pipetting up and down and microfuge briefly to collect reaction mixture at the bottom of the tube.
- 24 Incubate at **37 °C for 4 hours**.
- 25 Add 1 µL of TURBO DNase, mix well and incubate **37 °C for 15 minutes**.

## RNA product purification

- 26 Add the following components to a 1.5 mL tube.

Component	Volume in 81 µL reaction
Nuclease-free water	30 µL
Lithium Chloride	30 µL
T7 transcription product	21 µL

- 27 Mix and chill at **-20 °C for >30 minutes**.
- 28 Centrifuge at **4 °C for 15 minutes** at maximum speed to pellet the RNA.

- 29 Carefully remove supernatant. Wash pellet with 1 mL of 70% ethanol and re-centrifuge at **4 °C for 15 minutes** at maximum speed to maximize removal of unincorporated nucleotides.
- 30 Carefully remove the ethanol and resuspend in 20 µL of nuclease-free water.

#### Final quantification

- 31 Quantify the RNA concentration using the Qubit High Sensitivity RNA kit from 1 µL of each product, following manufacturer's protocol.



##### Qubit™ RNA HS Assay Kit

by Thermo Fisher

Catalog #: Q32852

- 32 Check RNA fragment distributions of the samples using the BioAnalyzer RNA 6000 pico kit, following manufacturer's protocol.



##### Bioanalyzer RNA 6000 Nano Kit

by Agilent Technologies

Catalog #: 5067

- 33 Use the following [RNA copy number calculator](#) to calculate genome equivalents per uL and prepare standards.



This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited