





Oct 05, 2022

# DNA extraction | CTAB-chloroform | 96 wells plate

 Forked from [96-well CTAB-chloroform DNA extraction](#)[Lila Fishman](#)<sup>1</sup>, [Simon Joly](#)<sup>2,3</sup>, [Andrea Corkal](#)<sup>3</sup>, [J  r  me Burkiewicz](#)<sup>3</sup><sup>1</sup>University of Montana; <sup>2</sup>Montreal Botanical Garden; <sup>3</sup>Universit   de Montr  al1 *Works for me* Share[dx.doi.org/10.17504/protocols.io.5qpvor5jdv4o/v1](https://dx.doi.org/10.17504/protocols.io.5qpvor5jdv4o/v1)[Centre sur la biodiversit  ](#)[Joly Laboratory](#) Simon Joly

This is a 96-well version of the classic CTAB-chloroform plant DNA extraction (Doyle & Doyle 1987), developed by John Willis and Lila Fishman in ~2000 and since optimized by multiple lab groups working on *Mimulus* (monkeyflowers). It has been used successfully for downstream applications in other plant taxa, including *Asclepias*, *Larix*, *Cynoglossum*, and *Populus*. In the Joly lab, we used it on *Impatiens* so far.

The protocol requires more specialized equipment (see below) than the single-tube version, but allows for inexpensive high-throughput DNA extraction. The yields are sufficient in quantity and quality for next-generation sequencing applications (e.g., ddrads, whole-genome Illumina sequencing), as well as PCR marker genotyping.

Genogrinder or similar bead-beater  
Benchtop centrifuge w/ deep-bucket plate rotor (e.g., Qiagen Sigma 4-15K)  
8-channel pipettor + tips  
Hood w/ waterbath  
96-well bead dispenser

dry ice (optional)  
liquid N<sub>2</sub>  
3mm grinding balls or "beads"  
chloroform  
CTAB DNA Extraction Buffer  
>95% isopropanol  
70% ethanol  
96-well "plates" of Costar 8-tube strips (Corning#4412, #4418) and strip-caps (Corning#4408)

[dx.doi.org/10.17504/protocols.io.5qpvor5jdv4o/v1](https://doi.org/10.17504/protocols.io.5qpvor5jdv4o/v1)

<https://doi.org/10.1371/journal.pgen.1009095>

## PROTOCOL CITATION

Lila Fishman, Simon Joly, Andrea Corkal, Jérôme Burkiewicz 2022. DNA extraction | CTAB-chloroform | 96 wells plate . **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.5qpvor5jdv4o/v1>



## FORK NOTE

## FORK FROM

Forked from 96-well CTAB-chloroform DNA extraction, Lila Fishman

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Oct 05, 2022

## LAST MODIFIED

Oct 05, 2022

## PROTOCOL INTEGER ID

70840

## MATERIALS TEXT

## STEP MATERIALS

[CTAB DNA Extraction buffer](#) **Contributed by users**

## Sample drying

- 1 Put the samples to be freeze-dried at -80 and let them sit for at least one hour.
- 2 Turn on the freeze dryer and the blue pump attached to the side. The freeze dryer will take <sup>30m</sup> minutes to reach the proper conditions.
- 3 When the machine is ready, it should say "Load/Unload" in the top left corner, load the <sup>2d</sup> samples that have been stored at -80C. Click "Start" and the machine. Let the machine run for 48 hours.

Do not stack samples as they may not dry properly. The best formation we found is staggering the samples in the box to minimize how much they touch each other.

- 4 After 48 hours, click "Aerate" to stop the machine and remove the red cap from the tube on the left side to drain the moisture collected. Unload the samples.

## Tissue collection

- ## 5 Pre-label the tubes

Use this labeling scheme throughout. Please verify with PI or project lead that a projectname is unique, concise, and informative for long-term storage, and that numbers correspond exactly to individual phenotypic data.

- 6 Pre-load cleaned and autoclaved beads into plate of tubes (could use a bead dispenser if available).

- 7 Add 10mg of leaves (max. 30 mg) to each tube. Try to be consistent to have a similar yield for each extraction.

The dried leaves can be a bit tricky to add to the tubes. Rolling and folding the leaves to fit them into the tube mouth works well. Be sure to cap the tubes that you are not currently loading to avoid cross-contamination and clean off the sides of the tubes after loading.

Keep all tubes capped that you are not currently loading, you can use a second strip to do this to have tubes on both sides of the current tube you are loading capped.

Try to add a similar standard amount of leaf to each sample.

This protocol has been used with individual seedlings, silica-dried leaf tissue, herbarium specimens, senesced (brown) wild plants, and fungus-containing roots. You may need to

adjust the amount of tissue and the grinding conditions (e.g. add sand or additional small beads) for other tissues.

## Tissue grinding

- 8 Put the samples in the tissue lyser for 2 min at a frequency of 30 (1/s).

Make sure the samples are finely ground. If not, do another round of grinding.

- 9 Centrifuge the samples for 2 min @2000rpm to get powder off caps and down into bottom of the tubes.

It is very important that tissue powder not spray all over the place when you open the tubes in next step.

## DNA extraction

10

This protocol works best (less time spent balancing the centrifuge...) if done on 2 entire 96-well plates OR an even number of 8-tube strips spread across 2 plates. If you are just doing 1 plate (recommended first time), spread the 12 tube-strips across 2 plates before starting.

For all reagent-dispensing steps, you can re-use pipet tips if they absolutely do not touch sample tubes (if in doubt, get fresh ones).

For all sample-transfer steps, use a fresh box of (filtered) tips for each plate. If you need to use extra tips, pull from a separate box so that the tip row ALWAYS corresponds to the tube-strip number.

Turn on water bath in fume hood (60°C) and clear workspace in hood.

Make sure you have sufficient CTAB buffer and reagents/supplies for all steps before starting.

- $\beta$ -mercaptoethanol is nasty, so we keep a dedicated L200 in the hood to get it from the bottle (stored in hood) rather than contaminating one of the main lab ones.

- Re-cap tightly with FRESH caps.

- o powder left at the bottom of the tube.

- Place brick + paper towel on top of samples to keep caps from popping off in waterbath. Invert at least once if trying to maximize yield -- put caps back on very tightly before inverting...

- centrifuge plates for 2 min to pellet solids @4000 rpm.

- Use FILTERED tips.
- If the tips clog, bounce gently on the beads to dislodge. We use filtered tips since the multichannel pipettor can spew grindate upward into shaft if tips clog and then unclog.

Reserve tube-strips with beads out of the way in hood for later bead recovery.

- 17

Do not MEASURE chloroform into the trough unless you have a plan to clean the graduated cylinder. Decant a reasonable amount and then get more if you need it.

Chloroform melts standard troughs. Use a hard-plastic trough dedicated to it.

Rinse fresh pipette tips 1-2x with chloroform (back into trough) before attempting to dispense (prevents dripping).

If you have not double-dipped tips, return excess chloroform to bottle carefully. If you have, put it into chloroform waste bottle immediately and re-start 16.

Chloroform melts standard troughs. Use a hard-plastic trough dedicated to it.

Rinse fresh pipette tips 1-2x with chloroform (back into trough) before attempting to dispense (prevents dripping).

If you have not double-dipped tips, return excess chloroform to bottle carefully. If you have, put it into chloroform waste bottle immediately and re-start 16.

Try not to drip any chloroform outside of the tubes. If you do, wipe it off before step 18!

- 18

- 19

- 20

Set pipettor to 300  $\mu$ L and SLOW speed (if automated), and get a fresh box of tips for each plate.

21 Bring plate of centrifuged samples to hood. Place on right (if right-handed).  
One strip at a time, put a sample strip and the corresponding empty strip in the handy-dandy-cut-down-rack.  
Un-cap sample (can do before, during, or after moving, depending on preference).

**Carefully pipet off the aqueous layer (up to 300mL) of sample strip and transfer it to the matched empty strip.**

Move chloroform waste to left waste plate and return filled strip to new base.

Repeat until all samples are transferred.

The handy-dandy-cut-down rack is a costar plate that has had the bottom several cm removed (or cut away from 1/2 of it, so a tipbox or costar lid can be slid beneath). When the tube-strips are placed in this contraption, the aqueous-chloroform interface should be completely visible, but the tube-strips should be stable enough to stay upright without being held.

You can buy/improvise a different striptube-holder -- the key is that you can see the interface and pipet efficiently.

Tips: Brace the pipette tips against the inner wall of the strip-tubes above the interface, while tilting the strip somewhat in the opposite direction (to create a deeper pool on the pipette side). Don't be greedy.

Try to find an arrangement of plates/hood height (if adjustable)/arm position that lets you stably and uni-directionally transfer the samples in a series of relaxed-but-robotic motions, without crossing sample tips over either sample or recipient plates (to avoid cross-contamination/mix-ups)

If this step is super-difficult or physically uncomfortable, experiment with the ergonomics until it is not (it doesn't have to be).

IMPORTANT: There is generally a solid/chalky/slimy interface between the chloroform (dark, green or orange-y, oily) and aqueous (clear, very light green or yellow depending on tissue) layers. Avoid drawing up interface debris or chloroform. Check the tips as you pipet.



Because the next step is the precipitation of DNA, anything denser than water will stay with your DNA and potentially inhibit downstream applications. So check one last time that all is clear before proceeding.

- The pellets are MUCH more slippery during the ethanol pour than the isopropanol step. Keep an eye on them and pour slowly but continuously (i.e., don't tilt back 1/2 way

through and then try to pour again).

- 26 Air-dry the pellet for 1-2 hours at room temp. with a Kimwipe placed over the open tubes. If you have to leave for longer, re-cap loosely and dry later. Don't overdry.
- 27 Once pellets are dry (no ethanol smell), add 50µL TE buffer to each sample.  
Cap with retained or new caps.  
Flick tube-strips to mix well and/or leave out at room temperature for at least an hour to fully re-suspend.  
Refrigerate if diluting/using within a few weeks, otherwise store in -20 or -80. (Do not freeze before pellet has fully resuspended!)

If you know that your yields will be unusually high or low, re-suspend pellet in more or less water (e.g. 20μL for dead field samples of *M. guttatus* or seedlings, or 100μL for silica dried buds).

The fluorometer that we use to roughly quantify DNA is optimized for concentrations 10-100 ng/μL and 50μL generally gets values on the high end of that range for good fresh bud tissue.

If you know that your yields will be unusually high or low, re-suspend pellet in more or less water (e.g. 20μL for dead field samples of *M. guttatus* or seedlings, or 100μL for silica dried buds).

The fluorometer that we use to roughly quantify DNA is optimized for concentrations 10-100 ng/μL and 50μL generally gets values on the high end of that range for good fresh bud tissue.