




MAR 21, 2024

🌐 HuBMAP | Digestion and scRNA Analysis of Skin

Arivarasan Karunamurthy¹, HJ Park¹

¹University of Pittsburgh Medical Center

GE Research



Liz McDonough
GE HealthCare

ABSTRACT

This method details digestion and scRNA analysis of the HuBMAP tissue skin specimens.

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.bp2l62krzgqe/v1

Protocol Citation: Arivarasan Karunamurthy, HJ Park 2024. HuBMAP | Digestion and scRNA Analysis of Skin. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bp2l62krzgqe/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Mar 21, 2024

Last Modified: Mar 21, 2024

PROTOCOL integer ID: 97072

Keywords: digestion, scRNA

- 1** The fresh skin sample is digested enzymatically using [Miltenyi Biotec Whole Skin Dissociation Kit, human](#) for 2 hours and further dispersed using the Miltenyi gentleMACS Octo Dissociator.
- 2** The resulting cell suspensions are filtered through 70 µm cell strainers twice and resuspended in phosphate-buffered saline (PBS) containing 0.04% BSA.
- 3** Resulting cell suspensions are then loaded into 10× Genomics Chromium instrument (Pleasanton, CA) for library preparation. V1 and V2 single-cell chemistries were used per manufacturer's protocol.
- 4** The resultant libraries are sequenced (~200 million reads/sample), using the Illumina NextSeq-500 platform.
- 5** The output sequencing reads were then examined by quality metrics, transcripts mapped to reference human genome (GRCh38) and assigned to individual cells according to cell barcodes, using Cell Ranger (10× Genomics).