



Sep 28, 2022

Ligation

Lizhen Zhu¹¹Shenzhen Polytechnic1 *Works for me*dx.doi.org/10.17504/protocols.io.bm34k8qw

2020 IGEM Shenzhen Polytechnic China

Lizhen Zhu
Shenzhen Polytechnic

ABSTRACT

Method of plasmid construction connection (T4 DNA ligase)

DOI

dx.doi.org/10.17504/protocols.io.bm34k8qw

PROTOCOL CITATION

Lizhen Zhu 2022. Ligation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bm34k8qw>

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

lizhen,ruixin,start use this ways

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 07, 2020

LAST MODIFIED

Sep 28, 2022

PROTOCOL INTEGER ID

42844

GUIDELINES

teacher ZHANG

MATERIALS TEXT

MATERIALS






☒ T4 DNA Ligase,

500u Promega Catalog #M1804

☒ T4 DNA Ligase Buffer Thermo

Fisher Catalog #46300018

system

- 1  1.5 μ L VECTOR DNA
 4.5 μ L Insert DNA
 0.5 μ L T4 DNA ligase
 2 μ L 5X T4 DNA Ligase Buffer
 12.5 μ L ddH2O

step1

- 2 Set up a system as required, all operations on ice Mix ligation system with 100 μ L of competent and ice bath for at least 30min (thorough contact of DNA and cell, enabling DNA to get inside cell walls)

step2

- 3 42°C heat shock for 90s (competent cell swell under high temperature, cell membrane and cell wall fully contact, enabling DNA to flow inside the cell)
Ice bath for 5min (Cell shrinks and DNA is dragged through the cell membrane)

step3

- 4 Pipette 200 μ L of full nutrient culture (in clean bench), place in a shaker at 37°C for 1h
- 5 Spread the entire sample on the plate, observe after 16h~20h

6