

Oct 13, 2021

Recombinant protein expression and purification of Taq DNA polymerase

Maira Rivera¹, Javiera Reyes¹, Javiera A Avilés², Amparo Núñez¹,
Fernan Federici², Cesar A Ramirez-Sarmiento¹

¹Institute for Biological and Medical Engineering, Pontificia Universidad Católica de Chile;

²Millennium Institute for Integrative Biology (iBio), Santiago, Chile

1



dx.doi.org/10.17504/protocols.io.bya3psgn

Reclone.org (The Reagent Collaboration Network)

Tech. support email: protocols@recode.org

[Click here to message tech. support](#)



Cesar A Ramirez-Sarmiento

Pontificia Universidad Catolica de Chile

This is a slightly modified and simplified version of a protocol by Thomas G.W. Graham et al, which is available at <https://gitlab.com/tjian-darzacq-lab/bearmix> and has been described in depth in the article [10.1371/journal.pone.0246647](https://doi.org/10.1371/journal.pone.0246647), for the recombinant expression of a E602D mutant of Taq DNA polymerase in pET-28a that is available in Addgene (Addgene plasmid # 166944 ; <http://n2t.net/addgene:166944> ; RRID:Addgene_166944).

The main goal of this protocol is to eliminate the use of large volumes for dialysis and potential issues with the protein crashing out of the solution due to the use of concentrators for buffer exchange of this enzyme into storage conditions.

DOI

dx.doi.org/10.17504/protocols.io.bya3psgn

Maira Rivera, Javiera Reyes , Javiera A Avilés , Amparo Núñez, Fernan Federici, Cesar A Ramirez-Sarmiento 2021. Recombinant protein expression and purification of Taq DNA polymerase. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bya3psgn>

ANID Millennium Science Initiative Program

Grant ID: ICN17_022

ANID CONCYTEC

Grant ID: covbio0012

COVID-19, SARS-CoV-2, PCR

_____ protocol ,

Sep 16, 2021

Oct 13, 2021

53307

MATERIALS

 **PMSF Sigma**

Aldrich Catalog #P7626

 **Imidazole Sigma Catalog #I5513**

 **NaCl Sigma**

Aldrich Catalog #53014

 **HisTrap FF Crude Column Ge**

Healthcare Catalog #17528601

 **Lysozyme Thermo Fisher**

Scientific Catalog #89833

 **Glycerol Merck**

Millipore Catalog #104092

 **DTT Millipore**

Sigma Catalog #DTT-RO

 **Triton X-100 Sigma**

Aldrich Catalog #X100-100ML

 **Trizma® base Merck Millipore**

Sigma Catalog #93362

 **EDTA Sigma**

Aldrich Catalog #ED2SS

 **Nonidet P40 Substitute Sigma**

Aldrich Catalog # 74385

 **2-mercaptoethanol Sigma**

Aldrich Catalog #M6250

Buffer A, pH 8.0

[M]**50 Milimolar (mM)** Tris-HCl pH 8.0
[M]**500 Milimolar (mM)** NaCl
[M]**0.1 % volume** Nonidet P-40
[M]**0.1 % volume** Triton X-100
[M]**10 Milimolar (mM)** Imidazole, pH 8.0
[M]**5 Milimolar (mM)** 2-mercaptoethanol (BME)
[M]**5 % volume** Glycerol

Buffer B, pH 8.0

[M]**50 Milimolar (mM)** Tris-HCl pH 8.0
[M]**500 Milimolar (mM)** NaCl
[M]**0.05 % volume** Nonidet P-40
[M]**10 Milimolar (mM)** Imidazole, pH 8.0
[M]**5 Milimolar (mM)** BME
[M]**5 % volume** Glycerol

Buffer C, pH 8.0

[M]**50 Milimolar (mM)** Tris-HCl pH 8.0
[M]**100 Milimolar (mM)** NaCl
[M]**0.05 % volume** Nonidet P-40
[M]**10 Milimolar (mM)** Imidazole, pH 8.0
[M]**5 Milimolar (mM)** BME
[M]**5 % volume** Glycerol

Buffer D, pH 8.0

[M]**50 Milimolar (mM)** Tris-HCl pH 8.0
[M]**100 Milimolar (mM)** NaCl
[M]**0.05 % volume** Nonidet P-40
[M]**300 Milimolar (mM)** Imidazole, pH 8.0
[M]**5 Milimolar (mM)** BME
[M]**5 % volume** Glycerol

Buffer HA, pH 8.0

[M]**50 Milimolar (mM)** Tris-HCl pH 8.0
[M]**100 Milimolar (mM)** NaCl
[M]**0.05 % volume** Nonidet P-40
[M]**5 Milimolar (mM)** BME

[M] **5 % volume** Glycerol

Buffer HB, pH 8.0

[M] **50 Milimolar (mM)** Tris-HCl pH 8.0

[M] **1000 Milimolar (mM)** NaCl

[M] **0.05 % volume** Nonidet P-40

[M] **5 Milimolar (mM)** BME

[M] **5 % volume** Glycerol

Storage conditions

[M] **50 Milimolar (mM)** Tris-HCl, pH 8.0

[M] **100 Milimolar (mM)** NaCl



[M] **0.1 Milimolar (mM)** EDTA

[M] **50 % volume** Glycerol

[M] **3 Milimolar (mM)** DTT



DAY 1 – Plasmid transformation

1d

- 1 Transform  **100 ng** of plasmid containing Taq DNA polymerase into *E. coli* C41 competent cells using either heat shock or electroporation. ^{2h}
- 2 Spread transformed cells in LB Agar plates supplemented with [M] **0.05 mg/mL** Kan. Grow plate overnight at  **37 °C** . ^{12h}




DAY 2 – Preinoculum



















1d

- 3 Select a single colony from the LB agar plate to prepare a preinoculum in  **10 mL** LB media supplemented with [M] **0.05 mg/mL** Kan. Grow overnight at  **200 rpm, 37°C** . ^{1d}


DAY 3 – Protein Overexpression

1d

- 4 Use the full volume of the preinoculum to inoculate  **1 L** of LB (or TB) media supplemented with [M] **0.05 mg/mL** Kan (1% inoculation). Grow at  **160 rpm, 37°C** until reaching an optical density at 600 nm (OD₆₀₀) = 0.8. ^{4h}
- 5 Upon reaching OD₆₀₀ = 0.8, add IPTG to a final concentration of [M] **0.5 Milimolar (mM)** and incubate  **160 rpm, 37°C, 02:00:00** . ^{16h}

- 6 Centrifuge the cell culture  **4000 x g, 4°C, 00:20:00**. Then, resuspend the cell pellet in ^{30m}
 **50 mL** of **Buffer A** freshly supplemented with  **1.0 Milimolar (mM)** PMSF and
 **0.2 mg/mL** lysozyme.
- 7 Incubate the resuspended cells at  **80 rpm, Room temperature , 00:20:00** . ^{30m}
- 8 Sonicate on ice for  **00:08:00** using cycles of  **00:00:01** ON and  **00:00:06** OFF at ^{10m}
40% amplitude (Qsonica Q125, 125W).
- 9 On an ultracentrifugation tube, incubate the unclarified lysate at  **75 °C** for  **00:30:00** to ^{1h}
precipitate most of *E. coli* proteins, and then place on ice for  **00:05:00** . Centrifuge
 **20000 x g, 4°C, 00:20:00** and collect the supernatant. You might want to collect a small
sample for SDS-PAGE afterwards.
- 10 On a **1 mL HisTrap column (GE Healthcare)** pre-equilibrated with 10 column volumes ^{1h}
(c.v.) (here, 10 mL) of **Buffer A**, load the supernatant. Wash with 10-20 c.v. of **Buffer B**. Repeat
the wash step with **Buffer C**. Then, elute with 5 c.v. of **Buffer D**, collecting the eluted
fractions every  **0.5 mL** in 1.5 ml tubes.
- 11 To quickly pool the fractions containing the protein of interest, prepare a 96-well plate or 1.5 ^{5m}
mL tubes with  **40 µL** of 5X Bradford reagent and  **150 µL** of distilled water. Then, add
 **10 µL** of each protein fraction and compare against a blank reference sample
corresponding to  **10 µL** of **Buffer C**. You can determine your protein-containing fractions
either by absorbance at 595 nm on a plate reader or visually by comparing the blue coloration
of each fraction against the blank reference. Pool your fractions and collect a  **10 µL**
sample for SDS-PAGE.


- 12 This method was preferred over protein dialysis or Amicon protein concentration to avoid ^{5m}
using large buffer volumes and proteins crashing out of the solution.

Dilute the pooled fractions 6X in buffer containing 50 mM Tris-HCl pH 8.0 and 100 mM NaCl.
- 13 Next, load the diluted sample onto a **1 ml HiTrap Heparin column** previously equilibrated ^{30m}
with 10 c.v. (here, 10 mL) **Buffer HA**. Then, elute the protein using a 10 c.v. linear gradient
against **Buffer HB**, collecting the eluted fractions every  **0.5 mL** in 1.5 ml tubes.

This linear gradient can be achieved by connecting two containers, one with 5 c.v. Buffer HA and the other with 5 c.v. buffer HB, with a syphon or a tube, and withdrawing solution from the Buffer HA container to the column using a cheap peristaltic pump or by gravity.

- 14 Again, determine your protein-containing fractions using the Bradford assay. Pool your fractions and determine its protein concentration using the same method and collect a **10 µL** sample for SDS-PAGE. 5m
- 15 For storage, supplement your pooled fractions with **100 Milimolar (mM)** Tris-HCl pH 8.0, **0.2 Milimolar (mM)** EDTA and **6 Milimolar (mM)** DTT. Then, dilute the sample by adding glycerol up to **50 % volume** to reach final storage conditions: **50 mM Tris-HCl pH 8.0, ~100 mM NaCl, 0.1 mM EDTA, 3 mM DTT**. 5m
- 16 Generate **200 µL** aliquots of the enzyme and store it at **-20 °C** until required. 30m

IMAC SDS-PAGE Result



SDS-PAGE of all purification steps of Taq DNA polymerase. Pooled eluted fractions from IMAC and Heparin correspond to P1, HP1 and HP2.