



Oct 29, 2021



# Human\_Tissue\_Nuclei\_Isolation\_Protocol\_2021\_10\_18

Jose.Bras<sup>1</sup>, lee.marshall<sup>1</sup>, Kimberly Paquette<sup>1</sup><sup>1</sup>Van Andel Institute, Grand Rapids, MI, USA

1

[dx.doi.org/10.17504/protocols.io.bzmgp43w](https://dx.doi.org/10.17504/protocols.io.bzmgp43w)

Jose.Bras

Homebrew protocol to isolate nuclei from human frozen brain tissue

DOI

[dx.doi.org/10.17504/protocols.io.bzmgp43w](https://dx.doi.org/10.17504/protocols.io.bzmgp43w)

Jose.Bras , lee.marshall , Kimberly Paquette 2021.  
Human\_Tissue\_Nuclei\_Isolation\_Protocol\_2021\_10\_18. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bzmgp43w>



ASAPCRN

\_\_\_\_\_ protocol ,

Oct 29, 2021

Oct 29, 2021

54664

## **Stock Solutions**

*10% Triton X-100 (100mL)*

Solution	Final Conc	Volume	Notes
Triton X-100	10%	10mL	
Ultrapure Water		Fill to: 100mL	

*1M MgCl<sub>2</sub> (100mL)*

Solution	Final Conc	Volume	Notes
<b>MgCl<sub>2</sub> (FW 203.31)</b>	1M	20.331g	203.31g/1000mL = 1M 5.08275g/100mL = 1M
<b>Ultrapure Water</b>		Fill to: 100mL	

### ***Working Solutions***

\*Use PBS without Calcium Chloride and Magnesium Chloride

*PBSTA (3mL) without TritonX-100 – store at 4°C for up to 1 month*

Solution	Final Conc	1 RXN	5.2 RXN
<b>10x PBS</b>	1x	300uL	1560uL
<b>1M MgCl<sub>2</sub></b>	3mM	9uL	46.8uL
<b>D-Sucrose (342.29)</b>	0.3M	0.3081g	1.602g
<b>** RNasein (add on day)</b>	0.4 U/uL	30uL	156uL
<b>Ultrapure Water</b>		Fill to: 3mL	Fill to: 15.6mL

- Add most of the Ultrapure water, then PBS, then MgCl<sub>2</sub>, then Sucrose
- Vortex to dissolve
- Then add Ultrapure water to final volume
- Add RNasein before using

*1.4M PBS Cushion (8mL) – make fresh weekly*

Solution	Final Conc	1 RXN	5.2 RXN
<b>10x PBS</b>	1x	0.8mL	4.16mL
<b>1M MgCl<sub>2</sub></b>	3mM	24uL	124.8uL
<b>D-Sucrose</b>	1.4M	3.8336g	19.935g
<b>10% Triton X-100</b>	0.1%	80uL	416uL
<b>Ultrapure Water</b>		Fill to: 8mL	Fill to: 41.6mL

- Add half of the Ultrapure water, then PBS, then MgCl<sub>2</sub>, then Sucrose
- Vortex to dissolve
- Then add 10% Triton,
- Then add Ultrapure water to final volume

*Nuclei Wash and Resuspension Buffer – make fresh daily*

Solution	Final Conc	1 RXN	5.2 RXN
<b>10x PBS</b>	1x	400uL	2080uL
<b>50mg/uL Ultrapure BSA</b>	1%	800uL (40mg)	4160uL
<b>** RNasein (add on day)</b>	0.2 U/uL	20uL	104uL
<b>Ultrapure Water</b>		Fill to: 4mL	Fill to: 20.8mL

- Add most of the Ultrapure water, then PBS, then BSA
- Tilt to dissolve, do not vortex

- Add RNasein before using

- 1 Homogenize
  - a. Transfer tissue to facs tube
  - b. Add 2mL of PBSTA to the facs tube containing the tissue
  - c. Homogenize on ice with hand-held mixer at lowest setting, 5 sec on, 5 sec off, repeat for a total of 4 times, (do not create bubbles, check no debris after homogenizing), store on ice
  - d. Add 40uL of 10% TX100 (final conc 0.2%), pipette mix 10 times with p1000, incubate for 5min on ice. For each sample leave 2min space between
  - e. Repeat for each sample a-c for each sample
- 2 Extract nuclei and filter onto gradient
  - a. Transfer the 2mL of PBSTA/tissue to a Dounce on ice
  - b. Dounce 10 times gently to release the nuclei, do a half turn at the top and bottom movement, do not create bubbles by never lifting the Dounce out of the volume at all
  - c. Filter the nuclei through a Miracloth (Cakbiochem, #475855) onto the 8mL of 1.4M PBS Cushion
  - d. Add 1mL of PBSTA to the Dounce to resuspend any remaining nuclei from the sides of the dounce, transfer these through the same Miracloth onto the 1.4M PBS Cushion
  - e. Repeat for each sample a-e for each sample
- 3 Sucrose Gradient Centrifuge
  - a. Do not disturb the gradient
  - b. Centrifuge nuclei through the gradient @ 3000 x g for 30 min at 4°C
  - c. Take out of centrifuge immediately when finished, should see a condensed white interphase
- 4 Resuspend Nuclei
  - a. On ice, remove the top layer first using a p1000
  - b. Then remove the interphase (cell debris) using a p1000
  - c. Once the cell debris is removed, continue to use a p1000 to remove all supernatant, switch to a p200 tip to remove the final 200uL supernatant, make sure to remove all the supernatant without touching the pellet
  - d. Resuspend the pellet in 210uL of Nuclei Wash and Resuspension Buffer, pipette 10 times using a p200, then transfer to LoBind tube on ice.
  - e. Further resuspend the nuclei into singlets by pipette 5 times using a p200 against the wall of the epi
  - f. Repeat for each sample a-e for each sample
- 5 Count Nuclei
  - a. Transfer 10uL of nuclei into a fresh epi on ice
  - b. Add 10uL of Trypan Blue into the 1:2 nuclei dilution, resuspend 10 times and count nuclei numbers using a Countess
- 6 Wash Nuclei 1
  - a. After incubation, add 900mL of Nuclei Wash Buffer to the 200uL nuclei, pipette 5 times with p1000
  - b. Centrifuge nuclei @ 500 x g for 5min at 4°C
  - c. Discard most of the supernatant using a p200 tip, leaving 20uL of supernatant and pellet, being

careful not to disturb the pellet

d. Pipette 10 times using a p200 to resuspend the pellet

## 7 Wash Nuclei 2

a. Add 1000uL of Nuclei Wash Buffer to the 100uL nuclei, pipette 5 times with p1000

b. Centrifuge nuclei @ 500 x g for 5min at 4°C

c. Discard most of the supernatant using a p200 tip, leaving 20uL of supernatant and pellet, being careful not to disturb the pellet

d. Pipette 10 times using a p200 to resuspend the pellet

## 8 Wash Nuclei 3

a. Add 1000uL of Nuclei Wash Buffer to the 100uL nuclei, pipette 5 times with p1000

b. Centrifuge nuclei @ 500 x g for 5min at 4°C

c. Discard most of the supernatant using a p200 tip, leaving 20uL of supernatant and pellet, being careful not to disturb the pellet

d. Pipette 10 times using a p200 to resuspend the pellet

e. Add 100uL of Nuclei Wash Buffer, pipette 10 times using a p200 to resuspend the pellet

## 9 Count Nuclei

a. Add 5uL of nuclei to 5uL of Nuclei Wash Buffer in a fresh epi on ice

b. Add 10uL of Trypan Blue resuspend 10 times and count nuclei numbers using a Countess and take photos on the EVOS level 3

## 10 Flow cytometry

a. Add 20uL of nuclei to 180uL of Nuclei Wash Buffer in a fresh epi on ice

b. Flow QC and Count Nuclei Dilution Protocol

## 11 Dilute Nuclei to 1000 nuclei/uL

a. Add 100,000 nuclei into 100uL of Nuclei Wash Buffer Snap Freeze leftover Nuclei Protocol

## 12 Snap Freeze leftover nuclei

b. Use liquid nitrogen to snap freeze leftover nuclei and store in -80 freezer