


JUN 28, 2023

## 🌐 Differentiation of mature neurons from mouse NPCs

Laura

Sonia Cruciani<sup>1</sup>, Jacqueline Severo<sup>1</sup>, Batlle<sup>1</sup>,  
Eva Maria Novoa<sup>1</sup>

<sup>1</sup>Center for Genomic Regulation (CRG)



sonia.cruciani

### ABSTRACT

This is a step by step protocol used to differentiate mouse NPCs to mature neurons.

### MATERIALS

- Mouse NPCs: vial of mNPCs derived from mES mettl3#6Cre from TEF
- Neurobasal (Thermo scientific, cat. no. 21103-049)
- DMEM/F12 (Thermo scientific, cat. no. 11320-033 )
- N2 supplement (Gibco, cat. no. 17502-048)
- Phosphate-buffered saline (PBS; Life technologies, cat. no. 14190-169)
- EDTA 0.5M pH 8
- B27 minus vitamin A (Gibco, cat. no. 11500446)
- Non-essential amino acids (Thermo scientific, cat no.11140-035)
- 2-mercaptoethanol (Thermo scientific, cat no.31350-010)
- Glutamax (life technologies, cat no. 35050-061)
- P/S (Thermo scientific, cat no. 15140-122)
- Laminin (Sigma cat no. L2020)
- Poli-L-Ornithine solution (Sigma cat no. A004-C)
- EGF (Gibco, cat no. PMG3043)
- FGF2 (aka bFGF) (Stem cell technologies, cat no. 78003)
- DAPT (Selleck,cat no. S2215)
- BDNF (Sigma,cat no. B3795)
- FBS
- DMSO

#### ● N2B27 medium:

Component	Volume for 50ml
Neurobasal	24 mL
DMEM/F12	24 mL
P/S	0.5 mL
MEM Non-Essential Amino Acids Solution	0.5 mL

### OPEN ACCESS

#### DOI:

[dx.doi.org/10.17504/protocols.io.6qpvr4oo2gmk/v1](https://dx.doi.org/10.17504/protocols.io.6qpvr4oo2gmk/v1)

**Protocol Citation:** Sonia Cruciani, Jacqueline Severo, Laura Batlle, Eva Maria Novoa 2023. Differentiation of mature neurons from mouse NPCs. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.6qpvr4oo2gmk/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Sep 21, 2022

**Last Modified:** Jun 28, 2023

**PROTOCOL integer ID:**  
70334

GlutaMAX™- 100X	0.5 mL
B-27 without vitamin A	0.5 mL
N2	0.25 mL
β-mercaptoethanol, 1000X. Sigma M-7522	50 uL

*This medium could be stored 1 month at 4 C*

●**NPC medium:**

N2B27 + EGF 10 ng/ml + FGF2 10 ng/ml

●**differentiation medium:**

N2B27 + BDNF 40 ng/ml + DAPT 10 uM

## Coating

- 1 -Dilute Poly-L-Ornithine (PLO) 1:6 in H2O
- 2 -Fully coat the wells with diluted PLO (i.e. 500 ul per well in a 12WP)
- 3 -Leave in the incubator for minimum 2h (also overnight works)
- 4 -Aspirate the PLO and wash three times with sterile H2O
- 5 -Dilute Laminin 1:500 in sterile H2O

6 -Fully coat the wells with diluted Laminin (i.e. 500 ul per well in a 12WP)

7 -Leave in the incubator for minimum 2h (also overnight works)

## Plating - Day -1

8 -count the desired number of cells: Plate ~100.000 cells on a coverslip, ~150.000 cells in a 24WP well, ~250.000 cells in a 12WP well.

9 -seed NPCs for the final differentiation in NPC medium

## Day 0

10 -aspirate all NPC medium

11 -immediately add differentiation medium not directly on the bottom of the well but making it go through the walls to avoid disturbing the cells.

## Day 2-14

12 -Change medium (remove half volume and add half volume of new medium) avoiding the cells' contact with air every 2-3 days. If at diff day 2-3 you notice overgrowth of glial cells, add AraC 1uM to the media in order to stop proliferation.

Cells can be kept in culture for >14 days, but they are already MAP2+ and TUJ1+ at day 7.

## Tamoxifen treatment

- 13 -Add 2.5  $\mu$ M Tamoxifen to the medium, and the corresponding volume of metOH as control. Incubate for at least 6 days.