

May 27, 2024 Version 4

## Tissue processing and freezing after surgery V.4

This protocol is a draft, published without a DOI.

Bettina Ergün<sup>1</sup>

<sup>1</sup>Charité



Bettina Ergün

Charité

OPEN  ACCESS



**Protocol Citation:** Bettina Ergün 2024. Tissue processing and freezing after surgery. protocols.io <https://protocols.io/view/tissue-processing-and-freezing-after-surgery-dein3cde>Version created by [Bettina Ergün](#)

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** January 18, 2024

**Last Modified:** May 27, 2024

**Protocol Integer ID:** 100654

## Abstract

The aim of this protocol is to document the processing of fresh tissue after surgery.

The collected tissue must also be documented in the files which are linked down below

The protocol describes in detail how to process the resulting tissue from following protocols:

Gewebesammlung Frischgewebe Zystektomie

Gewebesammlung Frischgewebe Prostatektomie

Gewebesammlung Frischgewebe Nephrektomie

in the current version.



## Materials

### Transport medium

**MEM** : Prozent: 98 %

**Zellshield** : Prozent: 2 %

**MEM mit Earle's Salzen, mit 2,2 g/l NaHCO<sub>3</sub>, mit stabilem Glutamin**

**Hersteller/Lieferant:**Biosell

**Katalog #:**BS.FG0325

### **ZellShield**

**Hersteller/Lieferant:**Minerva BioLabs

**Katalog #:**13-0050

### Freezing medium

**FKS** : Prozent: 90 %

**DMSO** : Prozent: 10 %

### **FBS Superior stabil**

**Hersteller/Lieferant:**Biosell

**Katalog #:**FBS. S 0615

sterile forceps, scapels and petri dishes


1 ml Cryotubes



## Sample processing in the lab

1h

- 1 Fresh tissue should be processed on the same day if possible, but can be stored at 4°C for a maximum of overnight.

To reduce pathogens, the tissue should be incubated for  01:00:00 at

 Room temperature

2

Patient#	Organ	Visit

3

### Note

*Comments/remarks on the tissue/processing:*

- 4 The tissue is divided in the lab into pieces with an edge length of 1-2 mm.  
Mix the fragments of a tumor piece so that each aliquot is as representative as possible of the tumor area.

- 5 16 (decision from June 29, 2023; previously: 8) fragments are stored in 1 ml of freezing medium in a cryotube.  
If there are enough fragments for more than 4 cryotubes, more fragments can be combined in one tube. A maximum of 5 cryotubes of each entity are frozen.

- 6 New procedure for storage and dokumentation started in May 2024 (change of storage tubes to LVL-Tubes with Barcode)

- 6.1 scan tube numbers and fill in the list mentioned in Step 9

- 6.2 Cryo tubes should be additionally labelled as follows:  
*Patient ID, Organ, T1/T2 resp. N1/N2...*  
*No. of Pieces*  
*Date*

- 7 Freeze samples in Mr. Frosty overnight.



- 8 On the next day transfer the tubes to N2 in the right order and document the location in the N2 binder and the list mentioned in Step 9

## Documentation

- 9 Store patient data :  
**BIOBANKING und EVEs\_ab 23052024.xlsx**
- 10 *Please add pictures of:*
  - *the origin of the sample in the organ (if available)*
  - *each tissue piece that has been taken before cutting the fragments. The picture should include the tissue ID (in the picture)*
  - *the documentation sheet*

### Note