VERSION 2
SEP 08, 2023

OPEN ACCESS

Protocol Citation: Bettina Ergün 2023. Prozedere Bacterial stabs in Kultur nehmen bis Übernachtkultur ernten. [protocols.io](https://protocols.io/view/prozedere-bacterial-stabs-in-kultur-nehmen-bis-ber-czrbx52n) <https://protocols.io/view/prozedere-bacterial-stabs-in-kultur-nehmen-bis-ber-czrbx52n> Version created by Bettina Ergün

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

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
GUIDELINES

Addgene ships plasmids as transformed bacteria in stab culture format. A stab culture is a type of Luria Broth (LB) Agar media, similar to a standard LB Agar plate. Unlike an LB Agar plate, a stab culture is created by piercing the LB agar with the bacteria instead of spreading it on the surface. The bacteria in a stab culture grow from the puncture site to spread across the surface of the stab culture.

Short Term Storage of Bacterial Stab

Stab cultures should be **stored at 4°C** upon receipt. The bacteria in the stab is guaranteed to live for at least 2 weeks when stored at 4°C.

Long Term Storage - Isolating and Verifying the Plasmid and Creating a Glycerol Stock

Within **2 weeks** of receiving your new Addgene plasmid you should verify that the plasmid is correct and create a glycerol stock. Please follow Addgene's recommended protocols below.

1. [Streak Bacteria for Single Colonies](#) - use your bacterial stab to streak bacteria onto a plate, grow overnight, and isolate single colonies.

2. Select an isolated, single colony to [inoculate an overnight culture](#). Using a liquid culture will allow you to grow enough bacteria for plasmid DNA purification and for [creating glycerol stocks](#).
1. Addgene recommends selecting and inoculating 2-4 separate colonies for plasmid screening and verification.
3. [Isolate your Plasmid DNA](#) - many companies sell miniprep or maxiprep kits for easy plasmid isolation. Addgene also provides a protocol for plasmid purification without a kit.
4. After isolating plasmid DNA, Addgene recommends verifying the plasmid by diagnostic digest or sequencing before beginning any experiments.
 1. [Perform a Diagnostic Digest](#) - verify backbone and insert sizes
 2. [Sequence your Plasmid](#) - verify key regions of the plasmid using DNA sequencing and compare these to sequences Addgene obtained from the plasmid, found on the plasmid page under "Sequences"
 3. After screening 2-4 colonies, if you find that your plasmid DNA is not correct, please view the instructions on how to [report a problem with a plasmid](#).
5. Once you have verified that the plasmid DNA is correct, keep the glycerol stocks corresponding to the colony that was verified- most plasmids are stable for years when stored as bacterial glycerol stocks.

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Protocol

NAME
E.Z.N.A.® Plasmid DNA Maxi Kit Centrifugation Protocol (modified)

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unter der Bench + AB; gießen (3-4 mm hoch)
mit AB einbringen