




MAR 07, 2024

Rab7 Phosphorylation reaction

 In 1 collection

Dan Tudorica¹

¹University of California, Berkeley



Dan Tudorica
Hurley Lab, QB3, UC Berkeley

ABSTRACT

Phosphorylation of Rab7 at S72 using TBK1

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.261ged92wv47/v1

Protocol Citation: Dan Tudorica 2024. Rab7 Phosphorylation reaction. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.261ged92wv47/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: Nov 09, 2023

Last Modified: Mar 07, 2024

Phosphorylation reaction

- 1 Mix phosphorylation reaction mixture. Combine 4 uM human Rab7 with 400 nM purified human strep-tagged TBK1 in a buffer consisting of 50 mM HEPES 7.5, 150 mM NaCl, 10 mM MgCl₂, and 200 uM ATP. Incubate at room temperature for 3 hours.
- 2 Pass reaction mixture over a strep-tactin sepharose gravity column. Pass reaction mixture over resin x4 in order to remove all strep tagged TBK1 from solution.
- 3 Buffer exchange via dilution and centrifugation or dialysis overnight in order to recover purified Rab7

PhosTag gel

- 4 In order to assess degree of phosphorylation, prepare a 15 ug sample of phosphorylated Rab7 for PhosTag gel electrophoresis. Buffer exchange the protein into pure water via dilution and concentration with centrifugal concentrator.

As a positive control, prepare a sample of purified Casein, purchasable in bulk as a lyophilized solid. This protein is natively phosphorylated. Incubate with lambda phosphatase in order to produce its unphosphorylated form. Run the native casein as a positive control for the PhosTag gel, and the unphosphorylated form as a negative control.
- 5 Follow guidelines for Fujifilm Wako's "SuperSep" precast PhosTag gel. Add loading buffer to 1x, and run gel at ~120 V until the dye front reaches the bottom of the gel.
- 6 Stain gel via Coomassie staining, and image.

