



MAR 25, 2024

bHEV genotyping RT-PCR

Siu Fung Ho¹

¹University of Hong Kong



Siu Fung Ho
University of Hong Kong

ABSTRACT

bHEV genotyping RT-PCR

PROTOCOL REFERENCES

<https://bmcbioinformatics.biomedcentral.com/articles/10.1186/s12859-022-04781-0#Sec11>

MATERIALS

- Superscript IV Reverse Transcriptase
- 2X Plus High Fidelity Fast PCR Master Mix
- Genotyping Primers
- Nuclease-free water
- Template RNA
- TBE Buffer
- Agarose
- 1 Kb ladder/Loading Dye

OPEN  ACCESS



Protocol Citation: Siu Fung Ho 2024. bHEV genotyping RT-PCR. protocols.io
<https://protocols.io/view/bhev-genotyping-rt-pcr-dabx2apn>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development
We are still developing and optimizing this protocol




Created: Mar 08, 2024



Last Modified: Mar 25, 2024

PROTOCOL integer ID: 96343

Reverse Transcription





- 1 Pre-warm the 5x SSIV buffer to RT
- 2 Make up the primer annealing master mix in a 2ml Eppendorf. Per sample volumes:

Component	Volume
dNTPs 10mM	 1 μL
Reverse Primer 2mM	 1 μL
Nuclease-free water	Up to  13 μL

- 3 Divide  13 μL of master mix between reaction tubes.
- 4 Add up to  5 μg RNA template to reaction tubes.
- 5 Mix and briefly centrifuge.
- 6 Heat the RNA-primer mix at 65°C for 5 minutes.
- 7 Incubate on ice for at least 1 minute.

8 Vortex and briefly centrifuge the 5x SSIV Buffer.

9 Make up the RT master mix in a 2ml Eppendorf. Per reaction volumes:

Component	Volume
5x SSIV Buffer	 4 µL
DTT 100 µM	 1 µL
RNase Inhibitor	 1 µL
Reverse Transcriptase (200 U/µL)	 1 µL

Note

Mastermix should be made up and aliquoted into PCR tubes in a mastermix cabinet. Cabinet and tubes should be cleaned with decontamination wipes/70% ethanol and UV sterilised before and after use.

10 Mix and briefly centrifuge RT master mix

11 Add  7 µL of RT master mix to annealed RNA

Note

This should be done in the pre-PCR template addition room.






12 Incubate the combined reaction mixture at 50°C for 10 minutes

- 13 Inactivate the reaction by incubating it at 80°C for 10 minutes. During this incubation, PCR mastermix can be setup (next step).

Polymerase Chain Reaction

18m

- 14 Setup mastermix in a 1.5 mL eppendorf tube. The following amounts are for a single reaction

Component	Volume
2x Reaction Mix	 25 µL
Forward Primer 10 µM	 2 µL
Reverse Primer 10 µM	 2 µL
Nuclease-free water	 16 µL
Total	 50 µL

Note

Mastermix should be made up and aliquoted into PCR tubes in a mastermix cabinet. Cabinet and tubes should be cleaned with decontamination wipes/70% ethanol and UV sterilised before and after use.



- 15 Add  5 µL template cDNA to the PCR tubes.

















Note

This should be done in the pre-PCR template addition room.

- 16 Gently mix contents of each tube by pipetting and spin down.

- 17 Program the following PCR cycles into the thermal cycler.

Step	Temperature	Time	Cycles
Reverse transcription	 50 °C	 00:10:00	1

Inactivation	 98 °C	 00:02:00	1
Denaturation	 98 °C	 00:00:10	10
Annealing	 50 °C	 00:00:10	10
Extension	 72 °C	 00:00:30	10
Denaturation	 98 °C	 00:00:10	30
Annealing	 60 °C	 00:00:10	30
Extension	 72 °C	 00:00:30	30
Final extension	 72 °C	 00:05:00	1