



1 ▼

Oct 27, 2021

# Biofilm growth with starch treatment V.1

Bjorn Bartholdy<sup>1</sup>, a.g.henry<sup>1</sup><sup>1</sup>Leiden University

protocol .

BYOC



Bjorn Bartholdy

Protocol to grow mature, mineralised oral biofilm with starches.

Modified from protocol by Sissons et al.:

Sissons, C. H., Cutress, T. W., Hoffman, M. P., & Wakefield, J. S. J. (1991). A Multi-station Dental Plaque Microcosm (Artificial Mouth) for the Study of Plaque Growth, Metabolism, pH, and Mineralization: *Journal of Dental Research*.

<https://doi.org/10.1177/00220345910700110301>

Bjorn Bartholdy, a.g.henry 2021. Biofilm growth with starch treatment.

**protocols.io**

<https://protocols.io/view/biofilm-growth-with-starch-treatment-bu7jnzkn>



European Union's Horizon 2020 research and innovation program

Grant ID: STG-677576 ("HARVEST")

protocol ,

May 20, 2021

Oct 27, 2021

50123

The preferred substratum for inoculation is glass or hydroxyapatite. Plastic substrata can be used but are less effective, so surface treatment of the plastic is recommended (e.g. heated HCl or acetone).

Substrata should be autoclaved prior to use, both before and during the experiment.

### Solutions

Artificial saliva



Artificial saliva  
by Bjorn Bartholdy

PREVIEW

RUN

CPMU



CPMU  
by Bjorn Bartholdy

PREVIEW

RUN

20% (v/v) sterile glycerol in dH<sub>2</sub>O

5% (w/v) sucrose in dH<sub>2</sub>O

0.25% (w/v) potato starch in dH<sub>2</sub>O

0.25% (w/v) wheat starch in dH<sub>2</sub>O

0.50% (w/v) equal parts wheat (0.25%) + potato (0.25%) in dH<sub>2</sub>O

### Equipment

24 deepwell polypropylene plates (w. lid containing pegs suspended from the lid)

Shaking incubator

### Saliva donor criteria

- Must have no/limited history of dental caries
- Must not have used antibiotics in the past 6 months
- Abstain from oral hygiene 24 hours prior to donation
- Refrain from eating and drinking (except water) 2 hours before donation

For experiments involving starches, donors should avoid eating starch-containing foods on the day of donation. To make this more bearable, saliva donation should take place in the morning before breakfast.

#### Saliva collection




- 1 Saliva donors rinse their mouth with water for 30 seconds.

- 2 Stimulate saliva production by chewing tasteless gum or parafilm.
- 3 Collect the saliva by spitting into 50 ml plastic centrifuge tubes.

Make sure donors wear gloves to avoid contamination with non-oral bacteria.

- 4 Make a 2-fold dilution of saliva in **sterile 20 % (v/v) glycerol** and vortex the solution.

#### Day 0: Inoculation and feeding

- 5 Before inoculation, vortex the saliva solution again.
- 6 Pipette the saliva solution into the wells, so approx. 1-2 cm of the substratum is submerged.
- 7 Place the plate in the incubator at 36°C for 4 hours for static inoculation  **04:00:00** <sup>4h</sup>
- 8 After inoculation, transfer the samples to a new plate containing the artificial saliva, and place <sup>4h</sup> in a shaking incubator at  **30 rpm, 36°C** for  **04:00:00**





Artificial saliva  
by Bjorn Bartholdy

PREVIEW

RUN



- 8.1 Add  **300 mL** distilled (or deionized) dH<sub>2</sub>O to a  **1000 mL** beaker, with stirring and heat

8 60 °C .

8.2 Add:

- ☐ ☐ 2.5 g ☐ Mucin from porcine stomach (Type III) Sigma Aldrich Catalog #M1778
- ☐ ☐ 5 g ☐ Trypticase™ Peptone Thermo Fisher Catalog #211921
- ☐ ☐ 10 g ☐ Oxoid™ Proteose Peptone Thermo Fisher Catalog #LP0085B
- ☐ ☐ 5 g ☐ Bacto Yeast Extract Becton-Dickinson

Let the reagents completely dissolve before continuing to the next step

8.3 Add:

- ☐ ☐ 2.5 g ☐ KCl Contributed by users
- ☐ ☐ 0.35 g ☐ NaCl Contributed by users
- ☐ ☐ 0.2 g ☐ CaCl<sub>2</sub> Contributed by users
- ☐ ☐ 0.74 g ☐ Sodium phosphate dibasic Sigma Aldrich Catalog #7558-79-4
- ☐ ☐ 0.54 g ☐ NaHCO<sub>3</sub> Contributed by users
- ☐ ☐ 2.5 mg ☐ Hemin Contributed by users

8.4 Add the remaining ☐ 700 mL distilled (or deionized) dH<sub>2</sub>O

8.5 Adjust to ☐ pH 7 with ☐ NaOH Contributed by users and stirring

8.6 Transfer to two 1000 ml bottles, so half of each bottle is filled.

15m

Autoclave at 8 121 °C , ☐ 1 Bar for ☐ 00:15:00 minutes



Do NOT screw bottle caps on tightly.

Loosely screw the caps on the bottles or cover the tops with foil

8.7 Once the solution has cooled, add:

- **1 mg** **Menadione Contributed by users**
- **0.3 g** **Urea Contributed by users**
- **0.17 g** **L-Arginine Contributed by users Catalog #A5006**

8.8 Store in fridge at ca. **4 °C**

Occasionally test the pH to ensure it stays around **pH 7**

9 Transfer the samples to a plate containing a 5% (m/v) **Sucrose Contributed by users** <sup>6m</sup> solution for **00:06:00** , then transfer back to the artificial saliva and leave **Overnight** .

Day 1-2: Feeding 8h 12m

10 First thing in the morning, transfer the samples to a new plate containing a 5% (m/v) **Sucrose Contributed by users** <sup>6m</sup> solution for **00:06:00** . While in the sucrose solution, add more artificial saliva to the wells on the original plate that have been partially depleted overnight.

11 After the 6 mins. return the samples to the plate with **artificial saliva**, and cover up the <sup>8h</sup> sucrose plate and leave for **08:00:00** .

12 After 8 hours, transfer the samples back to the plate with 5% (m/v) <sup>6m</sup> **Sucrose Contributed by users** solution for **00:06:00** . Transfer back to the artificial saliva and leave **Overnight** . Dispose of the sucrose.

Day 3: Inoculation and feeding 8h

13 Repeat steps from saliva collection and Day 0: Inoculation and feeding.

8h



A layer of clear plaque should be visible on the substrata

14 Prepare a new plate with artificial saliva. Transfer the samples from the inoculation plate to the **artificial saliva**.

Day 4: Feeding

8h

15 Repeat steps 10 through 12

8h

[go to step #10](#)

Prepare a new plate of artificial saliva every third day throughout the experiment. Every other morning, top up the wells with artificial saliva (so ca. 1-2 mm of the substratum is submerged).

Day 5: Inoculation

16 Repeat steps 1 through 9

[go to step #1](#)

Day 6-8: Feeding

17 Repeat steps 10 through 12

[go to step #10](#)

Day 9-14: Starch treatment

8h 6m

18 Transfer the samples to a plate with the starch treatment(s) for [00:06:00](#) at

6m

[60 rpm, 36°C](#)



If you are using multiple starch treatments within a single plate, take care to avoid cross-contamination.

19 Transfer the samples back to the artificial saliva plate for ⌚ 08:00:00 at 🌀 30 rpm, 36°C<sup>8h</sup>

20 Transfer the samples to a plate with the starch treatment(s) for ⌚ 00:06:00 at 🌀 60 rpm, 36°C<sup>6m</sup>

21 Transfer the samples back to the artificial saliva plate at 🌀 30 rpm, 36°C and leave ⌚ Overnight

Day 15-24: Mineralisation 8h 30m

22 Transfer the samples to a plate containing the calcium phosphate monofluorophosphate uréa (CPMU) solution for ⌚ 00:06:00 at 🌀 60 rpm, 36°C<sup>6m</sup>



CPMU

by Bjorn Bartholdy

PREVIEW

RUN



22.1 Add 🧴 300 mL distilled (or deionized) dH<sub>2</sub>O to a 🧴 1000 mL beaker, with stirring and heat 🔥 60 °C

22.2 Add:

- 🧴 1.55 g ☒ Calcium Chloride Contributed by users
- 🧴 1.44 g ☒ Sodium Phosphate monobasic Contributed by users
- ☒ Sodium Fluorophosphate Sigma –
- 🧴 0.72 g Aldrich Catalog #344443
- ☒ Magnesium Chloride Fisher
- 🧴 0.08 g Scientific Catalog #AC223210010
- 🧴 30 g ☒ Urea P212121

22.3 Add the remaining 🧴 700 mL and keep stirring until precipitate has completely dissolved

Store in fridge at **4 °C**

- 23 Transfer the samples back to the **artificial saliva** for **02:00:00** at **30 rpm, 36°C** 2h

Put a lid on the plate with CPMU (or cover with foil) to prevent evaporation.

Repeat step 22 and 23, four more times every two hours.

- 24 Transfer the samples to a plate with the **starch treatment(s)** for **00:06:00** at **60 rpm, 36°C** 6m

- 25 Transfer the samples back to the **artificial saliva plate** at **30 rpm, 36°C** and leave **Overnight**

#### Analysis

- 26 Samples should be dried before sampling.
- Transfer to a new plate with no liquid and leave in the incubator.
- 27 Once dried, samples are processed like archaeological samples.