

1 ▼

Dec 02, 2021

RT-qPCR Detection of SARS-CoV-2 from Wastewater Using the AB 7500 V.1

Jacquelina.Woods¹, rachel.rodriquez¹¹FDA

Jacquelina.Woods : Method Developer

1

dx.doi.org/10.17504/protocols.io.bypnpvme**GenomeTrakr**Tech. support email: genomeTrakr@fda.hhs.gov

Jessica Jones

US Food and Drug Administration

This method was developed at the FDA's Center for Food Safety and Applied Nutrition for GenomeTrakr's pandemic response project, monitoring SARS-CoV-2 variants in wastewater. Protocols developed for this project cover wastewater collection, concentration, RNA extraction, RT-qPCR detection, library prep, genome sequencing, quality control checks, and data submission to NCBI. This protocol describes triplex and duplex assays for the RT-qPCR detection of the nucleocapsid region of the SARS-CoV-2 genome. These assays, along with the murine norovirus (MNV; extraction control) and crAssphage (human indicator) RT-qPCR assay ([RT-qPCR Detection of Process Controls \(Murine norovirus and crAssphage\) from Wastewater \(protocols.io\)](https://dx.doi.org/10.17504/protocols.io.bypnpvme)), were developed for use on the AB 7500 platform using software version 2.0 or 2.3. All assays incorporate an internal amplification control (IC) to prevent the reporting of false negatives due to inhibition or failure of the RT-qPCR. These multiplexed detection assays were developed for the qualitative determination SARS-CoV-2 nucleocapsid gene extracted from wastewater. Valid sample results are contingent upon the detection of the MNV extraction control from the sample being tested.

DOI

dx.doi.org/10.17504/protocols.io.bypnpvme

Jacqueline.Woods , rachel.rodriquez 2021. RT-qPCR Detection of SARS-CoV-2 from Wastewater Using the AB 7500. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bypnpvme>



GenomeTrakr, wastewater, SARS-CoV-2, N gene, crAssphage, murine norovirus, process control, extraction control, endogenous control, RT-qPCR, AB 7500 Fast

_____ protocol ,

Sep 30, 2021

Dec 02, 2021

53710

Equipment and Supplies:

1. 7500 Fast Real-Time PCR System (ThermoFisher 4351106)
2. MicroAmp Fast Optical 96-Well Reaction Plate with Barcode, 0.1 mL (ThermoFisher 4346906), or MicroAmp Fast 8-Tube Strip, 0.1mL (ThermoFisher 4358293)
3. MicroAmp Optical Adhesive Film (ThermoFisher 4311971), or MicroAmp Optical 8-Cap Strips (ThermoFisher 4323032)
4. Platefuge (Fisher Scientific NC1823435), or Stripfuge (USA Scientific 2621-0016), or equivalent for either
5. 96-well cool block (USA Scientific 4051-0525, or equivalent)
6. Reagent cool block (USA Scientific 2312-2721, or equivalent) or ice bucket with ice
7. Adjustable calibrated micropipettes (0.2 – 1000 µl), two separate sets; one set dedicated for master mix setup and the other for template addition
8. Aerosol resistant pipet tips (0.2 – 1000 µl)
9. Personal microfuge (Labsource C90-044, or equivalent)
10. Hype-Wipe Disinfecting Towelettes (Fisher Scientific 14-412-56, or equivalent)
11. Vortex mixer (Labsource S16-200 or equivalent)

Reagents:

 [SARS-CoV-2 Genomic](#)

[RNA ATCC Catalog #1992D](#)

 [Internal Control](#)

[RNA BioGX Catalog #750-0001](#)

(Contact sales@biogx.com for ordering)

 [One-Step RT-qPCR](#)

[Kit Qiagen Catalog #210210 or 210212](#)

 [50mM MgCl₂ Bio-rad](#)

[Laboratories Catalog #17088872](#)

or

 [25 mM MgCl₂ Thermo Fisher](#)

[Scientific Catalog #AB0359](#)

 [Nuclease-free water Life](#)

[Technologies Catalog #AM9937](#)

 [Superase-In Life](#)

[Technologies Catalog #AM2696](#)

 [FAM reference dye Bio-rad](#)

[Laboratories Catalog #1708780](#)

, or equivalent, if running the

N1-N2-IC triplex assay

[ROX reference dye](#) **Fisher**

Scientific Catalog #12-223-012

, if running the N1-IC and N2-IC

duplex assays

[Tris \(1 M\) pH 8.0 RNase-free](#) **Thermo Fisher**

Scientific Catalog #AM9856

[EDTA \(0.5 M\) pH 8.0 RNase-free](#) **Thermo Fisher**

Scientific Catalog #AM9261

Primers and probe sequences in following table.

[Custom Primers and Probes](#) **IDT Technologies**

| Identification | Primers (Sequence 5' -> 3') | Location [#] |
|---------------------|--|-----------------------|
| N1F ^a | ACCCCAAATCAGCGAAAT | 28288-28306 |
| N1R | CTGGTACTGCCAGTTGAATCTG | 28357-28335 |
| N2F | TTACAAACATTGGCCGCAA | 29164-29182 |
| N2R | CGCGACATTCCGAAGAA | 29229-29213 |
| IC46F ^b | GACATCGATATGGGTGCCG | N/A |
| IC194R ^b | AATATTCGCGAGACGATGCAG | N/A |
| Identification | Probes | Location [#] |
| N1P ^c | FAM or JOE- ACCCCGCATTACGTTTGGTGGACC -IB-RQ* | 28309-28332 |
| N2P ^c | FAM or Cy5- ACAATTTGCCCCCAGCGCTTCAG -IB-RQ* | 29188-29210 |
| IACP ^b | TxRed- CCTTCCCGACCGATGGCATC -IB-RQ* | N/A |

[#]based on accession nos. JF320650, MT006214.1, and NC_045512

^aLuet al., 2020

^bDepaola, Jones, Woods et. al. 2010, U.S. Patent Application 0060166232

^cFluorophore for N1 and N2 should be FAM when running duplex assays, or JOE and Cy5, respectively, if running the triplex assay

*IB RQ- Iowa Black RQ

Always wear gloves during this procedure and never wear the same gloves when going between master mix and samples.

Always use aerosol resistant pipette tips for PCR.

Detection of SARS-CoV-2 can be conducted as a triplex assay (N1-N2-IAC), or as two

duplex assays (N1-IAC and N2-IAC).



Assembly of master mix should be done in a designated Master Mix PCR hood or BSC that has been decontaminated with 10% Bleach solution or HypeWipes followed by 70% Ethanol, or similar product and UV irradiated for 20 minutes prior to use. RNA sample template should be added in a separate designated area, physically separated from the Master Mix hood/area. Equipment should not be shared between the two areas.

Master Mix Preparation

5s

- 1 Prepare Master Mix for all sample and control reactions based on the **volumes per 25 µl reaction** in this table. Composition of mixes are listed here: [Reagent Mixes for RT-qPCR Detection of SARS-CoV-2 from Wastewater \(protocols.io\)](#) and should be prepared in advance and stored appropriately. Alternatively, Master Mixes can be prepared from individual components as described here: [Master Mix Tables for SARS-CoV-2 Assays .pdf](#).

| Reagent | Volume per 25 µl reaction |
|--|---------------------------|
| Buffer Mix | 15.55 µl |
| Primer Mix | 2 µl |
| Probe Mix | 1 µl |
| Enzyme Mix | 1.25 µl |
| FAM ^A or ROX ^B dye | 2 µl |
| Internal Control (IC) RNA* | 0.2 µl |
| | |
| RNA | 3µl |

^A Use a 1:1000 dilution (made in Primer TE) of FAM reference dye in the N1-N2-IC triplex assay.

^B Use a 1:10 dilution (made in Primer TE) of ROX reference dye in the N1-IC and N2-IC duplex assays.

*Amount varies with concentration of IC RNA. The amount of IC RNA template needs to be adjusted based on the prepared stock concentration to report a Cycle threshold (Ct) of 20-25 when no inhibition is present in the reaction (i.e., the negative control reaction).



Do not add IC or sample RNA at this step!

- 1.1 Thaw Master Mix reagents in bench top cool block (chilled at **2-8 °C**) or **On ice** in master mix preparation hood.



Keep Enzyme chilled continually; these enzymes are in glycerol and do not need to be thawed.

- 1.2 Vortex reagent tubes for **00:00:03 +/- 1 sec** at setting medium high to high (if vortex has settings).^{3s}

- 1.3 Briefly centrifuge all reagents **00:00:05 +/- 2 sec** in a personal microcentrifuge to bring liquid to the bottom of tube.^{5s}

- 1.4 Return all reagents to bench top cool block (chilled at **2-8 °C**) or **On ice**.

- 2 Proceed to hood/area or room where the template is added and thaw IC RNA and sample RNA in this hood/area.



RNA templates should be added to reaction tubes in a designated area separate from location where master mixes are prepared. A negative (water) and positive PCR control should be added to each reaction set-up.

- 2.1 Briefly centrifuge IC RNA **00:00:05 +/- 2 sec** in a personal microcentrifuge to bring liquid to the bottom of tube.^{5s}

- 2.2 Add appropriate volume of IC RNA (**0.2 µL** per reaction) to Master Mix from Step 1.4 in cold block/on ice.

2.3 Vortex briefly and centrifuge 🕒 **00:00:05 +/- 2 sec** in a personal microcentrifuge.

5s

Reaction Set-Up

5s

3 Add 🧴 **22 µL** of Master Mix to each designated reaction tube or sample wells.

4 Briefly centrifuge sample RNA 🕒 **00:00:05 +/- 2 sec** in a personal microcentrifuge to bring liquid to the bottom of tube.

5s

5 Add 🧴 **3 µL** of sample RNA template to each of three reaction tubes or wells.

6 ⚠️

Ensure each plate or run has appropriate controls (positive and negative controls) included.

Positive control prepared as described here: [Positive Control Material for RT-qPCR of SARS-CoV-2 and Process Controls \(protocols.io\)](#).

7 Seal sample plate or strip tubes. Then, briefly spin 🕒 **00:00:05 +/- 2 sec**.


5s

8 Start run on Applied Biosystems 7500 Fast instrument.


Assay parameters were optimized using the AB 7500 software versions 2.0-2.3. If other instruments or software versions are used, additional optimization may be needed.

8.1 Use the following settings for the Experiment Properties:
"7500Fast (96 wells)"

"Quantitation Standard Curve"
 "TaqMan Reagents"
 "Standard (~2 hours to complete run)"



8.2 Identify the appropriate target reporters and leave all quenchers as "NFQ-MGB".



Example in image is for the N1-N2-IAC triplex assay (N1-Joe, N2-Cy5, IC-Texas Red). If using the duplex assays, adjust appropriately (N1-FAM, N2-FAM, IC-Cy5).

8.3 Select appropriate passive reference dye (FAM for the triplex assay).

Experiment Menu << Experiment: Untitled

Setup

Experiment Properties

Plate Setup

Run Method

Reaction Setup

Materials List

Run

Analysis

Define Targets and Samples Assign Targets and Samples

Instructions:

To set up standards: Click "Define and Set Up Standards."

To set up unknowns: Select wells, assign target(s), select "U" (Unknown) as the task for each target assignment, then a

To set up negative controls: Select wells, assign target(s), then select "N" (Negative Control) as the task for each target :

Assign target(s) to the selected wells.

| Assign | Target | Task | Quantity |
|--------------------------|---------------|--|----------|
| <input type="checkbox"/> | SARS-CoV-2 N1 | <input type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> | |
| <input type="checkbox"/> | SARS-CoV-2 N2 | <input type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> | |
| <input type="checkbox"/> | IAC | <input type="checkbox"/> <input checked="" type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> | |

☒ Mixed ☒ Unknown ☒ Standard ☒ Negative Control

Define and Set Up Standards

Assign sample(s) to the selected wells.

| Assign | Sample |
|--------------------------|----------|
| <input type="checkbox"/> | Sample 1 |

Assign sample(s) of selected well(s) to biological group.

| Assign | Biological Group |
|--------|------------------|
|--------|------------------|

Select the dye to use as the passive reference.

FAM

View Plate Layout

Show in Wells

1

A

B

C

D


E

F

G

H

Wells: ☒ 0 Unknown ☒ 0 Standard ☒ 0 Negative Control

 Example in image is for the N1-N2-IAC triplex assay. If using the duplex assays, adjust appropriately (ROX).

8.4 Assign targets and samples.

Experiment Menu << Experiment: Untitled Type: S

Setup

Experiment Properties

Plate Setup

Run Method

Reaction Setup

Materials List

Run

Analysis

Define Targets and Samples Assign Targets and Samples

Instructions:
 To set up standards: Click "Define and Set Up Standards."
 To set up unknowns: Select wells, assign target(s), select "U" (Unknown) as the task for each target assignment, then assign a sample.
 To set up negative controls: Select wells, assign target(s), then select "N" (Negative Control) as the task for each target assignment.

Assign target(s) to the selected wells.

| Assign | Target | Task | Quantity |
|--------------------------|---------------|--|----------|
| <input type="checkbox"/> | SARS-CoV-2 N1 | <input type="checkbox"/> U <input type="checkbox"/> S <input type="checkbox"/> N | |
| <input type="checkbox"/> | SARS-CoV-2 N2 | <input type="checkbox"/> U <input type="checkbox"/> S <input type="checkbox"/> N | |
| <input type="checkbox"/> | IAC | <input type="checkbox"/> U <input type="checkbox"/> S <input type="checkbox"/> N | |

☐ Mixed ☐ Unknown ☐ Standard ☐ Negative Control

Define and Set Up Standards

Assign sample(s) to the selected wells.

| Assign | Sample |
|--------------------------|----------|
| <input type="checkbox"/> | Pilot_A |
| <input type="checkbox"/> | Pilot_B |
| <input type="checkbox"/> | Pilot_C |
| <input type="checkbox"/> | Positive |
| <input type="checkbox"/> | Negative |

Assign sample(s) of selected well(s) to biological group.

| Assign | Biological Group |
|--------|------------------|
|--------|------------------|

Select the dye to use as the passive reference.

FAM

View Plate Layout View Well Table

Show in Wells View Legend

| | 1 | 2 |
|---|--|---|
| A | Pilot_A IAC SARS-CoV-2 N1 SARS-CoV-2 N2 | Pilot_C IAC SARS-CoV-2 N1 SARS-CoV-2 N2 |
| B | Pilot_A IAC SARS-CoV-2 N1 SARS-CoV-2 N2 | Positive IAC SARS-CoV-2 N1 SARS-CoV-2 N2 |
| C | Pilot_A IAC SARS-CoV-2 N1 SARS-CoV-2 N2 | Negative IAC SARS-CoV-2 N1 SARS-CoV-2 N2 |
| D | Pilot_B IAC SARS-CoV-2 N1 SARS-CoV-2 N2 | |
| E | Pilot_B IAC SARS-CoV-2 N1 SARS-CoV-2 N2 | |
| F | Pilot_B IAC SARS-CoV-2 N1 SARS-CoV-2 N2 | |
| G | Pilot_C IAC SARS-CoV-2 N1 SARS-CoV-2 N2 | |
| H | Pilot_C IAC SARS-CoV-2 N1 SARS-CoV-2 N2 | |

8.5 Use the following settings for Run Method:

1h 6m 35s


☐ **25 µL** reaction volume

Holding stage 1: **50 °C** for **00:50:00**

Holding stage 2: **95 °C** for **00:15:00**

Cycling stage: 45 cycles of **95 °C** for **00:00:10** , **56 °C** for **00:00:25** ,
62 °C for **00:01:00**

Enable data collection on Step 3 of Cycling stage



These are the cycling conditions for the N1-N2-IAC triplex as well as the N1-IAC and N2-IAC duplex assays.

Data Analysis

- 9 Adjust analysis settings to appropriate thresholds. For the triplex assay, N1 and N2 should be set at **0.01** and the IC set at **0.1**. Baseline start cycle should be set at 3 and baseline end cycle should be set at 10.

Analysis Settings for Untitled

Cr Settings | Elag Settings | Advanced Settings

Review the default settings for analysis of targets in this experiment. To edit the default settings, click "Edit Default Settings." To use different settings for a target, select the target from the table, deselect "Use Default Settings," then change the settings that are displayed.

Default Cr Settings
 Default Cr settings are used to calculate the Cr for targets without custom settings. To edit the default settings, click "Edit Default Settings."
 Threshold: AUTO Baseline Start Cycle: AUTO Baseline End Cycle: AUTO [Edit Default Settings](#)

| Target | Threshold | Baseline Start | Baseline End |
|---------------|-----------|----------------|--------------|
| IAC | 0.1 | AUTO | AUTO |
| SARS-CoV-2 N1 | 0.01 | AUTO | AUTO |
| SARS-CoV-2 N2 | 0.01 | AUTO | AUTO |

Cr Settings for IAC
 Cr Settings to Use: ☐ Use Default Settings
☐ Automatic Threshold
 Threshold: 0.1
☒ Automatic Baseline
 Baseline Start Cycle: 3 End Cycle: 15

[Revert to Default Analysis Settings](#) [Apply Analysis Settings](#) [Cancel](#)




Example in image is for the N1-N2-IAC triplex assay. If using the duplex assays, adjust appropriately (all thresholds set at **0.01**).

- 10 Verify positive and negative calls for each reaction using either linear or log amplification plots.


10.1

Linear Amplification Plot




10.2

Log Amplification Plot



10.3

Linear Amplification Plot




False Positive

4/9/2019

10.4

Amplification Plot



False Positive

4/9/2019

11 RT-qPCR run is **invalid** if *any* of the following are observed:

1. Negative RT-qPCR control is positive (Ct value indicated) for any of the expected SARS-CoV-2 targets;
2. Positive RT-qPCR control is negative (undetermined) for expected SARS-CoV-2 targets; *or*

3. IC is negative (undetermined) in the negative RT-qPCR control; *or*

Run is **invalid** and the RT-qPCR assay must be repeated.

12 Sample is **negative** if:

1. Negative and positive RT-qPCR control reactions give appropriate results;
2. Sample reaction is negative (undetermined) for expected SARS-CoV-2 target/s; *and*
3. Internal amplification control (IC) is positive.



Valid SARS-CoV-2 results require the successful detection of MNV extraction control from the sample as described in *link protocol*. If MNV is negative for a sample, that sample is **invalid** and should not be reported as negative.

13 Sample is **positive** if:

1. Negative and positive RT-qPCR control reactions give appropriate results; *and*
2. Sample reaction is positive (Ct value indicated) for expected SARS-CoV-2 target/s.

If at least one of the triplicate sample reactions are positive, the sample should be reported as **positive**.