

Oct 27, 2022

High-throughput workflow for the genotypic characterization of transposon library variants

Lorea Alejaldre¹, Ana Mariya Anhel¹, Ángel Goñi-Moreno¹

¹Centro de Biotecnología y Genómica de Plantas, Universidad Politécnica de Madrid (UPM)-Instituto Nacional de Investigación y Tecnología Agraria y Alimentaria (INIA/CSIC), Madrid, Spain

Ángel Goñi-Moreno: angel.goni@upm.es

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.kqdg394jzg25/v1

biocomp.cbpp

ABSTRACT

This is a workflow for the genotypic characterization of transposon library variants. It has been developed using an open-source Opentrons OT-2 robot, BLASTN for genomic annotations and modular sub-protocols (e.g., PCR sample preparation, OT-2 volume transfer, OT-2 counter selection, etc) that can be used for other tasks, thus providing a general-purpose pipeline.

All steps follow a 96-well plate format for high-throughput analysis. The protocol is described for the characterization of transposon library variants generated with SEVA-Sib pBAMD1-x and pBLAM1-x plasmid sets that follow Standard European Vector Architecture (SEVA, <https://seva-plasmids.com>) and can be amplified with the standard PS1-PS6 primers. After the description of the protocol we present the results of an example generated at our laboratory (<https://biocomputationlab.com>) using the soil bacterium *Pseudomonas putida* KT2440 as acceptor strain.

DOI

dx.doi.org/10.17504/protocols.io.kqdg394jzg25/v1

PROTOCOL CITATION

Lorea Alejaldre, Ana Mariya Anhel, Ángel Goñi-Moreno 2022. High-throughput workflow for the genotypic characterization of transposon library variants. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kqdg394jzg25/v1>

FUNDERS ACKNOWLEDGEMENT

Comunidad de Madrid

Grant ID: Y2020/TCS-6555, 2019-T1/BIO-14053

MCIN/AEI

Grant ID: CEX2020-000999-S, PID2020-117205GA-I00

European Research Council

Grant ID: 101044360

KEYWORDS

High-throughput, Transposon library, marC9, Tn5, SEVA, OT-2, Opentrons, genotypic characterization, genotyping, bacterial genome, 96-well

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 19, 2022

LAST MODIFIED

Oct 27, 2022

OWNERSHIP HISTORY

Sep 19, 2022 | Lorea Alejaldre

Oct 27, 2022 | biocomp.cbpp

PROTOCOL INTEGER ID

70220

GUIDELINES

This workflow comprises the following sections: 1) Colony picking in selective media 2) Counter-selection and glycerol stocks pre-cultures 3) Colony selection in OT-2 liquid handler robot 4) Master 96-well plate for PCR steps 5) Control PCRs (spurious plasmid integration control and cargo insertion control) 5) Arbitrary PCRs 6) Sequencing and annotation. There is an additional section with an example on how to run the script.

We recommend the use of an OT-2 protocol specially if more than 2 libraries are to be analyzed. However, we recommend to do counter-selection in the OT-2 liquid handling robot even for one plate to avoid human errors. Note that other pipettes can be used to run the workflow in the OT-2 but these were deemed the most appropriate for the overall workflow to minimize pipette changes.

MATERIALS TEXT

Equipment:

Incubating mini-shaker

Incubating shaker

Fisherbrand 15554070 [Link](#)

OT-2

Liquid handler

Opentrons OT-2

SPECTROstar Nano

plate reader

BMG SPECTROstar Nano [Link](#)



Mastercycler® nexus - PCR Thermal
Cycler
Thermocycler
Eppendorf 12304943 [↗](#)

Centrifuge Tube Mini-Cooler
Cold Block
BRAND 10141921 [↗](#)

Electrophoresis machine

Wet-lab requirements:

Material

- Non-treated flat bottom sterile 96-well plates (i.e [↗ 96-well plates flat bottom non-treated VWR](#)
Avantor Catalog #734-2781)
- Sterile breathable membrane for 96-well plate (i.e [↗ Greiner Bio-One Sellador BREATHseal™ Fisher](#)
Scientific Catalog #11920667)
- Storage membrane for 96-well plate (i.e [↗ Thermo Scientific™ Láminas de papel de aluminio adhesivas para placas de PCR Fisher](#)
Scientific Catalog #10130853)
- PCR plates (pink ones)
- Lids for plates

Enzymes

- DNA polymerase with green buffer
[↗ Phire Green Hot Start II PCR Master Mix Fisher](#)
Scientific Catalog #15391732
[↗ Thermo Scientific™ Phire Hot Start II PCR Master Mix Fisher](#)
- DNA polymerase (**Scientific Catalog #15361732**)
- Agarose

Oligonucleotides:


- Spurious integration control: PS3, PS4, PS5, PS6 (
- Arbitrary PCR: ARB2, ME-O-Km-Ext-F, ME-O-Km-Int-F, ME-O-Sm-Ext-F, ME-O-Sm-Int-F, ME-O-Gm-Ext-F, ME-O-Gm-Ext-R (<https://doi.org/10.3389/fbioe.2014.00046>)
- Optional insert control: PSMCS

Dry-lab requirements:

- Python 3
- Command-line BLASTN

Colony picking in selective media

1d

- 1 Dispense  **100 µL** of selective media (M9-citrate for *P. putida* or Luria-Bertani plus 20 ng/µL nalidixic acid for *DH5α E. coli*) plus **transposon cassette antibiotic** in a 96-well plate
- 2 Pick individual colonies into a 96-well plate with selective media
Tip: Keep tips inside of wells to keep track
- 3 Cover with a sterile breathable membrane







Grow  **Overnight** at  **30 °C** (*P. putida*) or  **37 °C** (*E. coli*) /  **500 rpm**

Counter-selection and glycerol stocks pre-cultures

1d



Measure OD_{600nm} of overnight culture grown in selective media from  **go to step #4** plus **transposon cassette antibiotic** in a plate reader

- 6 Inoculation of **counter-selection plate** in selective media: 10m
 - Dispense  **100 µL** of selective media (M9-citrate for *P. putida* or Luria-Bertani plus 20 ng/µL nalidixic acid for *DH5α E. coli*) plus **ampicillin** (backbone antibiotic) to select against spurious integration events.
 - Transfer  **5 µL** of overnight culture from  **go to step #4** to counter-selection (ampicillin) plate
 - Cover with a sterile breathable membrane

For steps 6 and 7, if two or more 96-well plates are used as input it is advised to use the **OT-2** protocol below to minimize human error. Dispensed volume and culture volume inoculated should be that described in these steps. Note that steps 6 and 7 could be completed together in a single run depending on the number of initial plates.




OT-2 Media dispensing and culture inoculation protocol
by biocomp.cbcp


PREVIEW

RUN



- 7 Inoculation of precultures for **glycerol stock** in rich media: 10m
- Dispense **100 µL** of Luria-Bertani media plus **transposon cassette antibiotic** in a 96-well plate
 - Transfer **5 µL** of overnight culture from [go to step #4](#) to counter-selection (ampicillin) plate
 - Cover with a sterile breathable membrane

- 8  16h
- Grow counter-selection and glycerol stock pre-culture plates [Overnight](#) at **30 °C** (*P. putida*) or **37 °C** (*E. coli*) / **500 rpm**

- 9 
- Measure OD_{600nm} of overnight culture grown in selective media plus **ampicillin** in a plate reader

Colony selection in OT-2 liquid handler robot

- 10 ■ Selection of colonies to store as glycerol stocks and do further PCR reactions by running the following OT-2 protocol with its corresponding template.csv



OT-2 Counter-Selection
by biocomp.cbcp

PREVIEW


RUN




The OT-2 protocol will prepare three plates (2 glycerol stock plates and a "PCR plate") and perform the following:

- Dispense **75 µL** of PCR-grade water to "PCR plate"
- Dispense **25 µL** of 30% glycerol to two glycerol stock plates
- Transfer **25 µL** of grown pre culture in Luria-Bertani media plus **transposon cassette antibiotic** from [go to step #7](#) to "PCR plate" and glycerol stock plates

- 11 Cover glycerol stock plates with a storage membrane and store at -80°C

- 12 
- If not proceeding to the next step right away: Store "PCR plate" at 4°C for a few days or cover with an storage membrane and store at -20°C for longer term

Master 96-well plate for PCR steps

- 13 
- Transfer **50 µL** of selected colonies from one or more libraries to a 96-well plate with the following OT-2 protocol:



OT-2 Protocol to transfer volume from several plates to a single plate
by **biocomp.cbgp**

PREVIEW

RUN



Control PCRs

14



Positive control (donor plasmid) and wild-type control (*P. putida* or *E. coli*) should be added to every reaction in this section.

Spurious integration control with SEVA primers pairs PS3/PS4 and PS5/PS6



OT-2 PCR sample preparation protocol
by **biocomp.cbgp**

PREVIEW

RUN



If <48 cfu are to be analyzed both PS3/PS4 and PS5/PS6 spurious integration controls can be done in a single 96-well plate

15



Optional: Cargo integration control with primers PSMCS and either ME-O-Km-R/ME-O-Sm-R or ME-O-Gm-R (depending on **transposon cassette antibiotic**)

Arbitrary PCRs

16



Arbitrary PCR#1 using primer pairs ARB6 and ME-O-Km-Ext-F/ME-O-Sm-Ext-F or ME-O-Gm-Ext-F depending on transposon antibiotic cassette



OT-2 PCR sample preparation protocol
by **biocomp.cbgp**

PREVIEW

RUN



The OT-2 protocol will perform the following steps:

- Prepare a PCR master mix
- Dispense 19 µL of PCR mastermix

- Transfer 2 µL of pre-culture from [go to step #12](#)



If different primer pairs are added to a single 96-well plate, the OT-2 script should be run separately for each primer pair

- 17 Seal 96-well plate, place it in thermocycler and run the following PCR program:

A	B	C
98°C	5 min	
98°C	10 s	x6 cycles
30°C	30 s	
72°C	1 min 30 s	
98°C	10 s	
45°C	30 s	x30 cycles
72°C	1 min 30 s	
72 °C	5 min	
4°C	hold	

- 18 

Select 8-12 Arbitrary PCR#1 reactions from the 96-well plate and run them on a 1% agarose gel to verify amplification.

Several bands will appear and even DNA smears even when the reaction has worked perfectly.

- 19 **Arbitrary PCR#2** using primers pairs ARB2 and ME-O-Km-Int-F/ME-O-Sm-Int-F or ME-O-Gm-Int-F depending on transposon antibiotic cassette



OT-2 PCR sample preparation protocol
by [biocomp.cbgp](#)

PREVIEW

RUN



The OT-2 protocol will perform the following steps:

- Prepare a PCR master mix
- Transfer 1 µL of PCR product from Arbitrary PCR#1
- Dispense 19 µL of PCR mastermix

20 Seal 96-well plate, place it in thermocycler and run the following PCR program:

A	B	C
98°C	30 s	x30 cycles
98°C	10 s	
52°C	30 s	
72°C	1 min 30 s	
72°C	5 min	
4°C	hold	



21 

Select 8-12 Arbitrary PCR#2 reactions from the 96-well plate and run them on a 1% agarose gel to verify amplification

Several bands will appear and even DNA smears even when the reaction has worked perfectly.

Sequencing and annotation

1m

22 Prepare a PCR plate to send to sequencing by mixing  10 µL of unpurified **Arbitrary PCR#2** reaction and  10 µL of 10 µM sequencing primer (ME-O-Km-Ext-F/ME-O-Sm-Ext-F or ME-O-Gm-Ext-F depending on the transposon antibiotic cassette)

These guidelines may vary depending on the sequencing service arranged for your laboratory.

23 

1m

Annotate sequencing results by running the following protocol:



Bacterial genome annotation script using BLASTN
by biocomp.cbgp

PREVIEW

RUN



The python script uses as input:

1. DNA sequencing results in .txt or .seq
2. Reference genome file in fasta format
3. Genome annotation file in fasta format

Example

24 In this section we show a specific example on how to run the workflow with the OT-2 liquid handler using a starting point three 96-well plates with *P. putida* KT2440 colonies picked from matings using cargos with either a kanamycin (pBLAM1-2), streptomycin (pBLAM1-4) or gentamicin (pBLAM1-6) resistance gene.

25 **Counter-selection and glycerol stocks pre-cultures** using the following OT-2 protocol in [go to step #6](#) in two runs

25.1 The following template. csv will dispense M9-citrate media containing ampicillin to three plates. Subsequently, each plate will be inoculated from each of the three starting plates.

 [Variables-AntibioticPlatesCreation-OT.csv](#)

25.2 The following template. csv will dispense LB media with either gentamicin, kanamycin or streptomycin to each plates. Subsequently, each plate will be inoculated from each of the three starting plates.

 [Variables-AntibioticPlatesCreation-OT.csv](#)

26 **Colony selection in OT-2 liquid handler robot.**

Each LB plus cargo antibiotic plate will be run separately with the script in [go to step #10](#) . For each library there is an input .csv file for the absorbance at 600nm in M9-citrate with either kanamycin, gentamicin or streptomycin and another .csv file for the absorbance at 600nm in M9-citrate ampicillin. The LB plate from step 25.2 will be used to inoculate the two glycerol stock plates and PCR plate.

26.1 pBLAM1-2: .csv files and template .csv to run script

 [Variables-ColonieScreening-OT.csv](#)

 [220920_M9Cit_ant_LibpBLAM12.csv](#)

 [220921_M9Cit_amp_LibpBLAM12.csv](#)

Output from the run (map with identifiers)

 [map_220921_LibpBLAM12.csv](#)

26.2 pBLAM1-4: .csv files and template .csv to run script

 [Variables-ColonieScreening-OT.csv](#)

 [220920_M9Cit_ant_LibpBLAM14.csv](#)

 [220921_M9Cit_amp_LibpBLAM14.csv](#)

Output from the run (map with identifiers)

 [map_220921_LibpBLAM14.csv](#)

26.3 pBLAM1-6: .csv files and template .csv to run script

 [Variables-ColonieScreening-OT.csv](#)

 [220920_M9Cit_ant_LibpBLAM16.csv](#)

 [220921_M9Cit_amp_LibpBLAM16.csv](#)

Output from the run (map with identifiers)

 [map_220921_LibpBLAM16.csv](#)

27 **Master 96-well plate for PCR steps** to combine cultures from different transposon libraries in a single 96-well plate

The following template.csv is to be run with the OT-2 script in [go to step #13](#)

 [Variables-SamplesMerging-OT.csv](#)

Output from the run (map with identifiers)

 [maps_master_pcr-220923_5.csv](#)


28 **Control PCRs** to account for spurious integrations and the correct integration of the cargo

28.1 **Spurious integration control** with PS3/PS4 and PS5/PS6 SEVA oligonucleotides


The following template.csv was used for the OT-2 script in [go to step #14](#) to prepare each mastermix with either PS3/PS4 or PS5/PS6 primer pairs and transfer cultures from Master 96-well plate

 [Variables-PCRs-OT.csv](#)


Gels after PCR protocol in thermocycler




Gel of the first 7 rows with the primers ps3 and ps4



Gel of the last 5 rows with the primers ps3 and ps4



Gel of the first 7 columns with the primers ps5 and ps6




Gel of the last 5 columns with the primers ps5 and ps6

28.2 Cargo integration control with PSMCS and either ME-O-Km-R/ME-O-Sm-R or ME-O-Gm-R


The following template.csv was used for the OT-2 script in [go to step #14](#) to prepare the mastermix and transfer cultures from Master 96-well plate

We have discarded the first 3 columns because they came as negative (without integration) in the spurious PCRs.


Gels after PCR protocol in thermocycler



Gel of the columns 4 to 9 of the control of integration primers



Gel of the columns 10 and half of the 11 of the control of integration primers



Gel of the columns 11 and half of the 12 of the control of integration primers

29 Arbitrary PCRs

29.1 Arbitrary PCR#1

The following three different template.csv were used (one for each primer pair) in three independent runs with the OT-2 script in [go to step #16](#) using the same input and output plate to amplify the genomic context around each transposon cassette

Variables for source plate 1 arbitrary PCR 1: [Variables-PCRs-OT.csv](#)

Variables for source plate 2 arbitrary PCR 1: [Variables-PCRs-OT.csv](#)

Variables for source plate 3 arbitrary PCR 1: [Variables-PCRs-OT.csv](#)

29.2 Arbitrary PCR#2


The following three different template.csv were used (one for each primer pair) in three independent runs with the OT-2 script in [go to step #19](#) using the Arbitrary PCR#1 plate as input and the same output plate to amplify the genomic context around each transposon cassette

Variables for source plate 1 arbitrary PCR 2: [Variables-PCRs-OT.csv](#)

Variables for source plate 2 arbitrary PCR 2: [Variables-PCRs-OT.csv](#)

Variables for source plate 3 arbitrary PCR 2: [Variables-PCRs-OT.csv](#)

Gel results for arbitrary PCRs



Gel for several wells in the final plate of the arbitrary PCRs. .1 is result of the first PCR and .2 is the second PCR

30 Sequencing and annotation of arbitrary PCR#2

The following files were used:

[alignment_and_annotation_blastn.py](#) [Pseudomonas_putida_KT2440_110.fna](#)
[Pseudomonas_putida_KT2440_110.csv](#) [sequencing_results.zip](#)

We only sent 61 samples, that is why sequencing_results.zip has 61 files