

May 13, 2020

Cell-free lysate (E. coli) preparation with sonication

Nadanai Laohakunakorn¹¹University of Edinburgh

1

Works for me

dx.doi.org/10.17504/protocols.io.8nthvenNadanai Laohakunakorn
University of Edinburgh

ABSTRACT

Production of cell-free lysate from E. coli BL21 Star (DE3) with optional induction of T7 RNAP. Adapted from Kwon and Jewett 2015. Features:

- Variable starting culture sizes (10 mL - 1 L)
- Constant energy sonication
- S12 centrifugation
- Run-off and optional dialysis



Kwon YC, Jewett MC (2015). High-throughput preparation methods of crude extract for robust cell-free protein synthesis.. Scientific reports.



<https://doi.org/10.1038/srep08663>

Protocol successfully used at the University of Edinburgh by Nadanai Laohakunakorn, and at EPFL by Barbora Lavickova and the 2017 iGEM team.

MATERIALS TEXT

- Big centrifuge: Eppendorf 5810R with A-4-62 swing-bucket rotor (holds 50 mL tubes)
- Small centrifuge: Eppendorf 5424R with FA-45-24-11 rotor (holds 2 mL tubes)
- Vibra Cell 75186 sonicator / Qsonica Q125 + CL-18 probe
- spectrophotometer
- incubator with shaker
- autoclave
- magnetic stirrer
- LB/2xYTP/2xYTPG medium autoclaved
- E. coli strain of interest
- Tris base (Sigma T1503-100G)
- magnesium glutamate (L-glutamic acid hemimagnesium salt tetrahydrate) (Sigma T1503-100G)
- potassium glutamate (L-glutamic acid potassium salt monohydrate) (Sigma 49601-500G)
- DTT (1,4-dithiothreitol) (Sigma 10708984001)
- acetic acid
- deionized water
- autoclaved tips and petri dishes
- 500 mL Erlenmeyer flasks sterile
- culture tubes sterile
- 50 mL falcon tubes
- 2 mL eppendorf tubes
- 1 L beakers
- 10k MWCO dialysis cassettes (Slide-A-Lyzer, 3mL, Life Technologies)
- magnetic stir bar
- spectrophotometer cuvettes
- liquid nitrogen, dewar, -80 storage

1 Prepare all materials for bacterial culture

This protocol will make  **200 ml** of culture to yield around  **1 ml** of lysate with a protein concentration of around 40 - 80 mg/mL.


- 1.1 Reconstitute LB media from premix as per instructions on box.



Protocol has also been successfully carried out with 2x YTP and 2x YTPG (34°C, original). Cells grow slower on YTP media (~5-6h to OD 1.5-1.8) compared to ~4h in LB.

- 1.2 Autoclave LB media for  **00:20:00** at  **121 °C** along with tips and dishes

2 Grow overnight mini-culture

- 2.1 Add  **5 ml** of LB medium to a culture tube and inoculate with small amount of bacteria from glycerol stock.





Strains used successfully with this protocol are

- BL21 (DE3) and BL21 star (DE3) (lacking RNaseE), T7 and E. coli RNAP
- Rosetta (DE3) for oscillators (contains rare tRNA encoding plasmid, good for eukaryotic protein expression)
- Top10, Top10-GamS (GamS inhibits RecBCD degradation of linear DNA), E. coli RNAP, alpha-complementation
- M15 (E. coli RNAP, alpha-complementation)
- DH5alpha (E. coli RNAP, alpha-complementation)

- 2.2 Grow overnight culture in an incubator at  **37 °C** and 200 RPM. Make sure cap is in loose position to allow for air exchange

3 Grow lysis culture (the next morning)

- 3.1 **Measure and record** OD600 of overnight culture at **10x dilution**: pipette  **900 µl** LB and  **100 µl** overnight culture into a cuvette and take absorbance readings at 600 nm.



OD600 should be around 4 (0.4 at 10x dilution)

3.2 Add **1 ml** of overnight culture to **200 ml** of LB medium in a 500 mL Erlenmeyer flask.




3.3 Incubate the culture for **04:00:00**, **37 °C**, and 200 RPM.

4h



At the end of the step, the OD₆₀₀ should be ~1.5-2. While Kwon and Jewett grow to OD₆₀₀ = 3, we grow instead for a fixed time. They observe varying robustness of yield to final OD which is strain dependent. **In the end we should probably grow to fixed OD.**



Activity of cells harvested at different densities, grown in 1L 2xYTPG and sonicated at 1.5 mL/556 J. From Kwon and Jewett 2015, Figure 3, CC-BY license.

3.4 (If required) induce after 2 hours with 0.4 mM IPTG (**800 µl of 100 mM stock** in **200 ml** culture) to express T7 polymerase in the BL21 (DE3) strains.

4 A few minutes before incubation step finishes, prepare centrifuge and tubes

4.1 Cool down big centrifuge (Eppendorf 5810R) to **4 °C** with fast temp mode; this takes around 10 mins.

4.2 Weigh one 50 mL Falcon tube along with its cap, and record. Clearly label the tube as well as its cap; this tube will be used to spin the final pellet.

- 5 As soon as the incubation finishes, put the Erlenmeyer flask **On ice** to arrest growth.

Centrifugation and cleaning

20m

- 6 Put **50 ml** of culture into each of four 50 mL Falcon tubes (one of which has been labelled), and spin at **4000 rpm**, = **3220g** **4 °C** **00:20:00**. Tubes in swing-bucket rotor A-4-62.

- 7 Put tubes back **On ice**; keep bacterial pellet cold as much as possible.

- 8 Carefully discard supernatant using a pipette, then add **10 ml** of **Buffer A (no DTT)** to each pellet. Resuspend by carefully pipetting up and down. Finally, transfer all four parts into the single labelled tube.

10m

- 9 Balance centrifuge with a second tube containing **40 ml** water. Then spin at **4000 rpm**, = **3220g** **4 °C** **00:10:00**.

10m

- 10 Carefully discard supernatant using a pipette, then add **10 ml** of Buffer A to the pellet. Then spin at **4000 rpm**, = **3220g** **4 °C** **00:10:00**.

- 11 **go to step #10** and repeat spin once more.

- 12 Finally, carefully discard as much supernatant as possible.

- 13 Weigh the tube containing the pellet (make sure cap is on), and record:

Final weight of tube + pellet

The **wet mass** of cells is the difference between this weight and the dry weight measured earlier. This is typically around ~1 g.

- 14 Flash freeze the cells using liquid nitrogen, respecting all safety procedures (wear protective glasses as there is danger of tube explosion etc.). Store the cells at **-80 °C**.



Flash-freezing is optional (the protocol will work without) but it serves two purposes:

- Convenient pause point
- Frozen cells, once rethawed, lyse slightly more easily than fresh cells. So this step must be kept consistent

to ensure reproducible protocol.

Sonication


15 Keeping the pellet **On ice**, add **Buffer A with DTT**: **1 ml** buffer per **1 g wet mass** as determined earlier. **This ratio is critical**. The DTT must be prepared fresh as it will degrade otherwise.

16 Resuspend by vortexing, and put back **On ice**

17 Put **1 ml exactly** of resuspended cells into a new 2 mL Eppendorf tube. This volume is critical.



The larger the volume the more robust the lysis:



Yield as a function of volume and sonication energy. From Kwon and Jewett Figure 4, CC-BY license.

So if possible, scale up the production to achieve (ideally) at least 5 mL lysis volumes (= 5 g wet mass).

- 18 Place the tube in an **ice bath** - ice mixed with water - and place the sonicator tip inside the tube so that it is immersed as much as possible but **no part of the tip touches the tube surface**. This is one of the most variable parts of the protocol. A suggested technique is to drop the tip to the bottom of the tube, and then raise it slightly. Holding the tube using a cool rack (e.g. Fisher 15592801) is recommended.
- 19 Sonicate until **400J** has been achieved, using 50% amplitude and pulses of 10s : 10s (i.e. energy for 10s followed by pause for 10s). This typically takes around ⌚ **00:01:24** but the time is variable. **The total sonication energy is a critical parameter.**

Separation and Clarification


- 20 Cool down the small centrifuge (Eppendorf 5424R) to ⚡ **4 °C** using the fast temp mode: this takes around 20 minutes.
- 21 Balance and centrifuge the lysate at 🌀 **12000 x g** , = **11304 rpm** ⚡ **4 °C** ⌚ **00:10:00** 10m
- 22 Transfer the supernatant to a new tube. To prevent any transfer of bacterial debris, it is critical that you **do not take all the lysate**.
- 23 Place the lysate in an incubator at ⚡ **37 °C** , 200 RPM for ⌚ **01:30:00** , to degrade remaining DNA and RNA. This is called the run-off reaction. 1h 30m



Even though Kwon and Jewett observe deterioration in BL21 production yield with run-off, we do not see the same behaviour. If E. coli RNAP will be used, do dialysis as well as run-off.

- 24 Making sure the centrifuge is still cool, centrifuge the lysate at 🌀 **12000 x g** , = **11304 rpm** ⚡ **4 °C** ⌚ **00:10:00** 10m
- 25 Transfer the supernatant to a new tube. Again, **do not take all the lysate**. Aliquot around 📏 **3 µl** into separate tube for a Bradford assay.
- 26 **(IF REQUIRED) carry out dialysis**

- 26.1 Prepare **Buffer A + DTT** and add **900 ml** to a 1 L beaker at **4 °C** . Add a magnetic stir bar.
- 26.2 Rehydrate the required number of 10k MCWO cassettes in **Buffer A+DTT** in a beaker at **4 °C** for **00:02:00**
- 26.3 Load cassettes with **2.5 ml max (partial loading possible)** of extract and dialyze for **03:00:00** **4 °C** with stirring. 3h



- 26.4 Finally load extract into centrifuge and centrifuge the lysate at **12000 x g , = 11304 rpm** **4 °C** for **00:10:00**
- 26.5 Transfer the supernatant to a new tube. Again, **do not take all the lysate**. Proceed to step 27.

- 27 Keeping everything **On ice** , aliquot the lysate into small tubes as required (around **25 µl** per tube is recommended).
- 28 Flash-freeze the remaining tubes in liquid nitrogen and store at **-80 °C**
- 29 To perform the Bradford assay, dilute **1 µl** of the lysate in **99 µl** of Buffer A. Mix **5 µl** of this diluted lysate with **250 µl** of Bradford reagent. Incubate for **00:05:00** before measuring the absorbance using a Nanodrop. The final protein concentration is typically 40 - 80 mg/mL.

30 Preparation of 1M stock solutions of tris acetate, magnesium glutamate, and potassium glutamate

- 30.1 Weigh **12.11 g** of tris base (121.14 g/mol), add it to **100 ml** of deionized water, and adjust to pH 8.2 with acetic acid, to make 1M tris acetate stock.
- 30.2 Weigh **38.86 g** of magnesium glutamate (388.61 g/mol), and add it to **100 ml** of deionized water, to make 1M magnesium glutamate stock.
- 30.3 Weigh **20.32 g** of potassium glutamate (203.23 g/mol), and add it to **100 ml** of deionized water, to make 1M potassium glutamate stock.

31 Preparation of 1L of buffer A

- 31.1 Add **10 ml** of 1M tris acetate pH 8.2, **14 ml** of 1M magnesium glutamate, and **60 ml** of 1M potassium glutamate, and fill to **1 L** with deionized water to make buffer A (10 mM tris acetate, 14 mM magnesium glutamate, 60 mM potassium glutamate).
- 31.2 Store at **4 °C**
- 31.3 Add **2 µl** of 1M DTT to **1000 µl** Buffer A (final concentration 2 mM) for the lysis solution, or **1.8 ml** of 1M DTT to **900 ml** of Buffer A for the dialysis solution. This should be made fresh every time but if short-term storage is required, put at **-20 °C**