



SEP 20, 2023

OPEN ACCESS



Protocol Citation: Carlos Costas Insua 2023. Preparation of DEPC-treated water. **protocols.io** <https://protocols.io/view/preparation-of-depc-treated-water-cz9qx95w>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Sep 20, 2023

Last Modified: Sep 20, 2023

PROTOCOL integer ID:
88080

🌐 Preparation of DEPC-treated water

Carlos Costas Insua¹

¹Center for Molecular Biology of Heidelberg University (ZMBH), DKFZ-ZMBH-Alliance, 69120 Heidelberg, Germany.



Carlos Costas Insua

ABSTRACT

Diethyl pyrocarbonate (DEPC) is a wide-spectrum RNase inhibitor. DEPC-treated water is generally used when working with RNA

MATERIALS

Diethyl pyrocarbonate (Sigma-Aldrich #D5758)

SAFETY WARNINGS

⚠ Trace amounts of DEPC will modify purine residues in RNA by carboxymethylation. Carboxymethylated RNA is translated with very low efficiency in cell-free systems. However, its ability to form DNA:RNA or RNA:RNA hybrids is not seriously affected unless a large fraction of the purine residues have been modified. Residual DEPC must always be eliminated from solutions or vessels by autoclaving or heating to 100°C for 15 minutes.

1 Add 1 mL of fresh DEPC to 1 L of H₂O. Shake well to disperse the DEPC through the H₂O.

1d

Incubate at 37°C for at least 12 h and/or autoclave at 15 psi on liquid cycle for 20 min to inactivate the remaining DEPC.