



Upload image

Version 2

Apr 30, 2020

Ab Lysis Buffer V.2

Colin Daniel¹¹Oregon Health and Science University

1

Works for me

dx.doi.org/10.17504/protocols.io.bft7jnm

Colin Daniel

Oregon Health and Science University

MATERIALS

NAME	CATALOG #	VENDOR
100g Deoxycholate (Deoxycholic acid sodium salt)	DG090	G-Biosciences
NaCl	53014	Sigma Aldrich
Ethylenediaminetetraacetic acid	E9884	Sigma – Aldrich
Triton X-100	BP151-100	Fisher Scientific
Tris Base	604204	Fisher Scientific
Sodium dodecyl sulfate solution	71736-100ML	Sigma Aldrich

1 [M]20 Milimolar (mM) Tris pH7.5

[M]50 Milimolar (mM) NaCl

[M]0.5 % (v/v) Triton-X100

[M]0.5 % (v/v) DOC

[M]0.5 % (v/v) SDS

[M]1 Milimolar (mM) EDTA

2 On ice add protease and phosphatase inhibitors (Roche)