

Dec 16, 2021

# Estimate phospholipids from microalgae

Yingyu Hu<sup>1</sup><sup>1</sup>Dalhousie University

1



protocol .

Marine Microbial Macroecology Lab

Tech. support email: [ruby.hu@dal.ca](mailto:ruby.hu@dal.ca)

Yingyu Hu

Dalhousie University

Here we describe a protocol to estimate phospholipids from microalgae. After the extraction and measurement of total lipids from microalgae, the rest of lipids extract is dried by nitrogen flow and then dried with magnesium sulphate at 90 °C. A 650 °C dry combustion is used to decompose 80% phospholipids. The ash is digested by 0.5 mL 0.2 M HCl for 30 minutes at 90°C. The resulting orthophosphate is detected by mixing the digested sample with a mixture of molybdate and ascorbic acid to produce molybdenum blue (Chen 1956).

P.S. Chen, T.Y. Toribara and Huber Warner. Microdetermination of Phosphorus. Anal. Chem..

<https://doi.org/10.1021/ac60119a033>

Yingyu Hu 2021. Estimate phospholipids from microalgae. **protocols.io**

<https://protocols.io/view/estimate-phospholipids-from-microalgae-b2yxqfxn>



phospholipids, high temperature dry combustion

protocol ,

Dec 16, 2021

Dec 16, 2021

## Prepare phospholipids sample

- 1 Dry remaining organic phase extract of total lipids at  $\delta$  **37 °C** under a stream of N<sub>2</sub> gas (<2 psi)

## High temperature dry combust

9h

- 2 Use diamond pen to engrave the sample vials with numbers. Log number and sample code.
- 3 **[M]0.17 M** MgSO<sub>4</sub> reagent:  
Dissolve **1.023 g** MgSO<sub>4</sub> in 50 mL MilliQ water  
**☒ Magnesium sulfate anhydrous Fisher**  
**Scientific Catalog #M65500**
- 4 Add **200 µL** **[M]0.17 M** MgSO<sub>4</sub> to the dry extract.

Sing-use pipet tip to avoid cross-contamination.

- 5 Cover the uncapped vials with foil and place in the oven at  $\delta$  **90 °C** until samples are completely dry.

Forced air oven

VWR 89511-410

Remove samples out of the oven as soon as they are dried. If muffle furnace is not available, keep samples in vacuum desiccator.

6 Combust dried samples at 🔥 **650 °C** for ⌚ **09:00:00**

9h

Muffle furnace  
F30428C  
Thermo 10-505-13

Only place glass vials in the muffle furnace. Foil can be combusted under 🔥 **650 °C**

7 Allow samples to gradually cool down in the muffle furnace.

#### Digestion

8 [M]**0.2 M** HCl reagent:

In a reagent bottle, dissolve one part of [M]**12 N** HCl in 59 parts of MilliQ water

🔗 **12 N Hydrochloric acid Contributed by users**

Volume of HCl\_0.2M\_mL = (0.5\_mL) X (#Sample + #Blank)

9 Preheat oven to 🔥 **90 °C**

10 Add 📄 **0.5 mL** [M]**0.2 M** HCl to each vial.

11 Tightly cap the vial and vortex.

12 Place vials in the oven for 🕒 00:30:00

30m

13 Cool samples down to 🌡 Room temperature

Preparing standard working solutions

2h

14 Standard working solutions and reagents can be prepared during sample digestion.

15  $\text{KH}_2\text{PO}_4$  primary standard stock solution ( $\approx 1 \text{ mM}$ )

🔗 Potassium dihydrogen orthophosphate ACP

Chemicals Catalog #P-4550

15.1 Transfer about 1 g  $\text{KH}_2\text{PO}_4$  into a beaker, cover the beaker with foil

15.2 Place the beaker into an oven, dry  $\text{KH}_2\text{PO}_4$  at 🌡 110 °C for at least  
🕒 02:00:00

2h

15.3 Move  $\text{KH}_2\text{PO}_4$  into a vacuum desiccator, allow  $\text{KH}_2\text{PO}_4$  to cool to room temperature

15.4 Dissolve around 📏 0.136 g dried  $\text{KH}_2\text{PO}_4$  in 📏 1 L MilliQ water.  
■ Use 1 L volumetric flask  
■ Take notes of the actual weight of  $\text{KH}_2\text{PO}_4$  for final concentration of standard stock solution

15.5 Transfer standard stock solution into a 1 L bottle and store in the fridge.

This stock solution lasts quite a long time, unless there is evidence for growth of algae or other extraneous biotic material.

## 16 Standard working solution


KH <sub>2</sub> PO <sub>4</sub>	Primary (ul)	MilliQ (ul)
S1	0	1000
S2	5	995
S3	10	990
S4	20	980
S5	50	950
S6	100	900
S7	150	850
S8	200	800

17 Transfer  500 µL of each standard working solution to 2 mL microtube.

### Preparing working reagents

## 18

All reagents are freshly prepared before colorimetric measurement.

19  6 N (3 M) sulfuric acid reagent:

Carefully add 1 part  18 M concentrated sulfuric acid into 5 part MilliQ water

 18M sulfuric acid Contributed by users

## 20 [M]2.5 % ammonium molybdate reagent:

Weigh **0.25 g** ammonium molybdate in a Falcon tube and top to **10 g** with MilliQ water.  
Cap and shake until totally dissolved.

 **Ammonium molybdate Sigma**

**Aldrich Catalog #09878-100G**

## 21 [M]10 % ascorbic acid reagent:

Weigh **1 g** ascorbic acid in a Falcon tube and top to **10 g** with MilliQ water;  
Cap and shake until all dissolved.

 **Ascorbic acid Sigma**

**Aldrich Catalog #A5960-100G**

Wrap the tube with foil if the reagent is not used right after prepared.

## 22 Calculate the volume of molybdate-ascorbic reagent:

Total volume of reagent\_mL = (0.5 mL) X (#standard working solution + #samples + #blanks)


## 23 Mix the reagents into Falcon tube:

A	B
Reagent	Part(s) as in volume
MilliQ	2
6N sulphuric acid	1
2.5% ammonium molybdate	1
10% ascorbic acid	1

Colorimetric measurement 3h

## 24 Preheat incubator/shaker to **37 °C**

SHAKING INCUBATOR  
71L  
Corning® LSE™ 6753

- 25 Add  **500 µL** reagent to each standard, sample (in the vial) and blank, starting from blanks, including blank for standards and blank for samples.

Finntip Stepper Tips  
5 mL  
Thermo Scientific 9404200

Before dispensing the reagent, wipe or dab the liquid drop on the outside of the tip, avoid wiping the open tip.

- 26 Vortex.

- 27 Incubate at  **37 °C** for  **03:00:00** while shaking at 150 rpm

3h

- 28 Load microplate with 250 ul reactant from each tube, duplicate.

	<u>1</u>	<u>2</u>	<u>3</u>	<u>4</u>	<u>5</u>	<u>6</u>	<u>7</u>	<u>8</u>	<u>9</u>	<u>10</u>	<u>11</u>	<u>12</u>
<u>A</u>	S1	S1	Samples and sample blanks: 40 with duplicate									
<u>B</u>	S2	S2										
<u>C</u>	S3	S3										
<u>D</u>	S4	S4										
<u>E</u>	S5	S5										
<u>F</u>	S6	S6										
<u>G</u>	S7	S7										
<u>H</u>	S8	S8										

Example of loading the microplate

96-Well Microplates, Polystyrene, Clear,  
Greiner Bio-One 655101

## 29 Read plate in microplate reader

A	B
Shake duration	00:00:05
Shaking type	Continuous
Shaking force	High
Shaking speed [rpm]	600
Wavelength [nm]	820
Use transmittance	No
Pathlength correction	No
Measurement Time [ms]	100


Varioskan LUX Multimode Microplate  
Reader  
Thermo Fisher VL0L00D0



- 30 Subtract the average absorbance at 820 nm of the blank standard replicates from the absorbance at 820 nm of all other standard working solutions.
- 31 Subtract the average absorbance at 820 nm of the blank sample (i.e. blank filter) replicates from the absorbance at 820 nm of all other individual samples.
- 32 Prepare a standard curve by plotting the average blank-corrected 820 nm absorbance for each standard working solution versus its concentration in uM.
- 33 Use the standard curve to determine the orthophosphate concentration of each unknown sample by using its blank-corrected 820 nm absorbance.
- 34  $(P_{\text{measured}})_{\text{ug/sample}} = (\text{orthophosphate})_{\text{uM}} \times (V_{\text{HCl}})_{\text{mL}} \times (0.001) \times (30.97)$

$$(P_{\text{corrected}})_{\text{ug/sample}} = (P_{\text{measured}}) / 0.8$$

Where, 0.8 is the average recovery of phospholipids after a high temperature dry combustion at **650 °C**.



35 (Phospholipids)<sub>ug/sample</sub> = (P<sub>corrected</sub>)/(0.01X4.3)