



VERSION 2

MAR 07, 2023

OPEN  ACCESS

**DOI:**  
[dx.doi.org/10.17504/protocols.io.ewov1ojd2lr2/v2](https://dx.doi.org/10.17504/protocols.io.ewov1ojd2lr2/v2)

**Protocol Citation:** Edwin Yoo 2023. Green Lab Nanoparticle For 6 Well Cardiomyocyte Transfection v2. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.ewov1ojd2lr2/v2> Version created by Edwin Yoo

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Mar 07, 2023

**Last Modified:** Mar 07, 2023

**PROTOCOL integer ID:**  
78272

## Green Lab Nanoparticle For 6 Well Cardiomyocyte Transfection v2 V.2

Edwin Yoo<sup>1</sup><sup>1</sup>johns hopkins

Edwin Yoo

### ABSTRACT

nanoparticle transfection

### MATERIALS

Sodium acetate 25mM pH5

Polymer stock (100 ug/ul)

Desired media

## Nanoparticle Synthesis

## 1 Defining Terms

- DNA refers to nucleic acid cargo encoding genes of interest, generally stored at a concentration of 1 ug/ul.
- Polymer refers to the [Poly\(beta-amino ester\)](#) (PBAE) that makes up the backbone of the nanoparticle.
- Nanoparticle refers to the final complexed polymer + DNA that is ready for transfection.

## 2 DNA Dilution

For 3 technical replicate wells of one condition in 6-well plate format, mix the following in microcentrifuge tube.

- 60 ul DNA at 1ug/ul
- 440 ul sodium acetate
- Total volume = 500 ul for 3 technical replicate wells per one condition

## 3 Polymer Dilution (7.2 mg.ml polymer concentration in NaAc)

For 3 technical replicate wells of one condition in 6-well plate format, mix the following in microcentrifuge tube.

- 40 ul polymer (e.g. 4-5-6)
- 460 ul sodium acetate
- Total volume = 500 ul for 3 technical replicate wells per one condition

## 4 Nanoparticle (DNA+polymer) Synthesis

Combine diluted DNA with diluted polymer

- Add 500 ul of diluted DNA to 500 ul diluted polymer (can freeze here at -80 at your own risk)
- Add 5000 ul of serum containing media to 1000 ul NP's for final volume master mix of 6000 ul per condition
- Incubate for 10 minutes

## 5 Transfection

- Change media from each well of 6 well plate with 2000 ul nanoparticle master mix (media + nanoparticles)
- Incubate for 24 hours or desired time
- Change media after desired transfection time.