

DEC 11, 2022

WORKS FOR ME

2

Expansion-assisted iterative fluorescence in situ hybridization (EASI-FISH) in *Drosophila* CNS

DOI

dx.doi.org/10.17504/protocols.io.5jyl8jmw7g2w/v1Mark Eddison¹, Gudrun Ihrke¹¹Project Technical Resources, Janelia Research Campus, Ashburn, VA, USA.Mark Eddison
[Janelia Research Campus](#)

COMMENTS 2

ABSTRACT

A straightforward, robust, and reliable protocol (EASI-FISH) that utilizes expansion microscopy and the hybridization chain reaction for multiplexed *in situ* hybridization in thick slices of the mouse brain has recently been described (Wang et al., 2021). Below details a modified version of the EASI-FISH protocol for adult *Drosophila* CNS, which includes antibody detection of fluorescent reporters. The protocol also works well for larval CNS and is expected to be applicable to other tissue types.

ATTACHMENTS

[Fly_EASI_FISH_publication_protocol.docx](#) [Materials_for_EASI_FISH.xlsx](#) [LIGHTSHEETHOLDER_V18R_and_V20M_DOCUMENTATION.zip](#)

DOI

dx.doi.org/10.17504/protocols.io.5jyl8jmw7g2w/v1

PROTOCOL CITATION

Mark Eddison, Gudrun Ihrke 2022. Expansion-assisted iterative fluorescence in situ hybridization (EASI-FISH) in *Drosophila* CNS. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5jyl8jmw7g2w/v1>

KEYWORDS

Drosophila CNS, In situ hybridization, HCR, Expansion Microscopy, EASI-FISH

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 28, 2022

LAST MODIFIED

Dec 11, 2022

OWNERSHIP HISTORY

Sep 28, 2022



maria.s

Oct 07, 2022



Mark Eddison Janelia Research Campus

PROTOCOL INTEGER ID

70604

GUIDELINES



Reference:

Wang et al., EASI-FISH for thick tissue defines lateral hypothalamus spatio-molecular organization. *Cell* 184 (2021) 6361-6377.e24.






<https://pubmed.ncbi.nlm.nih.gov/34875226/>





MATERIALS TEXT

Solutions for Fly brain dissection





- 2%  PARAFORMALDEHYDE 16% Aqueous SOL. EM GRADE **Fisher Scientific Catalog #15710** in S2 medium.
-  Schneider's Drosophila Medium **Thermo Fisher Catalog #21720024**
- PBT (0.5% Triton).
- 70% Ethanol (in nuclease-free water).

Solutions for Day 1: Labelling RNA


-  MOPS (Fine White Crystals/Molecular Biology) **Fisher Scientific Catalog #bp308100** . Store  Room temperature .
- Melphalan-X (Stock  2 mg/mL ; working  1 mg/mL). Store  -20 °C .

-  Acryloyl-X, SE **Thermo Fisher Scientific Catalog #A20770** , Stock $[1\text{M}]$ 10 mg/mL ; working $[1\text{M}]$ 0.1 mg/mL . Store $-20\text{ }^{\circ}\text{C}$.
-  Press-to-Seal™ Silicone Isolator with Adhesive, eight wells, 9 mm diameter, 0.5 mm deep **Thermo Fisher Catalog #P24743**
-  Poly-L-Lysine solution 10ml **TEDELLA Catalog #18026** , 1.6 mL + $3.2\text{ }\mu\text{L}$  PHOTO-FLO 200 SOLUTION **Electron Microscopy Sciences Catalog #74257** . Store $4\text{ }^{\circ}\text{C}$.
- PBS-Triton (0.1%)
- 0.2 ml PCR tubes (USA Scientific; 1402-4700).
- RNase Away (Thermo Scientific; 7003).



Solutions for Day 2: Gelation and Proteinase K digestion

- Stock-X. Store $-20\text{ }^{\circ}\text{C}$.
- 10%  Ammonium Persulfate **Sigma Catalog #A3678** . Store $-20\text{ }^{\circ}\text{C}$.
- 10%  NNN'N'-Tetramethylethylenediamine **Sigma Aldrich Catalog #T22500** . Store $-20\text{ }^{\circ}\text{C}$.
- 0.5%  4-Hydroxy-TEMPO **Sigma Catalog #176141**). Store $-20\text{ }^{\circ}\text{C}$.
-  Proteinase K, Molecular Biology Grade - 2 ml **New England Biolabs Catalog #P8107S** , $[1\text{M}]$ 800 U/ml). Store $-20\text{ }^{\circ}\text{C}$.
- $[1\text{M}]$ 50 millimolar (mM) Proteinase K Buffer: Store Room temperature .
- Small paint brush.

Solutions for Day 3: Hybridization

- Hybridization Buffer (Molecular Instruments, Store $-20\text{ }^{\circ}\text{C}$).
- Probe Wash Buffer (Molecular Instruments, Store $-20\text{ }^{\circ}\text{C}$).
- DNA oligo probes (Designed and made by Molecular Instruments, Store $-20\text{ }^{\circ}\text{C}$).
- DAPI (Sigma D9534)/  PBS **Fisher Scientific Catalog #BP24384** at $[1\text{M}]$ 500 ng/ml .

Solutions for Day 5: Hybridization Chain Reaction (HCR)

- Amplification Buffer (Molecular Instruments, Store 4 °C).
- Fluorescent Hairpins (448 or 546 from Molecular Instruments or 669 conjugated in lab, Store -20 °C).
- 5x SSCT (5x SSC, 0.1% Tween in Nuclease Free Water). Store Room temperature .
- 0.5x SSCT (0.5x SSC, 0.1% Tween in Nuclease Free Water). Store Room temperature .
-  Invitrogen GFP Polyclonal Antibody Alexa Fluor™ 488 **Thermo Fisher Scientific Catalog #A-21311** .
-  UltraPure™ BSA (50 mg/mL) **Thermo Fisher Scientific Catalog # AM2616** .

Solutions for Stripping Probes and Hairpins for Multiplexing

 DNase I, RNase-free **Qiagen Catalog #79254** .




Note

DNase1 buffer

A	B
Tris-HCL pH 8.0	10 mM
MgCl ₂	2.5 mM
CaCl ₂	0.5 mM

Recipes and Reagents



200mM MOPS Buffer (10X Stock):

 1046.5 mg in  25 mL in NFW, pH to 7.7 with  10 Normality (N) NaOH. Store -20 °C .

Melphalan stock (2.5mg/ml)

 Melphalan **Cayman Chemical Company Catalog #16665**

 DMSO, Anhydrous **Thermo Fisher Catalog #D12345**

- Dissolve  2.5 mg per ml in anhydrous DMSO (Invitrogen, D12345).
- Dissolve, heat to  37 °C and vortex vigorously and place on a shaker.
- May take an hour to dissolve.

- Aliquot in 800 μL batches.
- Store in a desiccated environment at $-20\text{ }^{\circ}\text{C}$.

Acryloyl-X (AcX) stock (10 mg/ml)

⊗ Acryloyl-X, SE Thermo Fisher Scientific Catalog #A20770

- Dissolve in anhydrous DMSO.
- Aliquot in 200 μL batches.
- Any extra aliquot in 5 μL batches (for extra AcX).
- Store in a desiccated environment at $-20\text{ }^{\circ}\text{C}$.
- Don't re-use AcX after thawing.

Melphalan-X (2mg/ml):

- Combine an equal concentration of Acryloyl-X (10 mg/mL) and Melphalan (2.5 mg/mL) (1-part AcX to 4-parts Melphalan (ie. 200 μL : 800 μL)).
- Incubate Overnight at Room temperature with shaking.
- Store in 50 μL aliquots in a desiccated environment at $-20\text{ }^{\circ}\text{C}$.
- Use at 1 mg/mL by 1:1 dilution with 20 millimolar (mM) MOPS.

Stock X:

4.04 Molarity (M) Sodium Acrylate* – made from Acrylic Acid as it is made at variable purity

⊗ 40% Acrylamide Solution Bio-rad Laboratories Catalog #1610140

⊗ 2% Bis Solution Bio-rad Laboratories Catalog #1610142

⊗ NaCl (5 M) RNase-free Thermo Fisher Scientific Catalog #AM9760G

⊗ Nuclease-free water Ambion Catalog #AM9932

10X PBS (ThermoFisher; AM9625)

4.04M Sodium Acrylate stock solution

- In a fume hood, place 5.5 mL of acrylic acid into a 50mL tube. Place tube in a Room temperature water bath (e.g. a beaker).
- Add 4.5 mL water.
- Add 7.2 mL 10 Molarity (M) NaOH gradually to prevent excessive heating and evaporation/boiling.
- Remove tube from hood (at this point most of the acrylic acid has been converted to non-volatile sodium acrylate).
- Add 1 Molarity (M) NaOH (nominally 1 mL) gradually until the pH is between 7.5 and 8 using a pH meter, at

🌡 Room temperature . Do not use pH test strips.

- Add water up to a final volume of 🧪 20 mL .

Note

For 9.4 mls Stock X (not including APS +TEMED + 4HT):

A	B
4.04M Sodium Acrylate	2275 ul
40% Acrylamide	625 ul
2% N,N MethylBisacrylamide	750 ul
5M NaCl	4000 ul
10x PBS	1000 ul
Nuclease Free Water	750 ul

10% APS:

A	B
APS	100 mg
H2O	900 ul

10% Temed:

A	B
TEMED	100 ul
H2O	900 ul

0.5% 4HT:

A	B
4HT	5 mg
H2O	995 ul

Note

50mM ProK/SDS Buffer (50ml)

⊗ EDTA (0.5 M, pH 8.0, nuclease-free) **Thermo Fisher Scientific Catalog #AM9260G**


⊗ 10% SDS solution **Thermo Fisher Scientific Catalog #15553027**

A	B
Tris-HCL- pH8	50 mM
EDTA	1 mM
TritonX	0.5 %
NaCl	50 mM
SDS	0.3 %

A	B
1 M Tris> 50 mM	2.5 ml
10 % Triton> 0.5 %	2.5 ml
5 M NaCl> 50 mM	0.5 ml
0.5 M EDTA	100 ul
10% SDS > 0.3% SDS	1.5 ml
Nuclease Free Water	42.9 ml

Note

RNase-Free DNase1

- Add  550 µL DNAase1 Buffer to DNase1 powder. Mix.

- Aliquot  50 µL per PCR tube, store @  -20 °C .

DNase1 buffer (50ml)

A	B
1 M Tris-HCL > 10 mM	500 ul
1 M MgCl ₂ > 2.5 mM	125 ul
1 M CaCl ₂ > 0.5 mM	25 ul
Nuclease Free Water	49.350 ml

HCR Wash Buffers (50 ml)

5x SSCT:







 SSC, RNase-free, 20× **Ambion Catalog #AM9763**

A	B
20x RNase free-SSC (ThermoFisher, AM 976)	12.5 ml
10% Tween	500 ul
Nuclease Free Water	37 ml

0.5x SSCT:

A	B
20x RNase free-SSC	1.25 ml
10% Tween	500 ul
Nuclease Free Water	48.25 ml

Reagents for JF-669 conjugation to unlabelled hairpins:

-  Janelia Fluor® 669 **Tocris Catalog #6420** Store  -20 °C .
- [M] 1 nanomolar (nM) unlabelled amine-modified hairpins (~70mer, Molecular Instruments, Store  -20 °C).
-  Acetonitrile anhydrous **Thermo Fisher Scientific Catalog #042311-K7**
- [M] 0.1 Molarity (M) Sodium Bicarbonate pH 8-9.
-  DMSO, Anhydrous **Thermo Fisher Catalog #D12345**
-  QIAquick Nucleotide Removal Kit (250) **Qiagen Catalog #28306**

-  SCREW CAP MICROCENTRIFUGE TUBES SELF-STANDING AMBER (0.5 ml) **USA Scientific Catalog #1405-9707**

Equipment



Equipment	
Savant™ SpeedVac™ SPD120 Vacuum Concentrator and Kits	NAME
Vacuum Concentrator	TYPE
Savant	BRAND
SPD120-230	SKU
https://www.thermofisher.com/order/catalog/product/SPD120-230	LINK

Fly brain dissection

1

Dissect fly CNS in S2 medium.
(for more details see attached protocol)

2

Fix up to 20 brains or 10 CNS in  2 mL of 2% PFA/S2 medium for  00:55:00 in the dark on a nutator. 55m



3

Rinse sample 1 x  2 mL PBST (0.5% Triton).



4

Wash sample 4 x  00:15:00  1 mL PBST(0.5%Triton) on a nutator. 15m





5

Rinse sample 1 x  2 mL 70% EtOH.



6


Store brains in  2 mL 70% EtOH @  4 °C for up to 6 months.





Day 1: Labelling RNA

7



Transfer brains to a 0.2ml PCR tube (2-4 brains per tube).

8 Rehydrate with 2 x 5 min wash in  150 μ L PBT (0.1%).



8.1 Rehydrate with 2 x  00:05:00 wash in  150 μ L PBT (0.1%) (1/2).

5m

8.2 Rehydrate with 2 x  00:05:00 wash in  150 μ L PBT (0.1%) (2/2).

5m

9 Incubate brains 1 x  00:30:00 in  150 μ L  20 millimolar (mM) MOPS Buffer.

30m



10 Thaw Melphalan-X (MelphX) and Acryloyl-X (Ac-X) solutions.

11 Using a P20 pipette, take off as much MOPS buffer from the brains as possible.



12 Dilute Melphalan-X stock 1:1 with MOPS Buffer (to 1mg/ml).

13 Add 1/100 Ac-X ( 10 mg/mL) to Melphalan-X working solution. Vortex, mix, and spin.




14 Add  30 μL of Melphalan-X/AcX solution to each PCR tube and gently mix.




15 Incubate  Overnight @  37 $^{\circ}\text{C}$.





2m

16 Prepare gel chambers for gelation the next day. Wipe a non-charged slide with RNase away. Adhere a gasket (4-6 wells max) and coat the glass surface of the chamber with  1 μL poly-Lysine using a P20 pipette tip. Air dry and repeat.



 Chambers.jpg

Day 2: Gelation and Proteinase K digestion

17 Wash brains 2 x 2 min  150 μL PBT and 1 x 2 mins  150 μL PBS.




4m

17.1 Wash brains 2 x  00:02:00  150 μL PBT and 1 x  00:02:00  150 μL PBS.

17.2 Wash brains 2 x  00:02:00  150 μL PBT and 1 x  00:02:00  150 μL PBS.

4m

18 Thaw Stock-X and 4HT, Temed and APS. Vortex well and keep  On ice.



19 Gently stick down brains in the center of the chamber, once stuck down, carefully add a drop of PBS to prevent dehydration.



Note

20

Mix together Stock-X and 4HT, Temed and APS at a ratio of 94:2:2:2. Vortex.






Note

21

Remove PBS from chamber with pipette tip and carefully wick away remaining PBS with a tissue.

22

Pipette  40 μ L of gel solution, on top of the brain, to each chamber. Incubate slide in the fridge ( 4 °C) for  00:10:00 .

10m






23

Take off gel solution and repeat step 22.

Note


24

Take off the gel solution and gasket surface adhesive. Add a final  40 μ L of gel solution and gently place a cover slip over the chamber. Gently press to seal the coverslip and incubate @  4 °C for  00:10:00 .


10m



Note

- 25 Polymerize the gel @ 37 °C for 01:30:00 to 02:00:00 . 3h 30m
- 26 Cool gels on the bench for a few minutes.
- 27 Take off the chamber lid and gasket with a razor blade. Trim the gels into a rectangle and nick the top right-hand corner to track the orientation of the sample.
- 28 Take off the gel from the slide with a paintbrush that has been wetted with a small amount of ProK Buffer and transfer it to a 2 mL Eppendorf.
- 29 Incubate each gel with 1 mL ProK Buffer and 10 µL (1/100) ProK Enzyme @ 37 °C Overnight . 2h
- 

Day 3: Hybridization

- 30 Wash gels 4 x 15 min with 1 mL PBS at Room temperature . 
- 30.1 Wash gels 4 x 00:15:00 in 1 mL PBS at Room temperature (1/4). 15m

30.2 Wash gels 4 x  00:15:00 in  1 mL PBS at  Room temperature (2/4).

15m

30.3 Wash gels 4 x  00:15:00 in  1 mL PBS at  Room temperature (3/4).

15m

30.4 Wash gels 4 x  00:15:00 in  1 mL PBS at  Room temperature (4/4).

15m

31 DAPI stain gels for  00:10:00 with  1 mL DAPI/PBS (500 ng/ml).

10m






32 Rinse with PBS.



33 Use a dissection scope with a UV bulb to neatly trim gel edges with a razor blade.




34 Thaw and mix hybridization (hyb) and probe wash buffer. Make sure hyb buffer is clear.



35 Incubate gel in  500 μ L hyb buffer for  00:30:00 @  37 °C .

30m



36 Dilute probes 1/100 ( 10 Nanomolar (nM)), in  300 μ L hyb buffer per gel. Vortex. Incubate @  37 °C .





37 Incubate gels with probes overnight @ 37 °C , no shaking necessary.



38 Put probe wash buffer and PBS @ 37 °C .

Day 4: Probe Washing

39 Wash 3 x 30 min 750 µL Probe Wash Buffer @ 37 °C .



39.1 Wash 3 x 00:30:00 750 µL Probe Wash Buffer @ 37 °C (1/3).

30m

39.2 Wash 3 x 00:30:00 750 µL Probe Wash Buffer @ 37 °C (2/3).

30m

39.3 Wash 3 x 00:30:00 750 µL Probe Wash Buffer @ 37 °C (3/3).

30m

40 Wash 3 x 30 min 1 mL PBS @ 37 °C .



40.1 Wash 3 x 00:30:00 1 mL PBS @ 37 °C (1/3).

30m

40.2 Wash 3 x  00:30:00  1 mL PBS @  37 °C (2/3).

30m

40.3 Wash 3 x  00:30:00  1 mL PBS @  37 °C (3/3).

30m

41 Wash 3 x 1hr  1 mL PBS @  37 °C .



41.1 Wash 3 x  01:00:00  1 mL PBS @  37 °C (1/3).

1h

41.2 Wash 3 x  01:00:00  1 mL PBS @  37 °C (2/3).

1h

41.3 Wash 3 x  01:00:00  1 mL PBS @  37 °C (3/3).

1h

42 Keep gels in PBS at  Room temperature  Overnight .



1h

Day 5: Hybridization Chain Reaction (HCR)

43 Incubate gels in  500 µL Amplification buffer for at least  00:30:00 @  Room temperature



30m

44 Snap cool hairpins with PCR machine @  95 °C for  00:01:30 and cool @  Room temperature for



31m 30s



00:30:00

45 Mix hairpins h1 and h2 @ 1/100 in 300 µL Amp Buffer per gel. Vortex.



46 Incubate gel with hairpins for 03:00:00 @ Room temperature in the dark.



47 Wash gels 2 x 20 min in 750 µL 5X SSCT @ Room temperature



47.1 Wash gels 2 x 00:20:00 in 750 µL 5X SSCT @ Room temperature (1/2).

20m

47.2 Wash gels 2 x 00:20:00 in 750 µL 5X SSCT @ Room temperature (2/2).

20m

48 Wash gels 2 x 40 min in 1 mL 0.5X SSCT @ Room temperature



48.1 Wash gels 2 x 00:40:00 in 1 mL 0.5X SSCT @ Room temperature (1/2).

40m

48.2 Wash gels 2 x 00:40:00 in 1 mL 0.5X SSCT @ Room temperature (2/2).

40m

49 Stain sample with 500 µL of anti-GFP-488 Ab @ (1/500) in PBT (0.1%) containing 5 mg/mL

15m





Ultrapure BSA and incubate Overnight (or the weekend) @ 4 °C .

Day 6: Mount and Image

50



Wash 2 x 30 min with 1 mL PBS-Triton (0.1%).

50.1

Wash 2 x 00:30:00 with 1 mL PBS-Triton (0.1%) (1/2).

30m

50.2

Wash 2 x 00:30:00 with 1 mL PBS-Triton (0.1%) (2/2).

30m

51



Wash 2 x 30 min and 1 x 1hr with 1 mL PBS.

51.1

Wash 2 x 00:30:00 with 1 mL PBS (1/2).

30m

51.2

Wash 1 x 01:00:00 with 1 mL PBS (2/2).

1h

52


DAPI stain gels for 00:15:00 with 1 mL PBS/DAPI (500 ng/ml).

15m

53



Mount gels (sample up) on an 8mm poly-lysine coated coverslip superglued to a Z.1 light-sheet sample holder and image.

54 For gel removal from the holder, incubate gels in  750 μL of 10% Dextran Sulphate for 20 mins. Gels will shrink and fall off the coverslip.

55 For long-term storage, keep gels in  750 μL 10% Dextran Sulphate @  4 $^{\circ}\text{C}$



2h 45m

Stripping Probes and Hairpins for Multiplexing



30m

56 Incubate gel for  00:30:00 in  1 mL of DNase1 Buffer @  37 $^{\circ}\text{C}$.




57 Add  450 μL of DNase Buffer to  50 μL DNase1. Mix.





58 Incubate gel in DNase1 for  02:00:00 @  37 $^{\circ}\text{C}$.





2h

59 Wash 4 x 15 min with  1 mL PBS.





59.1 Wash 4 x  00:15:00 with  1 mL PBS (1/4).



15m

59.2 Wash 4 x  00:15:00 with  1 mL PBS (2/4).

15m

59.3 Wash 4 x  00:15:00 with  1 mL PBS (3/4).

15m

59.4 Wash 4 x  00:15:00 with  1 mL PBS (4/4).

15m

60 Hybridize with next round of probes (Day 3, step 34).



JF-669 conjugation to unlabelled hairpins

1h 45m






61

30m



Note

Turn on speed vacuum (eg. Thermofisher, SPD120) and defrost dye at  Room temperature for  00:30:00 .

62








 Resuspend  2 mg of JF-669, SE in  630 μ L Acetonitrile, mix and vortex, and aliquot in  30 μ L into labelled skirted 0.5ml screw-cap centrifuge tubes > each will contain  0.1 mg of dye.

63

Evaporate acetonitrile in speed vacuum (in organic solvent mode) for  00:45:00 . Can store @  -20 °C .

45m

64

 In two 1.5ml Eppendorf tubes, evaporate  5 μ L ( 500 picomolar (pM) /  10 μ g) of unlabelled hairpins h1 and h2 using a speed vac, in aqueous mode, for  00:30:00 . If you have  10 μ L ( 1 nanomolar (nM)) use two tubes per hairpin. Check they have been fully evaporated.

30m

65 Add 3 μL of 0.1 Molarity (M) Sodium Bicarbonate pH 8-9 to each evaporated hairpin. Mix.



66 Add 2 μL of anhydrous DMSO to 0.1 mg of dye. Mix.



67 Add 2 μL dye mix (100 μg) to each 3 μL of hairpin (10 μg). Mix. It will change colour.



68 Leave hairpin-dye mixture to react Overnight at Room temperature in the dark.



69 The next morning, add 5 μL nuclease-free water to bring to 10 μL .



70 Remove excess dye with a QIAquick Nucleotide removal kit (add 100 μL PN1).



71 Elute dye-oligo conjugate in 50 μL nuclease-free water.



72 Check the hairpin concentration on a spectrophotometer (eg. a NanoDrop One) and dilute to 60 ng/ μL .



73 Separately store hairpins h1-669 and h2-669 in  25 µL aliquot's in a PCR tube@  -20 °C .



74 Test conjugation by HCR.