



Jun 22, 2022

Wireless electrochemical measurement of kanamycin in whole blood

Jun-Chau Chien¹, Alexandra E. Rangel¹, Hyongsok T. Soh¹¹Stanford University

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.j8nlk5dydl5r/v1

SPARC

Tech. support email: info@neuinfo.org

Alexandra Rangel

ABSTRACT

This work evaluates the life-time of the structural-switching aptamer sensors in the whole blood using in vitro blood circulation setup as the first step toward implantable usage. The device consists of a 3-electrode sensing probe with the tip of the working electrode (WE) functionalized with the antibiotics aptamers. Square-wave voltammetry (SWV) is used to continuously track both the signal peak and the background current at 5 sec temporal resolution. The sensing probe features a diameter of 800um and is inserted through a 18G catheter. Two samples of whole blood (EDTA treated), spiked with either 500uM kanamycin or SSC 1x (salin sodium citrate) buffer, are alternated toward the sensor to evaluate the degradation of the aptamer sensitivity and the degree of drift.

DOI

dx.doi.org/10.17504/protocols.io.j8nlk5dydl5r/v1

PROTOCOL CITATION

Jun-Chau Chien, Alexandra E. Rangel, Hyongsok T. Soh 2022. Wireless electrochemical measurement of kanamycin in whole blood. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.j8nlk5dydl5r/v1>



KEYWORDS

aptamer, implant, kanamycin, electrochemical

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jan 16, 2019

LAST MODIFIED

Jun 22, 2022

PROTOCOL INTEGER ID

19401

MATERIALS TEXT

MATERIALS

 [Kanamycin Monosulfate 50 mg/mL Solution](#) **Gold**

Biotechnology Catalog #K-120-SL

 [TCEP-HCl](#) **Gold**

Biotechnology Catalog #TCEP

 [Human whole blood with Na₂EDTA](#) **Contributed by users**

 [Kanamycin aptamer](#) **Biosearch Technologies**

Aptamer sensor probe overview

- 1 The aptamer-immobilized sensing probe will be approximately 6 ~ 10 cm in length and ~1mm in diameter in order to fit into the 18G catheter. The sensing probe consists of 3 wires:
(1) A working electrode (WE) made of a gold wire (Au, Alfa Aesar Inc.) with its tip functionalized with aptamers at approximately 2 mm in length, as defined by the heat shrink tubing (McMaster Carr Inc.). The diameter of the gold wire is approximately 250um.
(2) A reference electrode (RE) made of silver chloride (AgCl, A-M Systems) with a diameter approximately 100 um (by immersing silver wires in bleach overnight).
(3) A counter electrode (CE) made of platinum (Pt, A-M Systems).

Aptamer Immobilization

- 2 Use heat-shrink and bare gold wire (250 um) to define the sensing area at the wire tip.
- 3 Solvent clean with sonication in acetone, methanol, and DI-water.
- 4 Perform electrochemical cleaning using 3-wire electrochemical system
 - a. 3 CV cycles in 200mM EC4 H₂SO₄
 - b. 3 CV cycles in 50mM EC3 H₂SO₄CV parameters: -0.4V ~ 1.5V, step = 0.001V, rate = 0.001V/sec

- 5 Prepare 100 micromolar stock of kanamycin aptamer (5'-HS-(CH₂)₆-GGGACTTGTTTAGGTAATGAGTCCC-(CH₂)₇-NH-MB-3'). Aptamer is functionalized with a thiol and methylene blue on the 5' and 3' termini respectively.
- 6 Prepare fresh TCEP solution at a concentration of 28.7 mg/mL in water.
- 7 Pipette **2 μL** TCEP solution in to an aptamer aliquot (**1 μL** at 100 micromolar concentration, stored at **-20 °C** upon receiving in dry format) and incubate for **00:40:00** in dark conditions.
- 8 Add **97 μL** of 1X SSC buffer into the aliquot after **00:40:00** time period.
- 9 Immerse the gold wire into the TCEP-aptamer-SSC solution, incubate for **01:00:00** in dark conditions.
- 10 After **01:00:00** take out the gold wire, rinse with DI-water for **00:01:00** and immerse in DI-water for **00:02:00**
- 11 Prepare C6-solution: **1.245 mL** of DI-water + **1 μL** of C6 solution.
- 12 Immerse the gold wire into the C6-solution, incubate for **02:00:00** .
- 13 After 2 hours, rinse the gold wire with DI-water for **00:01:00** , and immerse in DI-water for

- 14 Store the gold wire in 1X SSC buffer at 4°C overnight (24:00:00) to stabilize the formation of the self-assembly monolayer (SAM).

Aptamer real-time sensor measurement

- 15 Blood circulation loop (Figure 1 and Figure 2) is created with pump and tubing to emulate actual condition.




Figure 1: Schematic illustration of the blood circulation loop. The sensing probe, with the tip functionalized with aptamer sensors, is inserted into the tube through a 18G catheter. The exposed wires are connected to the potentiostat for electrochemical sensing.









Figure 2: Snapshot of the blood circulation loop setup in the lab.

16 Human whole blood sample ( 9 mL , stored at  4 °C upon receiving) is split into two separate reservoir, each with  4.5 mL . In the first reservoir,  0.5 mL of 5mM kanamycin in SSC buffer is added. The final concentration is therefore 500uM. In the second reservoir, 0.5mL of SSC buffer without kanamycin is added. Each reservoir has two opening for connecting to the two ports of the circulation loop.

17 Insert the sensing probe and connect to the electronics.

18 After the insertion of the aptamer probe, the pump is activated to circulate the blood. A  00:30:00 stabilization period is needed. During this period, the potential are swept based on the parameters below.

	SWV frequency	# of steps	Potential step size	Start potential	End potential	SWV amplitude
1 st SWV sweep	400 Hz	400	1mV	-0.4 V	0 V	36 mV
2 nd SWV sweep	60 Hz	80	5mV	-0.4 V	0 V	36 mV

19 PalmSense software is used to perform square-wave voltammetry (SWV) continuously. The parameters are listed below. Two SWV sweeps are performed sequentially at two different frequencies. The results are subtracted to mitigate the drift.

Parameters are aptamer-specific, i.e. different aptamers will need different sweeping parameters for minimum drift.

20 After each sweep, the measured current is stored in .CSV file. The results are post-processed in Matlab to extract: (1) background current, and (2) signal current, as shown in Figure 3.




Figure3. Example of SWV sweep. The signal of interest is the difference between the peak and the background (or baseline) current.

- 21 After 30-min stabilization period, the two samples are alternated using valves. Due to difference in the antibiotics concentration, the measurements resemble a square wave. The rise and fall time indicates the achievable temporal resolution, which is currently governed by the diffusion rate of the target across the porous membrane. Figure 4 shows an example of the measured results

Example Measurements




Figure 4: Measured time-domain results.