



Jul 12, 2020

# TAP agar plate preparation

Joao Vitor Molino<sup>1</sup><sup>1</sup>University of Zürich

1

Works for me

This protocol is published without a DOI.

Joao Vitor Molino  
University of Zürich

## PROTOCOL CITATION

Joao Vitor Molino 2020. TAP agar plate preparation. **protocols.io**  
<https://protocols.io/view/tap-agar-plate-preparation-big8kbzw>



## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jul 12, 2020

## LAST MODIFIED

Jul 12, 2020

## PROTOCOL INTEGER ID

39168

## DISCLAIMER:

### DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

## 1 Always follow aseptic techniques

1. Autoclavate freshly prepared TAP media with **1.5 Mass / % volume bacteriological agar** with a magnetic bar inside
2. Mix the flask using a magnetic stirrer to evenly dissolve the agar in the solution
3. For antibiotic containing plates, cool down the flask temperature until **55 °C** - **60 °C** or until it is possible to hold it without hurting the skin.
4. Add antibiotics in the desired amount, and mix using a magnetic stirrer to evenly dissolve the antibiotic
5. Pour the molten media in sterile petri dishes. (Keep volume constant among plates).
6. Wait for media to solidify
7. Stock at **4 °C** - **8 °C** until use. (Storage for 2 months was used without noticeable difference)