


Aug 15, 2022

# DNA Extraction from Yeast

Brian Teague<sup>1</sup><sup>1</sup>University of Wisconsin - Stout

1 Works for me



Share

This protocol is published without a DOI.

Yeast ORFans CURE



Brian Teague

University of Wisconsin - Stout

## ABSTRACT

This is a "quick and dirty" way to get some genomic DNA out of yeast cells. Not pure enough for many things, but should be fine for a PCR template.

This protocol is adapted from

Blount BA, Driessen MR, Ellis T (2016). GC Preps: Fast and Easy Extraction of Stable Yeast Genomic DNA.. Scientific reports.

<https://doi.org/10.1038/srep26863>

who in turn adapted it from

Lööke M, Kristjuhan K, Kristjuhan A (2011). Extraction of genomic DNA from yeasts for PCR-based applications.. BioTechniques.

<https://doi.org/10.2144/000113672>

## PROTOCOL CITATION

Brian Teague 2022. DNA Extraction from Yeast. **protocols.io**

<https://protocols.io/view/dna-extraction-from-yeast-cfagtibw>



## KEYWORDS

dna, extraction, saccharomyces, pcr, lioac, sds

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

#### IMAGE ATTRIBUTION

By gskx via Flickr. <https://www.flickr.com/photos/gskx/89462961>

#### CREATED

Aug 15, 2022

#### LAST MODIFIED

Aug 15, 2022

#### PROTOCOL INTEGER ID

68648

#### GUIDELINES

The centrifugation steps all specify 🌀**21000 x g** . If your microcentrifuge doesn't go this high, spin at the fastest speed available.

## MATERIALS TEXT

### Equipment

- Dry bath at 72°C
- Dry bath at 42°C (optional)

### Materials

- Sterile water

 [Lithium Acetate Dihydrate](#) **Sigma**

- **Aldrich Catalog #L4158** Step 1

solution,

[M] **1 Molarity (M)**

 [Sodium dodecyl sulfate](#) **Sigma**

- **Aldrich Catalog #436143-25G** Step 1

solution,

[M] **10 Mass / % volume**

- 

 [Ethanol \(100%, Molecular Biology Grade\)](#) **Fisher**

**Scientific Catalog #BP2818500** In 2 steps

- 

 [Ethanol \(100%, Molecular Biology Grade\)](#) **Fisher**

**Scientific Catalog #BP2818500** In 2 steps

solution, [M] **70 % (v/v)**

-  [TE Buffer](#) **Contributed by users** Step 10

## SAFETY WARNINGS

Both lithium acetate and SDS are irritants, particularly in the eyes. Wear appropriate PPE, including safety glasses, lab coats and gloves.

SDS is particularly gnarly if it's inhaled. If you're making a solution from powdered SDS, use a dust mask and/or weigh it out in a hood.

1 Make  **100 µL** of yeast lysis solution by mixing the following in a 1.7 ml microcentrifuge tube:

-  **70 µL** H2O

 [Lithium Acetate Dihydrate](#) **Sigma**

-  **20 µL** **Aldrich Catalog #L4158**

solution,

[M]1 Molarity (M)

☒ Sodium dodecyl sulfate Sigma

- 10 µL Aldrich Catalog #436143-25G

solution,

[M]10 Mass / % volume

2 Choose a yeast colony to analyze and circle it on the bottom of the petri dish. A large one is best.

3 Using a micropipette tip, scrape some of the colony off and resuspend it in the lysis solution. Vortex vigorously until the colony is mixed completely into the lysis solution.

**THIS IS NOT A CASE WHERE MORE IS BETTER.** Your solution should be slightly cloudy. If it is quite "thick", then try again.

4 Incubate at 70 °C for 00:05:00 5m

5 Add 300 µL of 3m

☒ Ethanol (100%, Molecular Biology Grade) Fisher

Scientific Catalog #BP2818500

and

vortex briefly. Centrifuge 21000 x g, 00:03:00 .

Make sure the centrifuge is balanced!

6 Using a P-1000 micropipettor, carefully aspirate the supernatant and discard as biological waste. **Do not disturb the pellet.**

7 Add 500 µL of [M]70 % (v/v) 3m

 Ethanol (100%, Molecular Biology Grade) Fisher



Scientific Catalog #BP2818500


to the

microcentrifuge tube. Centrifuge  **21000 x g, 00:03:00**.


Make sure the centrifuge is balanced!

- 8 Using a P-1000 micropipettor, carefully aspirate the supernatant and discard as biological waste. **Do not disturb the pellet. Try to get as much of the ethanol off as you can.**


- 9 Let the pellet dry by leaving it in a  **42 °C** dry bath for  **00:15:00**. Leave the cap open. 15m

If you don't have a dry bath, just leave the tube (cap open) at room temperature. Extend the time to  **00:30:00**

- 10 Add  **100 µL** of  **TE Buffer Contributed by users** and vortex to resuspend the pellet.

- 11 Centrifuge  **21000 x g, 00:00:30** to collect the cellular debris at the bottom. The DNA remains suspended in the supernatant. 30s

- 12 Label the tube and store at  **-20 °C**, or proceed directly to PCR.

- 13 If you need to use this sample again: 30s
- Thaw the sample completely.
  - Vortex briefly to resuspend everything.
  - Centrifuge  **21000 x g, 00:00:30** to (re)collect the cellular debris at the bottom of the tube.