



MAR 07, 2024

## ssDNA2.0: Klenow mix

Sarah Nagel<sup>1</sup>, Anna Schmidt<sup>1</sup>, Matthias Meyer<sup>1</sup>

<sup>1</sup>Max Planck Institute for Evolutionary Anthropology



Anya Patova

The Max Planck Institute for Evolutionary Anthropology

### ABSTRACT

Protocol for the preparation of Klenow mix for automated single-stranded DNA library preparation using the ssDNA2.0 method (Gansauge et al. 2020).

### References

Gansauge, M.-T., Aximu-Petri, A., Nagel, S., & Meyer, M. (2020). Manual and automated preparation of single-stranded DNA libraries for the sequencing of DNA from ancient biological remains and other sources of highly degraded DNA. *Nature Protocols*, 15, 2279-2300.

OPEN ACCESS



DOI:

[dx.doi.org/10.17504/protocols.io.ewov1qkbygr2/v1](https://dx.doi.org/10.17504/protocols.io.ewov1qkbygr2/v1)

**Document Citation:** Sarah Nagel, Anna Schmidt, Matthias Meyer 2024. ssDNA2.0: Klenow mix. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.ewov1qkbygr2/v1>

**License:** This is an open access document distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Created:** Jan 09, 2024

**Last Modified:** Mar 07, 2024

**DOCUMENT integer ID:** 93169

## Funders Acknowledgement:

Max Planck Society

Grant ID: -

## Note

The volume of Klenow mix suffices for one 96-well library preparation plate (96 + 20 reactions to account for dead volumes and loss of reagent). It is advisable to prepare 10-20 mixes at once.

## Materials

Reagent/Consumable	Supplier	Catalogue number	Decontamination*
Reagents			
Water, HPLC-grade	Sigma Aldrich/Merck	1153332500	UV
Tween-20 †	Thermo Fisher Scientific	11417160	UV
Klenow reaction buffer	Thermo Fisher Scientific	EP0052	-
Consumables			
5 ml screw cap tubes (rack 2d Lp W / barcode)	VWR	NUNC374320-BR	-

\* Decontamination of reagents should be performed as detailed in the Appendix.

† Use to prepare a 2% (vol/vol) solution in water. NOTE: Tween-20 is highly viscous, pipette slowly and with care.

## Equipment

- Label printer (e.g. Brady M611, cat. no. M611-EU-LABS) and tube labels (e.g. Labels for TLS2200/TLS PC Link/Polyester, cat. no. PTL-82-499)

## Protocol

1. Prepare the Klenow mix in a 5 ml screw-cap tube by combining the following reagents. Mix thoroughly by vortexing. Spin tube briefly in a microcentrifuge.

Reagent	Volume (µl)	Final concentration in reaction
Water	945	
Klenow reaction buffer (10x)	140	1x

Reagent	Volume (μl)	Final concentration in reaction
2% Tween-20 (v/v)	35	0.05%
sum	1120	

#### Note

##### [Labeling]

Prepare tube labels using Brady printer including name of the mix, date (dd.mm.yyyy) and the name of the person who prepared the Klenow Mix.

2. Freeze at -20 °C until used.

#### Note

##### [Documentation]

Note the lot/batch numbers of the reagents used for master mix preparation in Labfolder (orange fields).

## Appendix

#### Document



NAME

UV decontamination of reagents/buffers

CREATED BY

Elena Essel

PREVIEW