



APR 04, 2024

## 🌐 Rotarod Test

Sabina Marciano<sup>1</sup>, Roberta Marongiu<sup>2</sup>

<sup>1</sup>Weill Cornell Medicine; <sup>2</sup>Weill Cornell Medical College

ASAP Collaborative Research Network



Eileen Ruth Torres  
Weill Cornell Medicine

### ABSTRACT

Behavioral assay to measure motor impairments, specifically balance and coordination, in mice.

### MATERIALS

- MedAssociates mouse rotarod equipment
- computer with rotarod software
- Papertowels/kimwipes
- ethanol
- water
- timer

OPEN  ACCESS



DOI:

[dx.doi.org/10.17504/protocols.io.8epv5r6qjg1b/v1](https://dx.doi.org/10.17504/protocols.io.8epv5r6qjg1b/v1)

**Protocol Citation:** Sabina Marciano, Roberta Marongiu 2024. Rotarod Test. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.8epv5r6qjg1b/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working  
We use this protocol and it's working

**Created:** Apr 04, 2024

Last Modified: Apr 04, 2024

PROTOCOL integer ID: 97805

**Keywords:** ASAPCRN, behavior,  
mouse

**Funders Acknowledgement:**

Aligning Science Across  
Parkinson's  
Grant ID: 020608

- 1 Set up rotarod equipment as needed and clean thoroughly.
- 2 Turn on rotarod using green button and log into computer. Open software directly from desktop once machine is ON and connected.
- 3 Specify settings:
  - a. Protocol: Increasing speed
  - b. Speed: 4.0- 40 RPM
  - c. File name [EXP NAME].txt
    - i. Important: save file as .txt to then convert to excel once saved.
  - d. Subject ID: animal number
- 4 Place mouse on rotarod—can place up to five animals at a time. Make sure mice are all facing the same direction (i.e. wall, AWAY from experimenter) before starting.
- 5 Press RESET to start. IMPORTANT. IF YOU DO NOT PRESS RESET, THE SPEED SETTINGS WON'T APPLY.
- 6 Record behavior for 5 minutes or until mouse falls. If mouse falls <5s, place back on rod, reset zone using button, and continue recording. If mouse stays static (not walking/running but gripping rod) and spins more

than three times, remove and mark as failed trial.

a. IMPORTANT: The sensors register when a mouse falls and stops experiment for that zone. You must remove the fallen mouse as quickly as possible so that it doesn't intercept another zone.

**7** Place mice back in home cage and leave to rest for one minute before starting next trial. Clean rod in the meantime (within cage, use only water and/or kimwipe. Between cages, clean rod with ethanol and water).

**8** Perform five 5-minute (300 second) trials per cage.